



Product Catalogue

2024

SUSPENSION SOLUTIONS FOR THE HEALTHCARE SECTOR

MADE OF NECESSITY AND DESIGNED TO LAST

JBMEDICO

Highlights

- + Suspension solutions, hygiene, drip and infusion stands
- + Various stainless steel dispensers
- + Design, development and production of new products
- + Assembly

CONTACTS

Production and Design
Jan Haurum Kristensen
Phone: +45 7020 5826
E-mail: jan@jbmedico.com

Purchasing and Warehouse:
Jesper A. Berg
Phone: +45 7020 5826
Mobile: +45 2868 5826
Area: Denmark and Exports
E-mail: jb@jbmedico.com

Marketing and Website:
Erik Them
Mobil: +45 4221 5826
E-mail: et@jbmedico.com



ENVIRONMENTAL OBJECTIVES

JB Medico develops and produces high quality products with strong materials to increase the lifetime of the products. In this way we spare the environment in the long term.

Having a Green Environmental Profile and constantly focusing on sustainability means a lot to us and is something we deal with on a daily basis and constantly try to improve. We think we have come a long way, but there are of course places where we can improve even more.

Packaging:

We pack a lot of parcels every day and therefore use a lot of packaging in the form of boxes and stuffing. So this is one of the areas where we focus a lot on making it as sustainable as possible. We like to recycle cardboard boxes and fillings that still have a certain quality. So when our customers receive a box that looks used or has a logo from someone other than us, it's simply an attempt to protect the environment as much as possible.

FSC Eco-certified cardboard and paper from responsible sources. (www.fsc.org)
Minimal use of Plastic for box filling. Otherwise it is recycling from received shipments.
Recycling of cardboard boxes and filling.

Waste sorting:

We recycle for example:

- Paper and Cardboard.
- Packaging.
- Soft Plastics.
- Hard Plastic.
- Metal.
- Wiring.

Our hope for our customers:

We would be delighted if you would help us to take care of our planet by sorting and returning the packaging you receive from us to the recycling centre. Thank you!

Thanks in advance!

With best regards
Team JB Medico

CONTACT INFO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

TABLE OF CONTENTS

CONTACTS	2	SHARPS CONTAINER HOLDER - ROUND	64 - 65
ENVIRONMENTAL OBJECTIVES	2	SHARPS CONTAINER HOLDER - OVAL	65 - 66
TABLE OF CONTENTS	3	SHARPS CONTAINER HOLDER - SQUARED	66 - 68
MEDICAL RAIL CLAMPS and SLIDING CLAMPS	4 - 7	SHARPS CONTAINER HOLDER - ROUND	68 - 69
RAIL CLAMPS FOR LONG AND HEAVY EXAMINATION LIGHTS	8	WIRE BASKETS AND WIRE PRODUCTS	70 - 71
RAIL CLAMPS - EXTRA SOLID - FOR LONG OG HEAVY EXAMINATION LIGHTS	9	WIRE BASKETS AND WIRE ACCESSORIES	71
RAIL CLAMPS FOR DRIP STANDS, IT AND MEDICAL EQUIPMENT	10 - 11	WALL BRACKETS FOR BED RACKS	71
MEDICAL RAIL CLAMPS FOR EQUIPMENT RAILS	12	BRACKET FOR SUCTION CANISTER AND CATHETER CONTAINER	72
RAIL CLAMPS	13 - 14	VARIOUS ACCESSORIES	73 - 77
LAMP RAIL CLAMPS - EXTRA SOLID - FOR LONG AND HEAVY LAMPS	15	POWER TOWERS / CHARGING TOWERS	75
SUSPENSION RAIL CLAMPS FOR DROP STANDS, IT AND MEDICAL EQUIPMENT	15 - 16	EMERGENCY PULL CORD + ACCESSORIES	78
ACCESSORIES FOR VARIOUS LAMP RAIL CLAMPS	17	TERMS OF SALE AND DELIVERY	79
MULTIBRACKET - ADJUSTABLE - 16-41MM.	18 - 19	FAVORITES	80
SUSPENSION HOUSING FOR RAIL CLAMPS AND MULTIBRACKET	20	NOTES	80
FITTINGS SUITABLE FOR GCX RAILS	21 - 23		
MEDICAL EQUIPMENT RAILS AND ACCESSORIES	24		
MEDICAL EQUIPMENT RAILS DK AND EU DIN STANDARD	25		
ACCESSORIES FOR MEDICAL RAILS - SCANDINAVIAN STANDARDS	26		
WALL BRACKETS FOR BOTH RAIL STANDARDS	27		
WALL BRACKET FOR PROTECTION OF MEDICAL EQUIPMENT RAILS ETC.	28		
MEDICAL LIGHTING FROM GLAMOX LUXO	29 - 31		
ELECTRIC PILL CRUSHER - SEVERO V3	32		
ACCESSORIES FOR SEVERO	33		
MONITOR ARMS AND ACCESSORIES	34 - 38		
WALL BRACKETS	35		
ROD SUSPENSION	39		
SOAP AND ALCOHOL DISPENSERS	40 - 42		
DISPENSERS WITH DRIP TRAY AND ADAPTER BRACKET	43 - 44		
STAINLESS STEEL COVER FOR WIRE DISPERSERS	45		
TABLE AND WALL DISPENSER FOR SQUARE AND ROUND BOTTLES	46		
GLOVE, WET WIPE AND SANITARY FACEMASK HOLDER	47		
WALL DISPENSER FOR WET WIPES	48		
ACCESSORIES FOR GLOVE BOX HOLDERS	49		
ALUMINIUM SLIDING BRACKET FOR T-SLOT RAILS	50		
GLOVE BOX HOLDER	51		
IV STANDS / IV POLES	52 - 55		
IV POLE XL	53		
IV POLE XXL	54		
IV STAND / IV POLE - SUPPORT COLUMNS AND BASES	56		
IV STAND / IV POLE - ACCESSORIES	57 - 61		
WALL BRACKET FOR CO2 BOTTLES, Ø183 MM. HOLE	62		
OXYGEN AND GAS CYLINDERS HOLDERS	63		

CONTACT INFO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

MEDICAL RAIL CLAMPS and SLIDING CLAMPS

Scandinavian standard for 10x30 mm. medical equipment rails



Medical rail clamps for 10x30 mm. rails

All our medical rail clamps fit only Scandinavian standard 10x30 mm. rails.

To secure hospital equipment, a rail is used which is placed on a wall or bed panel.

When planning and ordering, it is important for you to know which of our systems you would like to use.

JB-Medico offers several different rail variants. We have developed two own rails 10x25 mm. and 10x30 mm., which can handle most of the cleats available on the market.

JB-Medico customises solutions for all clinical departments.

CONTACT INFO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

MEDICAL RAIL CLAMPS - SLIDING CLAMPS

Scandinavian standard for 10x30 mm. medical equipment rails



JB 88-00-00

Rail clamp, wide model, locked using two hex allen screws.



JB 103-00-00

Rail clamp, wide model, locked with one ball lock.



JB 104-00-00

Rail clamp, wide model, locked with two ball locks.



JB 207-00-00

Rail clamp, wide model, locked with two ball locks. It has three countersunk holes Ø6,6 mm.

"Prepared for suspension."

Item no: JB 88-00-00
Material: Aluminium (Natural Anodised)
Dimensions: Length x Width x Depth:
61,5 mm. x 90 mm. x 25 mm.
Weight: 280 grams
MOQ: 5 pcs.

Item no: JB 103-00-00
Material: Aluminium (Natural Anodized)
Dimensions: Length x Width x Depth:
61,5 mm. x 110 mm. x 25 mm.
Weight: 295 grams
MOQ: 5 pcs.

Item no: JB 104-00-00
Material: Aluminium (Natural Anodized)
Dimensions: Length x Width x Depth:
61,5 mm. x 110 mm. x 25 mm.
Weight: 304 grams
MOQ: 5 pcs.

Item no: JB 207-00-00
Material: Aluminium (Natural Anodized)
Dimensions: Length x Width x Depth:
61,5 mm. x 110 mm. x 25 mm.
Weight: 299 grams
MOQ: 5 pcs.

CONTACT INFO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

MEDICAL RAIL CLAMPS - SLIDING CLAMPS

Scandinavian standard for 10x30 mm. medical equipment rails

JBMEDICO



JB 206-00-00

Rail clamp, wide model, locked using two hex allen screws. It has three countersunk holes Ø6,6 mm.

"Prepared for suspension."



JB 217-00-00

Rail clamp, wide model, locked with one ball lock. It has three countersunk holes Ø6,6 mm.

"Prepared for suspension."



JB 120-00-00

Rail clamp, narrow model, locked with one ball lock, and a T-slot slot.



JB 121-00-00

Rail clamp, narrow model, locked with one ball lock.

Item no: JB 206-00-00
Material: Aluminium (Natural Anodized)
Dimensions: Length x Width x Depth:
61,5 mm. x 90 mm. x 25 mm.
Weight: 275 grams
MOQ: 5 pcs.

Item no: JB 217-00-00
Material: Aluminium (Natural Anodized)
Dimensions: Length x Width x Depth:
61,5 mm. x 110 mm. x 25 mm.
Weight: 293 grams
MOQ: 5 pcs.

Item no: JB 120-00-00
Material: Aluminium (Natural Anodized)
Dimensions: Length x Width x Depth:
61,5 mm. x 64 mm. x 25 mm.
Weight: 118 grams
MOQ: 5 pcs.

Item no: JB 121-00-00
Material: Aluminium (Natural Anodized)
Dimensions: Length x Width x Depth:
61,5 mm. x 64 mm. x 25 mm.
Weight: 155 grams
MOQ: 5 pcs.

CONTACT INFO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

MEDICAL RAIL CLAMPS - SLIDING CLAMPS

Scandinavian standard for 10x30 mm. medical equipment rails



JB 121-03-00

Rail clamp, narrow model, locked with one ball lock. It has three countersunk Ø6,6 mm. holes.

"Prepared for suspension."

Item no.: JB 121-03-00
Material: Aluminium (Natural Anodized)
Dimensions: Length x Width x Depth:
61,5 mm. x 64 mm. x 25 mm.
Weight: 150 grams
MOQ: 5 pcs.



JB 216-03-00

Rail clamp, narrow model, locked with two socket screws. It has three countersunk Ø6,6 mm. holes.

"Prepared for suspension."

Item no.: JB 216-03-00
Material: Aluminium (Natural Anodized)
Dimensions: Length x Width x Depth:
61,5 mm. x 45 mm. x 25 mm.
Weight: 132 grams
MOQ: 5 pcs.



JB 176-00-180

Rail Clamp, an extra wide model, is locked using four hex allen screws.

Item no.: JB 176-00-180
Material: Aluminium (Natural Anodized)
Dimensions: Length x Width x Depth:
61,5 mm. x 180 mm. x 25 mm.
Weight: 560 grams
MOQ: 4 pcs.



JB 176-00-180 supporting a transformer

CONTACT INFO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361



RAIL CLAMPS FOR LONG AND HEAVY EXAMINATION LIGHTS

Scandinavian standard for 10x30 mm.

Rail clamps for long, heavy and standard examination lights and lamps.

All our lamp rail clamps fit Scandinavian standard 10x30 mm. rails.

Can be fitted and removed without the use of tools, can be washed and disinfected in an instrument dishwasher when required.

For attachment of hospital equipment, a rail is used, which is placed on the wall or bed panel.

CONTACT INFO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

**RAIL CLAMPS - EXTRA SOLID -
FOR LONG OG HEAVY EXAMINATION LIGHTS**
Scandinavian standard for 10x30 mm. medical equipment rails



JB 206-03-18

Rail clamp, wide model, locked using two hex allen screws and adapter fitting Ø18 mm.

Three bush sizes are available.
Handles can be purchased.

Item no.: JB 206-03-18
Material: Aluminium (Ceramic Coating)
Dimensions: Length x Width x Depth:
61,5 mm. x 90 mm. x 60,6 mm.
Weight: 468 grams
MOQ: 2 pcs.



JB 217-03-18

Rail clamp, wide model, 1 pcs. ball lock and adapter fitting Ø18 mm.

Three bush sizes are available.
Handles can be purchased.

Item no.: JB 217-03-18
Material: Aluminium (Ceramic Coating)
Dimensions: Length x Width x Depth:
61,5 mm. x 110 mm. x 60,6 mm.
Weight: 475 grams
MOQ: 2 pcs.



JB 207-03-18

Rail clamp, wide model, 2 pcs. ball locks and adapter fitting Ø18 mm.

Three bush sizes are available.
Handles can be purchased.

Item no.: JB 207-03-18
Material: Aluminium (Ceramic Coating)
Dimensions: Length x Width x Depth:
61,5 mm. x 110 mm. x 60,6 mm.
Weight: 482 grams
MOQ: 2 pcs.



JB 216-03-18

Rail clamp, narrow model, locks with 2 pcs. hex allen screws and adapter fitting Ø18 mm.

Three bush sizes are available.
Handles can be purchased.

Item no.: JB 216-03-18
Material: Aluminium (Ceramic Coating)
Dimensions: Length x Width x Depth:
61,5 mm. x 45 mm. x 60,6 mm.
Weight: 316 grams
MOQ: 4 pcs.

CONTACT INFO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

RAIL CLAMPS FOR DRIP STANDS, IT AND MEDICAL EQUIPMENT

Scandinavian standard 10x30 mm.

JBMEDICO



Rail clamps for drip stands, IT and medical equipment

Can be used for:

- Drip stand on bedside panels
- Flat screen arms
- Medical equipment
- Medical lighting
- Rod suspension

JB-Medico offers several options for attaching hospital equipment to bedside panels or walls.

Different types of hospital equipment are attached to the rail by means of rail clamps. This provides a quick, secure and flexible attachment.

JB-Medico adapts solutions for all clinical departments.



JB 206-03-20

Rail clamp, wide model, locks with two hex allen screws and adapter bracket Ø20 mm.

Three bush sizes available.
Handles can be purchased (JB 94-00-00)

Item no.:	JB 206-03-20
Material:	Aluminium (Ceramic Coating)
Dimensions:	Lenght x Width x Depth: 61,5 mm. x 90 mm. x 60,5 mm.
Weight:	507 grams
MOQ:	2 pcs.

CONTACT INFO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

RAIL CLAMPS FOR DRIP STANDS, IT AND MEDICAL EQUIPMENT

Scandinavian standard 10x30 mm.

JBMEDICO



JB 217-03-20

Rail clamp, wide model, with a one-ball lock, adapter fitting and brass bush Ø20 mm. hole.

Three bush sizes available.
Handles can be purchased (JB 94-00-00)

Item no.: JB 217-03-20
Material: Aluminium (Ceramic Coating)
Dimensions: Length x Width x Depth:
61,5 mm. x 110 mm. x 60,5 mm.
Weight: 475 grams
MOQ: 2 pcs.



JB 207-03-20

Rail clamp, wide model, with a two-ball lock, adapter fitting and brass bush Ø20 mm. hole.

Three bush sizes available.
Handles can be purchased (JB 94-00-00)

Item no.: JB 207-03-20
Material: Aluminium (Ceramic Coating)
Dimensions: Length x Width x Depth:
61,5 mm. x 110 mm. x 60,5 mm.
Weight: 482 grams
MOQ: 2 pcs.



JB 216-03-20

Rail clamp, narrow model, with two hex allen screws, adapter fitting and brass bush Ø20 mm. hole.

Three bush sizes available.
Handles can be purchased (JB 94-00-00)

Item no.: JB 216-03-20
Material: Aluminium (Ceramic Coating)
Dimensions: Length x Width x Depth:
61,5 mm. x 45 mm. x 60,5 mm.
Weight: 316 grams
MOQ: 4 pcs.



JB 121-03-20

Rail clamp, narrow model, with a one-ball lock, adapter fitting and brass bush Ø20 mm. hole.

Three bush sizes available.
Handles can be purchased (JB 94-00-00)

Item no.: JB 121-03-20
Material: Aluminium (Ceramic Coating)
Dimensions: Length x Width x Depth:
61,5 mm. x 64 mm. x 60,5 mm.
Weight: 383 grams
MOQ: 4 pcs.

CONTACT INFO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

MEDICAL RAIL CLAMPS FOR EQUIPMENT RAILS

EU DIN standard 10x25 mm.



JB 146-00-00

Rail clamp, narrow model, with one ball lock and T-slot.



JB 147-00-00

Rail clamp narrow model, with one ball lock.



JB 147-03-00

Rail clamp, narrow model, with one ball lock and three countersunk Ø6,6 holes.

"Prepared for suspension."



JB 143-00-00

Rail clamp, wide model, locked using two hex allen screws.

Item no.: JB 146-00-00
Material: Aluminium (Ceramic Coating)
Dimensions: Height x Width x Depth:
61,5 mm. x 64 mm. x 28,6 mm.
Weight: 132 grams
MOQ: 5 pcs.

Item no.: JB 147-00-00
Material: Aluminium (Ceramic Coating)
Dimensions: Height x Width x Depth:
61,5 mm. x 64 mm. x 28,6 mm.
Weight: 168 grams
MOQ: 5 pcs.

Item no.: JB 147-03-00
Material: Aluminium (Ceramic Coating)
Dimensions: Height x Width x Depth:
61,5 mm. x 64 mm. x 28,6 mm.
Weight: 161 grams
MOQ: 5 pcs.

Item no.: JB 143-00-00
Material: Aluminium (Ceramic Coating)
Dimensions: Height x Width x Depth:
61,5 mm. x 90 mm. x 28,6 mm.
Weight: 308 grams
MOQ: 5 pcs.

CONTACT INFO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

RAIL CLAMPS

EU DIN standard for 10x25 mm. clamp rails



JB 143-03-00

Rail clamp, wide model, locked using two hex allen screws. It has three countersunk holes Ø6,6 mm.

"Prepared for suspension."

Item no.: JB 143-03-00
Material: Aluminium (Ceramic Coating)
Dimensions: Height x Width x Depth:
61,5 mm. x 90 mm. x 28,6 mm.
Weight: 303 grams
MOQ: 5 pcs.



JB 144-00-00

Rail clamp, wide model, locked with two ball locks.

Item no.: JB 144-00-00
Material: Aluminium (Ceramic Coating)
Dimensions: Height x Width x Depth:
61,5 mm. x 110 mm. x 28,6 mm.
Weight: 328 grams
MOQ: 5 pcs.



JB 144-03-00

Rail clamp, wide model, locked with two ball locks. It has three countersunk holes Ø6,6 mm.

"Prepared for suspension."

Item no.: JB 144-03-00
Material: Aluminium (Ceramic Coating)
Dimensions: Height x Width x Depth:
61,5 mm. x 110 mm. x 28,6 mm.
Weight: 318 grams
MOQ: 5 pcs.



JB 145-00-00

Rail clamp, wide model, locked with one ball lock.

Item no.: JB 145-00-00
Material: Aluminium (Ceramic Coating)
Dimensions: Height x Width x Depth:
61,5 mm. x 110 mm. x 28,6 mm.
Weight: 325 grams
MOQ: 5 pcs.

CONTACT INFO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

RAIL CLAMPS

EU DIN standard for 10x25 mm. clamp rails



JB 145-03-00

Rail clamp, wide model, locked with one ball lock. It has three countersunk holes Ø6,6 mm.

"Prepared for suspension."

Item no.: JB 145-03-00
Material: Aluminium (Ceramic Coating)
Dimensions: Height x Width x Depth:
61,5 mm. x 110 mm. x 28,6 mm.
Weight: 317 grams
MOQ: 5 pcs.



JB 045-03-00

Rail clamp, narrow model, locked with two socket screws. It has three countersunk Ø6,6 mm. holes.

"Prepared for suspension."

Item no.: JB 045-03-00
Material: Aluminium (Ceramic Coating)
Dimensions: Height x Width x Depth:
61,5 mm. x 45 mm. x 28,6 mm.
Weight: 151 grams
MOQ: 5 pcs.



JB 143-03-18

Rail clamp, wide model, locked using two hex allen screws and adapter fitting Ø18 mm.

Three bush sizes are available.
Handles can be purchased.

Item no.: JB 143-03-18
Material: Aluminium (Ceramic Coating)
Dimensions: Length x Width x Depth:
61,5 mm. x 90 mm. x 64,2 mm.
Weight: 468 grams
MOQ: 2 pcs.



JB 145-03-18

Rail clamp, wide model, 1 pcs. ball lock and adapter fitting Ø18 mm.

Three bush sizes are available.
Handles can be purchased.

Item no.: JB 145-03-18
Material: Aluminium (Ceramic Coating)
Dimensions: Length x Width x Depth:
61,5 mm. x 110 mm. x 64,2 mm.
Weight: 501 grams
MOQ: 2 pcs.

CONTACT INFO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

**LAMP RAIL CLAMPS - EXTRA SOLID -
FOR LONG AND HEAVY LAMPS
EU DIN STANDARD FOR 10X25**

**SUSPENSION RAIL CLAMPS FOR
DROP STANDS, IT AND MEDICAL EQUIPMENT
EU DIN STANDARD FOR 10X25 MM.**

JB MEDICO



JB 144-03-18

Rail clamp, wide model, 2 pcs. ball locks and adapter fitting Ø18 mm.

Three bush sizes are available.
Handles can be purchased (JB 94-00-00)



JB 147-03-18

Rail clamp, narrow model, 1 pcs. ball locks and adapter fitting Ø18 mm.

Three bush sizes available.
Handles can be purchased (JB 94-00-00)



JB 143-03-20

Rail clamp, wide model, locked using two hex allen screws and adapter fitting Ø18 mm.

Handles can be purchased (JB 94-00-00)



JB 145-03-20

Rail clamp, wide model, with a one-ball lock, adapter fitting and brass bush Ø20 mm. hole.

Handles can be purchased (JB 94-00-00)

Item no.: JB 144-03-18
Material: Aluminium (Ceramic Coating)
Dimensions: Length x Width x Depth:
61,5 mm. x 110 mm. x 64,2 mm.
Weight: 502 grams
MOQ: 2 pcs.

Item no.: JB 147-03-18
Material: Aluminium (Ceramic Coating)
Dimensions: Length x Width x Depth:
61,5 mm. x 64 mm. x 64,2 mm.
Weight: 343 grams
MOQ: 4 pcs.

Item no.: JB 143-03-20
Material: Aluminium (Ceramic Coating)
Dimensions: Length x Width x Depth:
61,5 mm. x 90 mm. x 60,6 mm.
Weight: 536 grams
MOQ: 2 pcs.

Item no.: JB 145-03-20
Material: Aluminium (Ceramic Coating)
Dimensions: Length x Width x Depth:
61,5 mm. x 110 mm. x 64,2 mm.
Weight: 550 grams
MOQ: 2 pcs.

CONTACT INFO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

**SUSPENSION RAIL CLAMPS FOR
DROP STANDS, IT AND MEDICAL EQUIPMENT
EU DIN STANDARD FOR 10X25 MM.**



JB 144-03-20

Rail clamp, wide model, with a two-ball lock, adapter fitting and brass bush Ø20 mm. hole.

Handles can be purchased (JB 94-00-00)

Item no.: JB 144-03-20
Material: Aluminium (Ceramic Coating)
Dimensions: Length x Width x Depth:
61,5 mm. x 110 mm. x 64,2 mm.

Weight: 551 grams
MOQ: 2 pcs.



JB 147-03-20

Rail clamp, narrow model, with a one-ball lock, adapter fitting and brass bush Ø20 mm. hole.

Handles can be purchased (JB 94-00-00)

Item no.: JB 147-03-20
Material: Aluminium (Ceramic Coating)
Dimensions: Length x Width x Depth:
61,5 mm. x 64 mm. x 64,2 mm.

Weight: 550 grams
MOQ: 2 pcs.



JB 045-03-20

Rail clamp, narrow model, with two hex allen screws, adapter fitting and brass bush Ø20 mm. hole.

Handles can be purchased (JB 94-00-00)

Item no.: JB 045-03-20
Material: Aluminium (Ceramic Coating)
Dimensions: Length x Width x Depth:
61,5 mm. x 45 mm. x 64,2 mm.

Weight: 384 grams
MOQ: 2 pcs.



CONTACT INFO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

ACCESSORIES FOR VARIOUS RAIL CLAMPS

JB 128-00-09

POM bushing with bushing hole Ø 9.6 mm.

Item no.: JB 128-00-09
Material: Plastic (White)
Weight: 12 grams
MOQ: 10 pcs.



JB 22-00-05

Brass bushing with bushing hole Ø 20 mm.

Item no.: JB 22-00-05
Material: Brass
Weight: 88 grams
MOQ: 10 pcs.



JB 128-00-01

POM bushing with bushing hole Ø 13 mm.

Item no.: JB 128-00-01
Material: Plastic (White)
Weight: 8 grams
MOQ: 10 pcs.



JB 94-00-00

Handle with M6 x 20 mm. thread

Item no.: JB 94-00-00
Material: Stainless steel (AISI 304)
Weight: 44 grams
MOQ: 5 pcs.



JB 127-00-01

POM bushing with bushing hole Ø 16 mm.

Item no.: JB 127-00-01
Material: Plastic (White)
Weight: 3,6 grams
MOQ: 10 pcs.



JB 46-00-00

Handle with M10 x 40 mm. thread

Item no.: JB 46-00-00
Material: Stainless steel (AISI 304)
Weight: 120 grams
MOQ: 2 pcs.



CONTACT INFO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

MULTIBRACKET - ADJUSTABLE - 16-41 MM.



JB 158-00-00

Multibracket - The fitting can be connected anywhere. Fits Ø16-41 mm.

This multibracket makes it possible to assemble various equipment together.

Item no.: JB 158-00-00
Material: Aluminium (Natural Anodized), SS 304, POM og Nylon (PGA)
Dimensions: Length x Width x Depth: 36 mm. x 51 mm. x 134/96 mm.
Weight: 209 grams
MOQ: 4 pcs.



JB 158-00-05

Multibracket with 2xØ6 clearance holes. Fits up to. Ø16-41 mm.

This multibracket makes it possible to assemble various equipment together.

Item no.: JB 158-00-05
Material: Aluminium (Natural Anodized), SS 304, POM og Nylon (PGA)
Dimensions: Length x Width x Depth: 36 mm. x 51 mm. x 134/96 mm.
Weight: 204 grams
MOQ: 4 pcs.



JB 158-06-05

Double multibracket for mounting of for example bed and drop stand.

Includes 2 pcs. Unbrako bolt screws.

Item no.: JB 158-06-05
Material: Aluminium (Natural Anodized), SS 304, POM og Nylon (PGA)
Weight: 421 grams
MOQ: 2 pcs.



JB 158-00-163

Multibracket with T-slot fitting fits from Ø16-41 mm. Fitting can be placed on table, trolley, drop stand, DIN rail and columns.

Includes 2 pcs. hex allen screws.

Item no.: JB 158-00-163
Material: Aluminium (Natural Anodized), SS 304, POM og Nylon (PGA)
Dimensions: Length x Width x Depth: 61,5 mm. x 64 mm. x 60,5 mm.
Weight: 245 grams
MOQ: 4 pcs.

CONTACT INFO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361



JB 158-03-18

Multibracket with a fixing device, Ø18mm. hole.
Fits Ø16-41 mm.

Three bush sizes available.
Handle can be purchased separately.

Item no.: JB 158-03-18
Material: Aluminium (Natural Anodized),
SS 304, POM, Nylon (PA6)
Dimensions: Length x Width x Depth:
36 mm. x 51 mm. x 134/96 mm.
Weight: 288 grams
MOQ: 2 pcs.

ACCESSORIES: Handle can be purchased.



JB 158-03-20

Multibracket with a fixing
device and brass bushing.
Fits up to. Ø16-41 mm.

Handle can be purchased separately.

Item no.: JB 158-03-20
Material: Aluminium (Natural Anodized),
SS 304, Brass, POM, Nylon (PA6)
Dimensions: Length x Width x Depth:
36 mm. x 51 mm. x 134/96 mm.
Weight: 440 grams
MOQ: 2 pcs.

ACCESSORIES: Handle and T-slot fitting can be
purchased.



JB 163-00-01

T-slot fitting for multi-bracket.

The holes are only suitable for multibracket.

Item no.: JB 163-00-01
Material: Aluminium (Natural Anodized)
Dimensions: Height x Width x Depth:
60 mm. x 36,2 mm. x 12,5 mm.
Weight: 36 grams
MOQ: 1 bag of 20 pcs.

CONTACT INFO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

SUSPENSION HOUSING FOR RAIL CLAMPS AND MULTIBRACKET



JB 205-00-00

Fixing device for rail clamps and multi-brackets. Hole dimension: Ø25 mm.

Prepared for suspension.

Item no.: JB 205-00-00
Material: Aluminium (Natural Anodized)
Dimensions: Length x Width x Depth:
36 mm. x 35,6 mm. x 61 mm.
Weight: 133 grams
MOQ: 2 pcs.

ACCESSORIES: Handle and bushing is purchased separatly.



JB 205-00-18

Fixing device for rail clamps and multi brackets. Hole dimension: Ø18 mm.

Prepared for suspension.

Item no.: JB 205-00-18
Material: Aluminium (Natural Anodized)
Dimensions: Length x Width x Depth:
36 mm. x 35,6 mm. x 61 mm.
Weight: 171 grams
MOQ: 2 pcs.

ACCESSORIES: Handle and bushing is purchased separatly.



CONTACT INFO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361



JB 820-625-06

GCX rail fitting, prepared with mounting holes.

Item no.: JB 820-625-06
Material: Aluminium (Natural Anodized)
Colour: Metal
Dimensions: Height x Width x Depth:
82 mm. x 63 mm. x 6 mm.
Weight: 83 grams
MOQ: 4 pcs.

ACCESSORIES:
JB 163-00-01 T-slot fitting
JB 205-00-00 Adapter fitting/ Suspension housing Ø20
JB 205-00-18 Adapter fitting/ suspension housing Ø18



JB 820-625-163

GCX rail fitting with a T-slot.

Item no.: JB 820-625-163
Material: Aluminium (Natural Anodized)
Colour: Metal
Dimensions: Height x Width x Depth:
82 mm. x 36 mm. x 18,5 mm.
Weight: 123 grams
MOQ: 4 pcs.



JB 820-625-20

GCX rail fitting with suspension housing. Ø20 mm. bushing hole.

Item no.: JB 820-625-20
Material: Aluminium (Natural Anodized) and Brass
Colour: Metal
Dimensions: Height x Width x Depth:
82 mm. x 63 mm. x 42 mm.
Weight: 313 grams
MOQ: 2 pcs.

ACCESSORIES: Handle can be purchased.



POM bushing not included as standard

JB 820-625-18

GCX rail fitting with suspension housing Ø18 mm. hole.

Optional plastic bushing Ø9,6/13 or 16 mm.

Item no.: JB 820-625-18
Material: Aluminium (Natural Anodized) and Plastic
Colour: Metal
Dimensions: Height x Width x Depth:
82 mm. x 63 mm. x 42 mm.
Weight: 270 grams
MOQ: 2 pcs.

ACCESSORIES: Handle can be purchased.

CONTACT INFO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361



JB 820-625-20
For flat screen / monitor holder.



JB 820-625-163
For catheter holder.



JB 820-625-18
For light fixture.



JB 820-625-163
For sharpsafe holder.

CONTACT INFO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361



JB 820-625-20
For monitor holder.



JB 820-625-163
For wire basket.



JB 820-625-20
Fitting with suspension housing.



JB 820-625-163
For catheter holder.

CONTACT INFO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

MEDICAL EQUIPMENT RAILS AND ACCESSORIES

Rail Scandinavian 10x30 mm. standard



JBM 100-01-06

JBM medical equipment rail in hygienic design.

With increased Depth to avoid collision with cable trays.
This allows suspension of medical equipment.

Item no.:	JBM 100-01-06
Material:	Aluminium (Natural Anodized),
Dimensions:	Length - Standard 6 meters Depth: 81 mm. Height: 93 mm.
Weight:	3,8 kg pr. meter
Standard length:	6 meters
End caps:	Plastic, easily mounted with wetroom silicone or screws
Screw holes:	Unplugged with plugs

Benefits include:

1. Very easy to clean
2. Optimal hygiene, the rail has smooth surfaces, no cracks or transitions that can be dirt collectors.
3. A steep slope, 25 degrees, so the rail cannot be used as a shelf.
4. Mounting holes are pre-drilled pr. 100 mm.



CONTACT INFO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

MEDICAL EQUIPMENT RAILS DK AND EU DIN STANDARD -
10X30 MM. AND 10X25 MM.

JBMEDICO



JB 252-00-00

10x25 mm. medical equipment rail.

EU DIN standard.



JB 249-00-00

End cap.

For EU DIN 10x25 mm. rail profile.



JB 234-00-01

10x30 mm. medical equipment rail.

Scandinavian standard.



JB 248-00-00

End cap.

For 10x30 mm. Scandinavian
standard rail profile.

Item no.: JB 252-00-00
Material: Aluminium
Surface: Natural Anodized
Colour: Metal
Dimensions: Height x Depth:
25 mm. x 10 mm.
Weight: 700 grams pr. m.
Standard Længde: 6 meters

Cutting:
Drilling Holes:

Item no.: JB 249-00-00
Material: Plastic (PE)
Colour: Grey
Dimensions: Height x Width:
25 mm. x 10 mm.
MOQ: 50 pcs.

Item no.: JB 234-00-01
Material: Aluminium
Surface: Natural Anodized
Colour: Metal
Dimensions: Height x Depth:
30 mm. x 10 mm.
Weight: 800 grams pr. meter
Standard Længde: 6 meters

Cutting:
Drilling Holes:

Item no.: JB 248-00-00
Material: Plastic (PE)
Colour: Grey
Dimensions: Height x Width:
30 mm. x 10 mm.
MOQ: 50 pcs.

CONTACT INFO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

ACCESSORIES FOR MEDICAL RAILS - SCANDINAVIAN STANDARDS



JB 100-01-07

End cap for JBM rail, right side.

Item no.: JB 100-01-07
Material: ABS Plastic
Colour: Grey
Dimensions: Height x Width
93 mm. x 81 mm.
Weight: 6 grams
MOQ: 50 pcs.



JB 100-01-10

Screw plugs for JBM rail.

Item no.: JB 100-01-10
Material: ABS Plastic
Colour: Grey
Dimensions: Diameter
Ø12 mm.
MOQ: 50 pcs.



JB 100-01-08

End cap for JBM rail, left side.

Item no.: JB 100-01-08
Material: ABS Plastic
Colour: Grey
Dimensions: Height x Width
93 mm. x 81 mm.
Weight: 6 grams
MOQ: 50 pcs.



CONTACT INFO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

WALL BRACKETS FOR BOTH RAIL STANDARDS



JB 286-00-00

Wall or equipment fitting for two types
of medical equipment rails.

Item no.:	JB 286-00-00
Material:	Stainless steel (AISI 304)
Colour:	Metal
Dimensions:	Diameter x Depth Ø50 mm. x 10 mm.
Weight:	137 grams
MOQ:	10 pcs.

JB 400-00-00

Wall or equipment fitting for two types
of medical equipment rails.

Item no.:	JB 400-00-00
Material:	Aluminium
Colour:	Metal
Dimensions:	Diameter x Depth Ø50 mm. x 10 mm.
Weight:	48 grams
MOQ:	10 pcs.

JB 286-00-18 (small)

JB 286-00-36 (large)

Spacer to two types of wall rails.

Item no.:	JB 286-00-18 (small)
Material:	Stainless steel (AISI 304)
Colour:	Metal
Dimensions:	Height / Outer diameter / Hole diameter 18 mm. / Ø15 mm. / Ø6,5 mm.
Weight:	20 grams
MOQ:	40 pcs.

Item no.:	JB 286-00-36 (large)
Material:	Stainless steel (AISI 304)
Colour:	Metal
Dimensions:	Height / Outer diameter / Hole diameter 36 mm. / Ø15 mm. / Ø6,5 mm.
Weight:	40 grams
MOQ:	20 pcs.

JB 400-00-18 (small)

JB 400-00-36 (large)

Spacer to two types of wall rails.

Item no.:	JB 400-00-18 (small)
Material:	Aluminium
Colour:	Metal
Dimensions:	Height / Outer diameter / Hole diameter 18 mm. / Ø15 mm. / Ø6,5 mm.
Weight:	7 grams
MOQ:	40 pcs.

Item no.:	JB 400-00-36 (large)
Material:	Aluminium
Colour:	Metal
Dimensions:	Height / Outer diameter / Hole diameter 36 mm. / Ø15 mm. / Ø6,5 mm.
Weight:	14 grams
MOQ:	20 pcs.

CONTACT INFO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

WALL BRACKET FOR PROTECTION OF MEDICAL EQUIPMENT RAILS ETC.



JB 180-02-02

Bracket that protects
medical equipment rails etc.

It absorbs shocks from beds that accidentally collide
with the rails.

Item no.:	JB 180-02-02
Material:	5 mm. Stainless steel (AISI 304)
Colour:	Metal
Dimensions:	Width x Height x Depth 400 mm. x 3.380 mm. x 1.450 mm.
Weight:	1.000 grams
MOQ:	4 pcs.



CONTACT INFO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

MEDICAL LIGHTING FROM GLAMOX LUXO

Examination lamps

JBMEDICO



Carelite LED

Carelite's LEDs are extremely energy efficient and have a long lifetime of at least 50,000 hours.

The powerful LED module has an excellent light output providing approximately 2,000 lux at 0.6 meters distance.

Item nr.

CAG026487
CAG026488
CAG026489
CAG026497
CAG026498
CAG026499



LHH LED

Modern medical examination lamp
LED - the light source of the future.

LHH's LEDs are extremely energy efficient and have a long lifetime of at least 50,000 hours. The powerful LED module has an excellent light output, providing more than 24,000 lux at 1 meter.

Item nr.

LHG026534
LHG026533
LHG026535
LHG026547
LHG026545
LHG026546

CONTACT INFO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

MEDICAL LIGHTING FROM GLAMOX LUXO

Magnifiers



Circus LED

Circus LED is suitable for demanding work in beauty salons, clinics and laboratories, or as a tool for people with impaired vision.

Circus LED has a large lens, aluminium head and a very flexible self-balancing arm and screen mount. The large lens allows comfortable viewing with both eyes.

Standard version: 3.5/5 diopter glass lens.

Item nr.

CIL026695
CIL026694
CIL026692
CIL026693



KFM LED

Advanced LED magnifier for many purposes
The KFM LED has a robust metal construction with a slim design, a flexible arm and excellent light.

It is suitable for many industrial applications, clinics and laboratories, and as an aid for people with impaired vision. Standard version comes with a 3 or 5 diopter glass lens. An included fabric bag protects the lens from dust and prevents the lamp from becoming a burning glass via sunlight. Colours: White or light grey.

Item nr.

KFL026034
KFL026035
KFL026036
KFL026037



CONTACT INFO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

MEDICAL LIGHTING FROM GLAMOX LUXO

Magnifiers and work lamps



KFM LED ESD

Developed for the electronics industry.

Static discharges can be devastating to electronic components. It is designed to dissipate static discharges in a controlled manner.

Item nr.

KFL026050
KFL026099



LFM LED G2

LFM LED is a lightweight magnifier with a metal construction, LED light sources, and a flexible spring balanced arm in the Luxo tradition.

Voltage from (V): 100, Lumen out: 700.

Item nr.

LFG028194
LFG028195



Wave LED

Wave LED provides a shadow-free magnification with light coming from both sides and which three-dimensional magnification with light from the left or right.

The 3D feature is particularly suitable for working with printed circuit boards and similar tasks.

Item nr.

WAL025948
WAL025949
WAL025950
WAL025951



CONTACT INFO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

ELECTRIC PILL CRUSHER V3 - SEVERO



JB 0830-2017

Crushes several standard pills at a time.



Item no.:	JB 0830-2017	Weight:	
Brand:	Severo		
Guarantee:	2 years on the device (1 year on battery - but does not include user errors)		
Cup volume:	30 ml.	Accessories:	Type:
Milling time:	8 seconds for standard pills	JB MED-020	Cups
Volume:	52 dB (low volume)	JB MED-002	Lids
Fully charged:	2 hours	JB MED-003	Rondels
		JB 3000001	Metal stamps
Contacts:	On/Off		
Indicator:	Battery level		
Dimensions:	Height x Width x Depth:		
	240 mm. x 70 mm. x 110 mm.		
	1.5 kg.		

CONTACT INFO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

ACCESSORIES FOR SEVERO ELECTRIC PILL CRUSHER V3



JB MED-020

Severo cups 2.100 pcs.

Disposable cup in hard plastic. The bottom has a serrated edge, so the cup can stand firmly in the Severo pellet crusher.

Item no.: JB MED-020
Material: PP
Dimensions: Height x Upper diameter
29 mm. x 47mm.
Volume: 30 ml.
MOQ: 2.100 pcs.



JB MED-002

Severo lid 2.000 pcs. for Severo cups.

Item no.: JB MED-002
Material: PE (soft polyethylene)
Colour: White
MOQ: 2.000 pcs.



JB MED-003

Rondel/Pill crusher 2.100 pcs. pr. box.

Keeps the pill crusher clean and ready to use immediately.

Item no.: JB MED-003
Material: PE
Diameter: 8 mm.
Thickness: 50 My
Certificates: CE
MOQ: 2 pcs. boxes with
2.100 pieces in each.



JB 3000001

Severo, extra metal stamp
for tablet crushing.

Item no.: JB 3000001
Material: Stainless steel (AISI 304)
Colour: Metal
Diameter: Ø55 mm.
MOQ: 2 pcs.

CONTACT INFO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

MONITOR ARMS AND ACCESSORIES

ADJUSTABLE ARM WITH GAS SPRING - JUNIOR ARMSTRONG



JB 21-00-00

Gas arm 2-8 kg. with handle for external locking of height position.

Mounting on wall, table, ceiling, equipment rails, equipment etc. All arms can be potential equalised. Supplied with all necessary accessories.

Item no.:	JB 21-00-00
Material:	Stainless steel (AISI 304)
Lifting height:	382 mm.
Range:	334 mm.
Weight:	2.300 grams



JB 21-00-15

Gas arm 8-18 kg. with handle for external locking of height position.

Mounting on wall, table, ceiling, equipment rails, equipment etc. All arms can be potential equalised. Supplied with all necessary accessories.

Item no.:	JB 21-00-15
Material:	Stainless steel (AISI 304)
Lifting height:	382 mm.
Range:	334 mm.
Weight:	2.300 grams



CONTACT INFO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

MONITOR ARMS AND ACCESSORIES

Junior Armstrong continued

WALL BRACKETS

JBMEDICO



JB 22-00-00

Manual extension arm with Ø20 mm. shaft, concealed cable routing.

The extension arm can be mounted on e.g. Wall bracket, table-clamp, ceiling-tube, rail clamps and IV-Poles.

Item no.: JB 22-00-00
Material: Stainless steel (AISI 304)
Total length: 255 mm.
Height: 150 mm.
Range: 220 mm.
Weight: 1.100 grams

Accessories:
Item no.: JB 94-00-00



JB 27-00-00

Monitor bracket with Ø20 mm. shaft.

Can be potential equalised.

Item no.: JB 27-00-00
Material: Stainless steel (AISI 304)
VESA: 75 mm. x 75 mm.
VESA: 100 mm. x 100 mm.
Dimensions: Depth x Height x Width
95 mm. x 159 mm. x 114 mm.
Weight: 1.000 grams

Accessories:
Item no.: JB 46-00-00



JB 31-00-00

Circular wall bracket fits Ø20 mm. shaft.

Unbrako bolt is standard.

Item no.: JB 31-00-00
Material: Stainless steel (AISI 304)
Dimensions: Depth x Height x Width
58 mm. x 123 mm. x 132 mm.
Weight: 688 grams

Accessories:
Item no.: JB 94-00-00



JB 241-00-00

Extended wall bracket fits Ø20 mm. shaft.

Unbrako bolt is standard.

Item no.: JB 241-00-00
Material: Stainless steel (AISI 304)
Dimensions: Depth x Height x Width
120 mm. x 123 mm. x 130 mm.
Weight: 1.100 grams

Accessories:
Item no.: JB 94-00-00

CONTACT INFO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

ACCESSORIES FOR MONITOR ARMS



JB 43-01-00

Keyboard Tray with hand support

The plate can be tilted -25°/ +115°

Item no.: JB 43-01-00
Material: Stainless steel (AISI 304)
Dimensions: Depth (without handle)
x Height x Width
277 mm. x 119 mm. x 579 mm.
Weight: 3.000 grams

Accessories:
Item no.: JB 46-00-00



JB 209-00-00

Suspension Tray for laptops with Ø20mm shaft.

Item no.: JB 209-00-00
Material: Stainless steel (AISI 304)
Dimensions: Depth (without handle)
x Height x Width
44 mm. x 62 mm. x 135 mm.
Weight: 3.800 grams

Accessories:
Item no.: JB 46-00-00



JB 253-00-00

Suspension Tray for medical equipment.

The plate can be tilted -25°/ +115°

Item no.: JB 253-00-00
Material: Stainless steel (AISI 304)
Dimensions: Depth (without handle)
x Height x Width
44 mm. x 62 mm. x 140 mm.
Weight: 2.200 grams

Accessories:
Item no.: JB 46-00-00



CONTACT INFO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361



JB 90-00-00

Standard clamp desk mount for Ø20 mm. shaft.

Mounting on the edge of the table top.
Table thickness 19-51 mm.
Unbrako bolt is standard.

Item no.: JB 90-00-00
Material: Stainless steel (AISI 304)
Dimensions: Depth x Height x Width
80 mm. x 186,5 mm. x 60 mm.
Weight: 2.500 grams
Surface: Sand blasted
Bushing: Ø20 mm.
Cone height: 94 mm.

Accessories: Locking handle can
be purchased
Item no.: JB 94-00-00



JB 91-00-00

Through hole table mount for Ø20 mm. shaft.

Mounting in the table top.
Table thickness 0-41 mm.
Unbrako bolt is standard.

Item no.: JB 91-00-00
Material: Stainless steel (AISI 304)
Dimensions: Depth x Height x Width
80 mm. x 135,5 mm. x 60 mm.
Weight: 2.000 grams
Surface: Sand blasted
Bushing: Ø20 mm.
Cone height: 95 mm.

Accessories: Locking handle can
be purchased
Item no.: JB 94-00-00



JB 46-00-00

Locking handle M10x40 mm.

For e.g. screen laptops and keyboard brackets for
tilt/rotate function.

Item no.: JB 46-00-00
Material: Stainless steel (AISI 304)
MOQ: 2 pcs.



JB 49-00-00

Bracket for mounting IT equipment on anaesthesia
equipment or as wall bracket for Ø20 mm. tube.

Different models for different brands of
equipment, e.g. Dräger and Dameca.

Item no.: JB 49-00-00
Material: Stainless steel (AISI 304)
Dimensions: Depth x Height x Width
122 mm. x 162 mm. x 73 mm.
Weight: 1.000 grams
Surface: Sand blasted
Materiale Bushing: Brass

Accessories: Locking handle can be
purchased
Item no.: JB 94-00-00



JB 94-00-00

Locking handle M6x20 mm.

For e.g. clamping and manual extension arm.

Item no.: JB 94-00-00
Material: Stainless steel (AISI 304)
MOQ: 5 pcs.

CONTACT INFO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361



JB 62-00-00

Pole clamp mount Ø20/Ø30/Ø20 mm.
for mounting IT equipment.

Item no.: JB 62-00-00
Material: Stainless steel (AISI 304)
Dimensions: Depth (without handle) x
Height x Width
44 mm. x 62 mm. x 235 mm.
Weight: 1.400 grams

Accessories:
Item no.: Locking handle available
JB 94-00-00



JB 63-00-00

Pole clamp mount Ø20/Ø20 mm.
for mounting IT equipment.

Item no.: JB 63-00-00
Material: Stainless steel (AISI 304)
Dimensions: Depth (without handle) x
Height x Width
44 mm. x 62 mm. x 135 mm.
Weight: 800 grams

Accessories:
Item no.: Locking handle available
JB 94-00-00



JB 64-00-00

Pole clamp mount Ø20/Ø20 mm.
for mounting IT equipment.

Item no.: JB 63-00-00
Material: Stainless steel (AISI 304)
Dimensions: Depth (without handle) x
Height x Width
44 mm. x 62 mm. x 135 mm.
Weight: 800 grams

Accessories:
Item no.: Locking handle available
JB 94-00-00



JB 65-00-00

Pole clamp mount Ø30/Ø20 mm.
for mounting IT equipment.

Item no.: JB 65-00-00
Material: Stainless steel (AISI 304)
Dimensions: Depth (without handle) x
Height x Width
44 mm. x 62 mm. x 235 mm.
Weight: 1.300 grams

Accessories:
Item no.: Locking handle available
JB 94-00-00

CONTACT INFO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

ROD SUSPENSION



JB 29-00-04

Solid rod
Ø20 x 600 mm.
2 end plugs.

Item no.: JB 29-00-04
Material: Stainless steel (AISI 304)
Dimensions: Diameter x Length
Ø20 mm. x 600 mm.
Surface: Glassblown
Weight: 1.700 grams



JB 29-00-10

Solid rod
Ø20 x 500 mm.
2 end plugs.

Item no.: JB 29-00-10
Material: Stainless steel (AISI 304)
Dimensions: Diameter x Length
Ø20 mm. x 500 mm.
Surface: Glassblown
Weight: 1.400 grams



CONTACT INFO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

SOAP AND ALCOHOL DISPENSERS For bottles



JB 06-21-30

Wire dispenser with 14 cm. arm
for 1.000 ml. bottles.

Item no.: JB 06-21-30
Material: Stainless steel (AISI 304)
Surface: Electropolished
Dimensions: Depth (incl. arm operation) x
Height x Width
149 mm. x 175 mm. x 90 mm.
Weight: 165 grams
MOQ: 12 pcs.



JB 06-21-40

Wire dispenser with 10 cm. arm
for 1.000 ml. bottles.

Item no.: JB 06-21-40
Material: Stainless steel (AISI 304)
Surface: Electropolished
Dimensions: Depth (incl. arm operation) x
Height x Width
145 mm. x 175 mm. x 90 mm.
Weight: 160 grams
MOQ: 12 pcs.



JB 06-21-50

Wire dispenser with 6 cm. arm
for 1.000 ml. bottles.

Item no.: JB 06-21-50
Material: Stainless steel (AISI 304)
Surface: Electropolished
Dimensions: Depth (incl. arm operation) x
Height x Width
141 mm. x 175 mm. x 90 mm.
Weight: 150 grams
MOQ: 12 pcs.

CONTACT INFO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

SOAP AND ALCOHOL DISPENSERS

For bags

JBMEDICO



JB 42-79-00

Wire dispenser with 14 cm. arm
for 1.000 ml. bags.



JB 42-90-03

Wire dispenser with 10 cm. arm
for 1.000 ml. bags.



JB 42-98-03

Wire dispenser with 6 cm. arm
for 1.000 ml. bags.



JB 728-87-67

Spacer, 85 mm. for all wire dispensers.
Used to mount the dispenser a distance from the wall,
e.g. over a table or sink.

Item no.: JB 42-79-00
Material: Stainless steel (AISI 304)
Surface: Electropolished
Dimensions: Depth (incl. arm operation) x
Height x Width
165 mm. x 240 mm. x 100 mm.
Weight: 165 grams
MOQ: 10 pcs.

Item no.: JB 42-90-03
Material: Stainless steel (AISI 304)
Surface: Electropolished
Dimensions: Depth (incl. arm operation) x
Height x Width
161 mm. x 240 mm. x 100 mm.
Weight: 145 grams
MOQ: 12 pcs.

Item no.: JB 42-98-03
Material: Stainless steel (AISI 304)
Surface: Electropolished
Dimensions: Depth (incl. arm operation) x
Height x Width
157 mm. x 240 mm. x 100 mm.
Weight: 145 grams
MOQ: 12 pcs.

Item no.: JB 728-87-67
Material: Stainless steel (AISI 304)
Surface: Electropolished
Dimensions: Depth x Height x Width
86 mm. x 58 mm. x 90 mm.
Weight: 132 grams
MOQ: 12 pcs.

CONTACT INFO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361



JB 102-00-01

Drip tray to suit everyone soap and alcohol dispensers.

For hanging on equipment or directly on the wall.
Protects floor and wall.

Item no.: JB 102-00-01
Material: Stainless steel (AISI 304)
Surface: Electropolished
Dimensions: Depth x Height x Width
84,5 mm. x 189,5 mm. x 91,5 mm.
Weight: 209 grams
MOQ: 12 pcs.



JB 71-72-51

Drip tray, drip catcher Ø60 mm. for mounting on wire soap and alcohol dispensers.

Protects the floor.

Item no.: JB 71-72-51
Material: Stainless steel (AISI 304)
Surface: Electropolished
Dimensions: Depth x Height x Width
80 mm. x 180 mm. x 70 mm.
Weight: 67 grams
MOQ: 12 pcs.



JB 213-00-00

Adapter T-slot fitting for soap and alcohol detergent Wet Wipe dispenser and drip trays.

Enables removal and cleaning of dispensers without the use of too.

Item no.: JB 213-00-00
Material: Stainless steel (AISI 304)
Surface: Electropolished
Dimensions: Depth x Height x Width
14 mm. x 50 mm. x 92 mm.
Weight: 174 grams
MOQ: 10 pcs.

CONTACT INFO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

DISPENSERS WITH DRIP TRAY AND ADAPTER BRACKET

For bottles



JB 50-213-102

Steel dispenser with 6 cm. long arm for 1.000 ml. bottles with drip tray and adapter fitting.

Item no.: JB 50-213-102
Material: Stainless steel (AISI 304)
Colour: Metal
Surface: Polished
Weight: 537 grams
MOQ: 4 pcs.

Accessories:
Item no.: JB 47-00-00

Item no.: JB 120-00-00



JB 40-213-102

Steel dispenser with 10 cm. long arm for 1.000 ml. bottles with drip tray and adapter fitting.

Item no.: JB 40-213-102
Material: Stainless steel (AISI 304)
Colour: Metal
Surface: Polished
Weight: 547 grams
MOQ: 4 pcs.

Accessories:
Item no.: JB 47-00-00

Item no.: JB 120-00-00



JB 30-213-102

Steel dispenser with 14 cm. long arm for 1.000 ml. bottles with drip tray and adapter fitting.

Item no.: JB 30-213-102
Material: Stainless steel (AISI 304)
Colour: Metal
Surface: Polished
Weight: 551 grams
MOQ: 4 pcs.

Accessories:
Item no.: JB 47-00-00

Item no.: JB 120-00-00



CONTACT INFO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

DISPENSERS WITH DRIP TRAY AND ADAPTER BRACKET

For bags



JB 98-213-102

Steel dispenser with 6 cm. long arm for 1.000 ml. bags with drip tray and adapter fitting.

Item no.: JB 98-213-102
Material: Stainless steel (AISI 304)
Colour: Metal
Surface: Polished
Weight: 537 grams
MOQ: 4 pcs.

Accessories:
Item no.: JB 47-00-00

Item no.: JB 120-00-00



JB 90-213-102

Steel dispenser with 10 cm. long arm for 1.000 ml. bags with drip tray and adapter fitting.

Item no.: JB 90-213-102
Material: Stainless steel (AISI 304)
Colour: Metal
Surface: Polished
Weight: 532 grams
MOQ: 4 pcs.

Accessories:
Item no.: JB 47-00-00

Item no.: JB 120-00-00



JB 79-213-102

Steel dispenser with 14 cm. long arm for 1.000 ml. bags with drip tray and adapter fitting.

Item no.: JB 79-213-102
Material: Stainless steel (AISI 304)
Colour: Metal
Surface: Polished
Weight: 553 grams
MOQ: 4 pcs.

Accessories:
Item no.: JB 47-00-00

Item no.: JB 120-00-00



CONTACT INFO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

STAINLESS STEEL COVER FOR WIRE DISPENSERS



JB 25-06-20

Cover and lock for wire dispensers in Electropolished stainless steel with a back plate

Fire retardant for installation in e.g. escape routes.

It consists of 2 parts:
1. Back plate that screws into the wall, it has the same hole patterns as the wire dispensers
2. Cover that slides over the back plate and can be secured with a lock on both sides (padlock or strips)

It fits all thread dispensers with standbag bottles. It also fits all thread dispensers with bags from e.g. Plum

Item no.:	JB 25-06-20
Material:	Stainless steel (AISI 304)
Surface / Colour:	Electropolished or white powder coated
Dimensions:	Depth x Height x Width
Weight:	110 mm. x 338 mm. x 130 mm.
MOQ:	1.297 grams
	2 pcs.
Accessories:	
Dispenser 14 cm. arm for bags	JB 42-79-00
Dispenser 10 cm. arm for bags	JB 42-90-03
Dispenser 6 cm. arm for bags	JB 42-98-03
Dispenser 14 cm. arm for bottles	JB 06-21-30
Dispenser 10 cm. arm for bottles	JB 06-21-40
Dispenser 6 cm. arm for bottles	JB 06-21-50
Drip tray with back plate	JB 102-00-01
Drip tray, drip catcher Ø60 mm.	JB 71-72-51
Adapter T-slot fitting	JB 213-00-00



CONTACT INFO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

TABLE AND WALL DISPENSER FOR SQUARE AND ROUND BOTTLES



JB 162-67-87

Table and wall dispenser for round and square 500- 600 ml. bottles

With a handle for T-slot, stainless steel

It consists of 2 pieces:

1. Wire dispensers in 5 mm. wire
2. Table clamp

Fits most square alcohol and round soap bottles. (500 and 600 ml.) Square bottles 67x67 mm.
Round bottles Ø70 mm.

Item no.:

Material:

Colour:

Surface:

Dimensions:

Weight:

MOQ:

Accessories:

Wall fitting with T-slot

Rail clamp with T-slot (DK)

Rail clamp with T-slot (EU)

Multibracket with T-slot

GCX rail fitting with T-slot

JB 162-67-87

Stainless steel (AISI 304)

Metal

Electropolished

Depth x Height x Width
87 mm. x 162,5 mm. x 77 mm.

100 grams

10 pcs. pr. pcs.

JB 47-00-00

JB 120-00-00

JB 146-00-00

JB 158-00-163

JB 820-625-163



Squared bottle 67x67 mm.

Round bottle Ø70 mm.

CONTACT INFO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

GLOVE, WET WIPE AND SANITARY FACEMASK HOLDER

New HandsUp

JBMEDICO



JB 800-00-00

With "New HandsUp", staff can finally see the glove sizes, which are always written on the sides of the boxes.

Fits 99 percent of people glove boxes.

Fits Wet Wipes (the large MAXI)

Fits sanitary towels and mask boxes

"New HandsUp" can be mounted in both horizontal and vertical position.

Mounts to wall fitting, T-slot or wall rail

Autoclavable

Item no.:
Material:

Surface:
Dimensions:

Weight:
MOQ:

JB 800-00-00
Stainless steel (AISI 304) -
4mm. spring steel
Electropolished
Depth x Height x Width
105 mm. x 148 mm. x 260 mm.
362 grams
4 pcs.

Accessories:

Wall bracket with T-slot
Rail clamp with T-slot (DK)
Rail clamp with T-slot (EU)
Multibracket with T-slot
GCX rail bracket with T-slot

JB 47-00-00
JB 120-00-00
JB 146-00-00
JB 158-00-163
JB 820-625-163



CONTACT INFO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361



JB 16-05-41-02

Steel dispenser, Wet Wipe MAXI

Stainless steel wall holder for wall mounting of Wet Wipe products.

Item no.:	JB 16-05-41-02
Material:	Stainless steel (AISI 316)
Brand:	Wet Wipe
Colour:	Metal
Surface:	Electropolished
Weight:	230 grams
MOQ:	25 pcs.

Accessories:
JB 213-00-00



JB 16-05-40-02

Steel dispenser, Wet Wipe MINI

Stainless steel wall holder for wall mounting of Wet Wipe products.

Item no.:	JB 16-05-40-02
Material:	Stainless steel (AISI 316)
Brand:	Wet Wipe
Colour:	Metal
Surface:	Electropolished
Weight:	180 grams
MOQ:	25 pcs.

Accessories:
JB 213-00-00

CONTACT INFO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

ACCESSORIES FOR GLOVE BOX HOLDERS

JBMEDICO



JB 47-00-00

T-slot wall bracket for 1 pcs. glove box holder.

C-profile with 2 pcs. countersunk mounting holes and stop screw.



JB 47-00-02

T-slot wall bracket for 2 pcs. glove box holder.

C-profile with 4 pcs. countersunk mounting holes and stop screw.



JB 47-00-03

T-slot wall bracket for 3 pcs. glove box holder.

C-profile with 6 pcs. countersunk mounting holes and stop screw.



JB 47-00-04

T-slot wall bracket for 4 pcs. glove box holder.

C-profile with 6 pcs. countersunk mounting holes and stop screw.

Item no.: JB 47-00-00
Material: Aluminium
Surface: Anodized
Dimensions: Depth x Height x Width
12,5 mm. x 60 mm. x 36 mm.
Hole Size: Ø 6,4 mm.
Weight: 36 grams
MOQ: 10 pcs.

Item no.: JB 47-00-02
Material: Aluminium
Surface: Anodized
Dimensions: Depth x Height x Width
12,5 mm. x 210 mm. x 36 mm.
Hole Size: Ø 6,4 mm.
Weight: 129 grams
MOQ: 5 pcs.

Item no.: JB 47-00-03
Material: Aluminium
Surface: Anodized
Dimensions: Depth x Height x Width
12,5 mm. x 360 mm. x 36 mm.
Hole Size: Ø 6,4 mm.
Weight: 222 grams
MOQ: 5 pcs.

Item no.: JB 47-00-04
Material: Aluminium
Surface: Anodized
Dimensions: Depth x Height x Width
12,5 mm. x 520 mm. x 36 mm.
Hole Size: Ø 6,4 mm.
Weight: 320 grams
MOQ: 5 pcs.

CONTACT INFO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

ALUMINIUM SLIDING BRACKET FOR T-SLOT RAILS



JB 90-92-05-30

Aluminium sliding bracket,
lockable

Hole patterns in plate part
fit various accessories and
dispensers

Item no.:
Material:
Colour:
Dimensions:


Weight:
MOQ:

JB 90-92-05-30
Aluminium
Metal
Depth x Height x Width
87 mm. x 162,5 mm. x 77 mm.
150 grams
4 pcs.

Accessories:

Steel dispenser, Wet Wipe MAXI
Steel dispenser, Wet Wipe MINI
T-slot bracket length 210 mm.
T-slot bracket length 360 mm.
T-slot bracket length 520 mm.
CombiPlum wire dispenser 14 cm.
CombiPlum wire dispenser 10 cm.
CombiPlum thread dispenser 6 cm.
Dispenser for 1.000 ml. bottles 14 cm.
Dispenser for 1.000 ml. bottles 10 cm.
Dispenser for 1.000 ml. bottles 6 cm.
Drip tray for wire dispensers
Fixing device Ø25 bushing hole
Fixing device Ø18 socket hole

JB 16-05-41-02
JB 16-05-40-02
JB 47-00-02
JB 47-00-03
JB 47-00-04
JB 42-79-00
JB 42-90-03
JB 42-98-03
JB 06-21-50
JB 06-21-40
JB 06-21-30
JB 102-00-01
JB 205-00-00
JB 205-00-18



Accessories proposal No. 1



Accessories proposal No. 2



Accessories proposal No. 3



Accessories proposal No. 4

CONTACT INFO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

GLOVE BOX HOLDER



JB 110-00-00

Glove holder for 3 glove boxes wall mounted.

Suitable for glove boxes, face wipes, aprons and handkerchiefs, etc.

Item no.:	JB 110-00-00
Material:	Stainless steel (AISI 304)
Dimensions:	Depth x Height x Width (outside) 93 mm. x 406 mm. x 285 mm.
Opening:	183 mm.
Weight:	1.854 grams
MOQ:	4 pcs.



JB 110-00-04

Glove holder for 4 glove boxes wall mounted.

Suitable for glove boxes, face wipes, aprons and handkerchiefs, etc.

Item no.:	JB 110-00-04
Material:	Stainless steel (AISI 304)
Dimensions:	Depth x Height x Width (outside) 93 mm. x 520 mm. x 285 mm.
Opening:	183 mm.
Weight:	2.346 grams
MOQ:	4 pcs.

CONTACT INFO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

IV STANDS / IV POLES



JB 306-2-317-211

"One Hand"/"Small" straight legs, large wheels,
long straight rod, 2 strong angled hooks

The flexible IV stand can be modified as needed. The solid and rust-proof specially designed base allows the attachment of various items in a safe way without the risk of the stand tipping over.

Item no.:

Material / Surface
Color of Plastic Parts
Hooks
Load capacity pr. hook
Weight without equipment
Max. load capacity
Height with/without wheels (telescopic max.)
Height with/without wheels (telescopic min.)
Outer tube diameter and length
Inner tube diameter and length
Setting
Height adjustment
Base diameter, straight leg
Wheel diameter according to DIN EN 60601

JB 306-2-317-211

Stainless steel (AISI 304) / Electropolished
Black
2 strong angled hooks
15 kg.
15 kg.
40 kg.
2.369 mm. / 2.232 mm.
1.605 mm. / 1.470 mm.
Ø32/1.281 mm.
Ø20 x 940 mm.
"One Hand" adjustable
765 mm.
515 mm.
Ø100 mm.
1 x antistatic wheels
2 x lockable wheels
2 x standard wheels

Safety distance between both tubes when outer tube reaches its lowest height.

Integrated stop prevents separation of the adjustable tubes.

Stainless steel ball bearing: Yes
An integrated base-weight facilitates a safe of various pumps and other medical equipment.
CE Approved



Security

- Columns are equipped with plastic stops (no risk of damage)
- Sufficient safety distance in the collapsed state as a precaution against risk of injury (e.g. pinching of fingers)
- Increased stability due to conical connection between stand and foot section
- Safe control of the lowering of the inner tube
- No risk of injury due to rounded design of hooks
- Due to the geometry/angle of the hooks, the infusion bottles cannot fall out
- Flexible handle for drip stand can be adapted to the individual user.

Technology

- All plastic parts are made of lightweight materials
- Pipes in Stainless steel according to DIN EN 10217-7
- Material ultrasonic and eddy current tested (according to EN 10246-9 /3).
- Integrated pull-out limiter prevents unintentional separation of the adjustable tubes
- Safety bottle cross according to DIN ISO 15375
- Good running characteristics through high quality 100 mm. mono Rolls, 1 Rolls are antistatic apply (One Hand stands).
- The introduction of hygienically sensitive liquids is prevented by the sealing lid.
- Reinforced bottle cross allows connection of 5 liters bottles.

CONTACT INFO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

IV POLE XL ONE HAND - FOUR HOOKS

JB MEDICO



JB 332-2-320-190

Drop stand "One Hand" with curved legs, large wheels,
curved long rod, 4 strong angled hooks.

The flexible IV stand can be modified as needed. The solid and rust-proof specially designed base allows the attachment of various items in a safe way without the risk of the stand tipping over.

Item no.:

Material / Surface
Color of Plastic Parts
Hooks
Load capacity pr. hook
Weight without equipment
Max. load capacity
Height with/without wheels (telescopic max.)
Height with/without wheels (telescopic min.)
Outer tube diameter
Inner tube diameter
Setting
Height adjustment
Base diameter, straight leg
Wheel diameter according to DIN EN 60601

Safety distance between both tubes when outer tube reaches its lowest height

Integrated stop prevents separation of the adjustable tubes.

JB 332-2-320-190

Electropolished Stainless steel (AISI 304)
Black
Four strong hooks
15 kg.
15 kg.
50 kg.
2.281 mm. / 2.144 mm.
1.516 mm. / 1.379 mm.
Ø32 mm.
Ø20 mm.
"One Hand" adjustable
765 mm.
545 mm.
Ø100 mm.
1 x antistatic wheels
2 x lockable wheels
2 x standard wheels

Stainless steel ball bearing: Yes
An integrated base-weight facilitates a secure attachment of various pumps and other medical equipment.
CE Approved



Security

- Pillars are equipped with plastic stops (no risk of injury)
- Sufficient safety clearance in collapsed condition as a precaution against risk of injury (e.g. pinching of fingers)
- Increased stability due to conical connection between stand and foot section
- Safe control of the lowering of the inner tube
- No risk of injury due to rounded design of hooks
- Due to the geometry/angle of the hooks, the infusion bottles cannot fall out
- Flexible handle for drip stand can be adapted to the individual user.

Technology

- All plastic parts are made of lightweight materials
- Pipes in Stainless steel according to DIN EN 10217-7
- Material ultrasonic and eddy current tested (according to EN 10246-9 /3).
- Integrated pull-out limiter prevents unintentional separation of the adjustable tubes
- Safety bottle cross according to DIN ISO 15375
- Good running characteristics through high quality 100 mm. mono Rolls, 1 Rolls are antistatic apply (One Hand stands).
- The introduction of hygienically sensitive liquids is prevented by the sealing lid.
- Reinforced bottle cross allows connection of 5 liters bottles.

CONTACT INFO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

IV POLE XXL ONE HAND - SIX HOOKS

JB MEDICO



JB 332-2-321-216

Drop stand "One Hand" with curved legs, large wheels,
curved long rod, 6 strong straight hooks

The flexible IV stand can be modified as needed. The solid and stainless steel frame, allows the attachment of various items in a safe way without the risk of the stand tipping over.

Item no.:

Material / Surface
Color of Plastic Parts
Hooks
Load pr. hook
Weight without equipment
Max. load
Height with/without wheels (telescopic max.)
Height with/without wheels (telescopic min.)
Outer tube diameter
Inner tube diameter
Setting
Height adjustment
Base diameter, straight leg
Wheel diameter according to DIN EN 60601

Safety distance between both tubes when outer tube reaches its lowest height.

Integrated stop prevents separation of the adjustable tubes.

Security

- Columns are equipped with plastic stops (no risk of damage)
- Sufficient safety distance in the collapsed state as a precaution against risk of injury (e.g. pinching of fingers)
- Increased stability due to conical connection between stand and foot section
- Safe control of the lowering of the inner tube
- No risk of injury due to rounded design of hooks
- Due to the geometry/angle of the hooks, the infusion bottles cannot fall out
- Flexible handle for drip stand can be adapted to the individual user.

JB 332-2-321-216

Electropolished Stainless steel (AISI 304)
Black
Six strong straight hooks
15 kg.
15 kg.
50 kg.
2.281 mm. / 2.144 mm.
1.516 mm. / 1.379 mm.
Ø32 mm.
Ø20 mm.
"One Hand" adjustable
765 mm.
545 mm.
Ø100 mm.
1 x antistatic wheels
2 x lockable wheels
2 x standard wheels

Stainless steel ball bearing: Yes
An integrated base-weight facilitates a secure attachment of various pumps and other medical equipment
CE Approved



Technology

- All plastic parts are made of lightweight materials
- Pipes in Stainless steel according to DIN EN 10217-7
- Material ultrasonic and eddy current tested (according to EN 10246-9 /3).
- Integrated pull-out limiter prevents unintentional separation of the adjustable tubes
- Safety bottle cross according to DIN ISO 15375
- Good running characteristics through high quality 100 mm. mono Rolls, 1 Rolls are antistatic apply (One Hand stands).
- The introduction of hygienically sensitive liquids is prevented by the sealing lid.
- Reinforced bottle cross allows connection of 5 liters bottles.

CONTACT INFO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

IV STAND / IV POLE



JB 306-2-317-190

Drop stand "One Hand" with straight legs, large wheels, straight long rod, 4 strong angled hooks.

The flexible IV stand can be modified as needed. The solid and stainless steel frame, allows the attachment of various items in a safe way without the risk of the stand tipping over.

Item no.:	JB 306-2-317-190
Material / Surface	Electropolished Stainless steel (AISI 304)
Color of Plastic Parts	Black
Hooks	Four strong straight hooks
Load pr. hook	15 kg.
Weight without equipment	15 kg.
Max. load	40 kg.
Height with/without wheels (telescopic max.)	2.369 mm. / 2.232 mm.
Height with/without wheels (telescopic min.)	1.605 mm. / 1.470 mm.
Outer tube diameter	Ø32/1.281 mm.
Inner tube diameter	Ø20 x 940 mm.
Setting	"One Hand" adjustable
Height adjustment	765 mm.
Base diameter, straight leg	515 mm.
Wheel diameter according to DIN EN 60601	Ø100 mm.
	1 x antistatic wheels
	2 x lockable wheels
	2 x standard wheels
Safety distance between both tubes when outer tube reaches its lowest height.	Stainless steel ball bearing: Yes An integrated base-weight facilitates a secure attachment of various pumps and other medical equipment
Integrated stop prevents separation of the adjustable tubes.	CE Approved



Security

- Columns are equipped with plastic stops (no risk of damage)
- Sufficient safety distance in the collapsed state as a precaution against risk of injury (e.g. pinching of fingers)
- Increased stability due to conical connection between stand and foot section
- Safe control of the lowering of the inner tube
- No risk of injury due to rounded design of hooks
- Due to the geometry/angle of the hooks, the infusion bottles cannot fall out
- Flexible handle for drip stand can be adapted to the individual user.


Technology

- All plastic parts are made of lightweight materials
- Pipes in Stainless steel according to DIN EN 10217-7
- Material ultrasonic and eddy current tested (according to EN 10246-9 / 3).
- Integrated pull-out limiter prevents unintentional separation of the adjustable tubes
- Safety bottle cross according to DIN ISO 15375
- Good running characteristics through high quality 100 mm. mono Rolls, 1 Rolls are antistatic apply (One Hand stands).
- The introduction of hygienically sensitive liquids is prevented by the sealing lid.
- Reinforced bottle cross allows connection of 5 liters bottles.

CONTACT INFO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361


IV STAND / INFUSION POLE - SUPPORT COLUMNS



JB 320-00-02

Rod long, curved.


Item no.: JB 320-00-02
Dimensions: Length x Depth:
1.281 mm. x 32 mm.
Weight: 2.000 grams



JB 317-00-06

Rod long, straight.

Item no.: JB 317-00-06
Dimensions: Length x Width:
1.281 mm. x Ø32 mm.
Weight: 1.800 grams



JB 317-00-13

Telescopic pole
"One Hand".

Item no.: JB 317-00-13
Dimensions: Length x Depth:
940 mm. x Ø20 mm.
Weight: 1.000 grams



CONTACT INFO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

INFUSION AND IV POLE - ACCESSORIES

Hooks



JB 204-00-00

Hook top - straight 2 strong hooks including bolt.



JB 211-00-00

Hook top - angled 2 strong hooks including bolt.



JB 212-00-00

Hook top - angled 3 strong hooks including bolt.



JB 190-00-08

Hook top - 4 strong hooks including bolt.

Item no.: JB 204-00-00
Material: Stainless steel (AISI 304)
Surface: Electropolished
Dimensions: Height x Width
406 mm. x 285 mm.
Hook distance: 252 mm.
Weight: 248 grams

Item no.: JB 211-00-00
Material: Stainless steel (AISI 304)
Surface: Electropolished
Dimensions: Height x Width
406 mm. x 285 mm.
Hook distance: 204 mm.
Weight: 248 grams

Item no.: JB 212-00-00
Material: Stainless steel (AISI 304)
Surface: Electropolished
Dimensions: Height x Width
406 mm. x 285 mm.
Hook distance: 120 mm.
Weight: 363 grams

Item no.: JB 190-00-08
Material: Stainless steel (AISI 304)
Surface: Electropolished
Dimensions: Height x Width
406 mm. x 285 mm.
Hook distance: 126 mm.
Weight: 448 grams

CONTACT INFO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

IV STAND / IV POLE - ACCESSORIES

Hooks



JB 216-00-00

Hook top - 6 strong hooks including bolt.




JB 170-02-10

IV 2-Hooks for infusion bags, AISI 304



JB 280-03-10

IV 3-Hooks for infusion bags, AISI 304



JB 390-04-10

IV 4-Hooks for infusion bags, AISI 304

Item no.: JB 216-00-00
Material: Stainless steel (AISI 304)
Surface: Electropolished
Dimensions: Height x Width
406 mm. x 285 mm.
Hook distance: 120 mm.
Weight: 596 grams

Item no.: JB 280-02-10
Material: 10 mm. stainless steel (AISI 304)
Surface: Electrically polished
Dimensioner: Height x Width
44,88 mm. x 170 mm.

Total load capacity: 4 kg.
Distance between hooks: 110 mm.
Weight: 299 grams
MOQ: 2 pcs.

Item no.: JB 280-03-10
Material: 10 mm. stainless steel (AISI 304)
Surface: Electrically polished
Dimensioner: Height x Width
44,88 mm. x 280 mm.

Total load capacity: 6 kg.
Distance between hooks: 110 mm.
Weight: 331 grams
MOQ: 2 pcs.

Item no.: JB 390-04-10
Material: 10 mm. stainless steel (AISI 304)
Surface: Electrically polished
Dimensioner: Height x Width
44,88 mm. x 390 mm.

Total load capacity: 8 kg.
Distance between hooks: 110 mm.
Weight: 433 grams
MOQ: 2 pcs.

CONTACT INFO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

IV STAND / IV POLE - ACCESSORIES

Wheels and bases



5387PJC100P30-13

Swivel castor, 100 mm. total brake.
Stainless (Levina), Grey.

Item no.: 5387PJC100P30-13
Wheel Ø: Ø13 mm.
Centre hole Ø: Ø13 mm.
Overall height: 137 mm.
Load: 100 kg. (dynamic)
Standard: EN 12530
Weight: 430 grams
Stainless steel ball bearing: Yes
Electrical conductivity: No
MOQ: 18 pcs.



5380PJC100P30-13

Swivel castor without brake 100 mm.
Stainless (Levina), Grey.

Item no.: 5380PJC100P30-13
Wheel Ø: Ø13 mm.
Centre hole Ø: Ø13 mm.
Overall height: 137 mm.
Load: 100 kg (dynamisk)
Standard: EN 12530
Weight: 430 grams
Stainless steel ball bearing: Yes
Electrical conductivity: No
MOQ: 20 pcs.



5380XSC100P30-13

Swivel castor without brake 100 mm.
Stainless (Levina) Grey, Conductive.

Item no.: 5380XSC100P30-13
Wheel Ø: Ø13 mm.
Centre hole Ø: Ø13 mm.
Overall height: 137 mm.
Load: 100 kg (dynamisk)
Standard: EN 12530
Weight: 430 grams
Stainless steel ball bearing: Yes
Electrical conductivity: Yes
MOQ: 20 pcs.



JB 332-00-00

Main base, bent legs for drop stand.
(Tente wheels are included)

Item no.: JB 332-00-00
Base diameter: Ø545 mm.
Leg Profile: Ø25 mm.
Weight: 12,3 kg.



JB 306-00-00

Main base, straight leg for drop stand.
(Tente wheels are included)

Item no.: JB 306-00-00
Base diameter: Ø545 mm.
Leg Profile: Ø25 mm.
Weight: 11,6 kg.

CONTACT INFO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361



JB 294-00-05

Stainless steel handle with an ergonomic grip.

The handle is easy to mount and adjustable, can easily be set in relation to other equipment/user's height.

Item no.: JB 294-00-05
 Material: Stainless steel (AISI 304)
 Surface: Electropolished
 Dimensions: Height x Width x Length
 36 mm. x 194 mm. x 280 mm.
 Handle diameter: Ø20 mm.
 Weight: 733 grams

Accessories: JB 163-00-01
 Dimensions: H: 60 mm. x W: 36,2 mm.
 x D: 12,5 mm.
 Weight: 36 grams



JB 190-00-05

Telescopic rod 940 mm. on rail or tube-mounted.

Several combinations (Picture is with hook top - to be ordered separately).

Item no.: JB 190-00-05
 Material: Stainless steel (AISI 304)
 Surface: Electropolished
 Length: 940 mm.
 Diameter: Ø20 mm.
 Hook distance: 204 mm.
 Weight: 795 grams



JB 289-00-00

Transducer holder for 4 pcs. A needles.

A-cannulas made of aluminium.
 Can be washed, fits most transducers.

Item no.: JB 289-00-00
 Material: Stainless steel (AISI 304)
 Surface: Electropolished
 Dimensions: Height x Width x Depth
 90 mm. x 185 mm. x 19,5 mm.
 Weight: 241 grams
 MOQ: 4 pcs.

Accessories: JB 158-00-163
 (Multibracket med T-spors
 beslag)



JB 205-08-05

Drip pole for infusion bags for IV fluids.

Fits all our products with T-slot handle.

Item no.: JB 205-08-05
 Material: Stainless steel (AISI 304)
 Surface: Electropolished
 Dimensions: Height x Width
 50 mm. x 205 mm.
 Capacity: 5-6 drip bags
 Weight: 135 grams
 MOQ: 5 pcs.

CONTACT INFO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361



JB 25-38-BB

Lock for B. Braun Pump Rack

Securing against rotation of the "Docking Station" on the stand so that the weight of the pumps does not cause the Infusion Stand to be unstable

The lock consists of a front plate that slides into the "rack" and a rear plate that is clamped around the infusion stand with four bolts and four nuts.

+ an extra bolt to prevent unwanted removal without the use of tools

Item no.:

Material:

Surface:

Dimensions:

Weight:

MOQ:

Note:

JB 25-38-BB

Stainless steel (AISI 304)

Glassblown

Height x Width

130 mm. x 75 mm.

387 grams

2 pcs.

Suitable for pipe diameter from 25-38 mm.



CONTACT INFO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

WALL BRACKET FOR CO2 BOTTLES, Ø183 MM. HOLE



JB 283-01-00

Wall bracket for CO2 bottles fits for example Air Liquide 6 liters.

Bottle holder for wall for up to Ø183 mm. (Inside dimensions)

Item no.:	JB 283-01-00
Material:	Stainless steel (AISI 304) - 4 mm.
Surface:	Electropolished
Dimensions:	Height x Width 306 mm. x 215 mm.
Weight:	2.340 grams
Note:	Fits up to Ø183 mm. bottles



CONTACT INFO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361



Picture is with our multi fitting

JB 290-01-00

Oxygen and gas cylinder holder top
for Ø117 mm. bottles.

Can be ordered with a bottom ring JB 290-02-00

Item no.: JB 290-01-00
Material: Ø8 mm. SS, POM,
Nylon(PA6)
Width: 175 mm. Ext. w. handle
Weight: 255 grams
MOQ: 2 pcs.

Accessories: Multibracket
Item no.: JB 158-00-163



Picture is with our multi fitting

JB 290-02-00

Oxygen and gas cylinder holder-bottom
for Ø117 mm. cylinders.

Can be ordered with a bottom ring JB 290-02-00

Item no.: JB 290-02-00
Material: Ø8 mm. SS, POM,
Nylon(PA6)
Width: 189 mm. Ext. w. handle
Height: 128 mm.
Weight: 315 grams
MOQ: 2 pcs.

Accessories: Multibracket
Item no.: JB 158-00-163



Picture is with our multi fitting

JB 277-01-00

Oxygen and gas cylinder holder top
for Ø140 mm. cylinders.

Can be ordered with a bottom ring JB 290-02-00

Item no.: JB 277-01-00
Material: Ø8 mm. SS, POM,
Nylon(PA6)
Width: 189 mm. Ext. w. handle
Weight: 278 grams
MOQ: 2 pcs.

Accessories: Multibracket
Item no.: JB 158-00-163



Picture is with our multi fitting

JB 277-02-00

Oxygen and gas cylinder holder-bottom
for Ø140 mm. cylinders.

Can be ordered with a bottom ring JB 290-02-00

Item no.: JB 277-02-00
Material: Ø8 mm. SS, POM,
Nylon(PA6)
Width: 192 mm. Ext. w. handle
Height: 122 mm.
Weight: 341 grams
MOQ: 2 pcs.

Accessories: Multibracket
Item no.: JB 158-00-163

CONTACT INFO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

SHARPS CONTAINER HOLDER - ROUND FOR E.G. "Grathwol"

JBMEDICO



Sharps Container 0,7 and 1,5 liters

Sharps container 0,7 and 1,5 liters
AP Medical DISPO Line.

Yellow or red lid.

Accessories:

Wall bracket with T-slot
Rail clamp with T-slot (DK)
Rail clamp with T-slot (EU)
Multibracket with T-slot
GCX rail bracket with T-slot



JB 255-00-00

Sharps container holder til e.g.
AP Medical 0,7 and 1,5 liters Sharps Containers.

Item no.:	JB 255-00-00
Inner diameter:	Ca. Ø100 mm.
Material:	Stainless steel (AISI 304)
Surface:	Electropolished
Height:	50 mm.
Depth:	143 mm. (outer dimension)
Thread diameter:	8 mm.
Weight:	189 grams
MOQ:	5 pcs.

CONTACT INFO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

SHARPS CONTAINER HOLDER - ROUND FOR E.G. "AP MEDICAL"



Sharps Container 2 liters

Sharps container 2 liters
AP Medical DISPO Line.

Yellow or red lid.



JB 149-00-00

Sharps container holder for
1.0 - 1,3 - 1,5 - 2 and 3 liters Sharps Containers.

Fits perfectly to Klinion EasyCare "Mediq /UDSON" 1,3
and 2 liters AP Medical DISPO Line boxes.

Yellow or red lid.

Item no.:	JB 149-00-00
Inner diameter:	Ø 134 mm.
Material:	Stainless steel (AISI 304)
Surface:	Electropolished
Height:	50 mm.
Depth:	172 mm. (outer dimension)
Thread diameter:	8 mm.
Weight:	220 grams
MOQ:	5 pcs.

SHARPS CONTAINER HOLDER - OVAL FOR E.G. "AP MEDICAL"

JBMEDICO



Sharps Container 2, 4 and 7 liters

Sharps container 2 liters
AP Medical DISPO Line.

Yellow or red lid.

Accessories:
Wall bracket with T-slot
Rail clamp with T-slot (DK)
Rail clamp with T-slot (EU)
Multibracket with T-slot
GCX rail bracket with T-slot

JB 159-00-00

Sharps container holder for e.g. AP Medical 2, 4 and 7
liters oval boxes.

Fits perfectly to Klinion EasyCare "Mediq /UD-
SON" AP Medical DISPO Line boxes.
Yellow or red lid.

Item no.:	JB 159-00-00
Inner diameter:	133 mm. x 133 mm.
Material:	Stainless steel (AISI 304)
Surface:	Electropolished
Height:	50 mm.
Depth:	175 mm. (outer dimension)
Thread diameter:	8 mm.
Weight:	248 grams
MOQ:	5 pcs.

CONTACT INFO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

SHARPS CONTAINER HOLDER - OVAL

For E.G. "AP Medical"



Sharps Container 3 liters

Sharps container 3 liters AP Medical DISPO Line.

Yellow or red lid.



JB 169-00-00

Sharps container holder for e.g. AP Medical
3 liters, oval container.

Accessories:
Wall bracket with T-slot
Rail clamp with T-slot (DK)
Rail clamp with T-slot (EU)
Multibracket with T-slot
GCX rail bracket with T-slot

Item no.:	JB 169-00-00
Inner diameter:	143 mm. x 143 mm.
Material:	Stainless steel (AISI 304)
Surface:	Electropolished
Height:	50 mm.
Depth:	175 mm. (outer dimension)
Thread diameter:	8 mm.
Weight:	319 grams
MOQ:	5 pcs.

SHARPS CONTAINER HOLDER - SQUARED

For E.G. "Sharpsafe"

JBMEDICO



Sharps Container 4 and 7 liters

Sharps containers 4 / 7 liters
"Sharpsafe" container.

Yellow or red lid



JB 154-00-00

Sharps container holder for e.g. "Sharpsafe"
4 and 7 liters squared container.

Item no.:	JB 154-00-00
Inner diameter:	156 mm. x 156 mm.
Material:	Stainless steel (AISI 304)
Surface:	Electropolished
Height:	50 mm.
Depth:	208 mm. (outer dimension)
Thread diameter:	8 mm.
Weight:	280 grams
MOQ:	5 pcs.



CONTACT INFO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

SHARPS CONTAINER HOLDER - SQUARED For E.G. "Sharpsafe"



JB 279-00-00

Sharps container holder for e.g.
"Sharpsafe Quiver"
7,5 liters, squared container (long model)

Item no.: JB 279-00-00
Inner diameter: 101 mm. x 147 mm.
Material: Stainless steel (AISI 304)
Surface: Electropolished
Height: 50 mm.
Depth: 153 mm. (outer dimension)
Thread diameter: 8 mm.
Weight: 247 grams
MOQ: 5 pcs.

Sharps Container 7,5 Liters

Sharps container "Sharpsafe"
7.5 liters QUIVER.

Yellow or red lid.



SHARPS CONTAINER HOLDER - SQUARED For E.G. "Europlásticos"



JB 170-00-00

Sharps container holder for e.g. "Sharpsafe"
1 liter squared container.

Item no.: JB 170-00-00
Inner diameter: 96 mm. x 96 mm.
Material: Stainless steel (AISI 304)
Surface: Electropolished
Height: 50 mm.
Depth: 208 mm. (outer dimension)
Thread diameter: 8 mm.
Weight: 200 grams
MOQ: 5 pcs.

Sharps Container 1 Liters

Sharps container for 1 liter boks
from Europlásticos.

Yellow or red lid.



CONTACT INFO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

SHARPS CONTAINER HOLDER - SQUARED

Other Assortment - e.g. "Safebox"



Sharps Container 4L, 6L and 11 Liters

Sharps Container "Safebox" 4, 6 and 11 liters

Yellow or red lid.

Accessories:
Wall bracket with T-slot
Rail clamp with T-slot (DK)
Rail clamp with T-slot (EU)
Multibracket with T-slot
GCX rail bracket with T-slot



JB 279-00-00

Sharps Container Holder for e.g.
"Sharpsafe Quiver"

7,5 liters, squared container (long model)

Item no.:	JB 279-00-00
Inner diameter:	101 mm. x 147 mm.
Material:	Stainless steel (AISI 304)
Surface:	Electropolished
Height:	50 mm.
Depth:	153 mm. (outer dimension)
Thread diameter:	8 mm.
Weight:	247 grams
MOQ:	5 pcs.

SHARPS CONTAINER HOLDER - ROUND

Other Assortment

JBMEDICO



JB 150-00-00

Sharps Container Holder for
0,5- and 1 liter round container.

Fits e.g. Dispo 0,5 liter box from Grathwol.

Item no.:	JB 150-00-00
Inner diameter:	Ø73 mm.
Material:	Stainless steel (AISI 304)
Surface:	Electropolished
Height:	50 mm.
Depth:	112 mm. (outer dimension)
Thread diameter:	8 mm.
Weight:	150 grams
MOQ:	5 pcs.



JB 151-00-00

Sharps Container Holder for
0.6 liter - round container.

Fits e.g. Multisafe 1 liter and USON 500 ml.

Item no.:	JB 151-00-00
Inner diameter:	Ø89 mm.
Material:	Stainless steel (AISI 304)
Surface:	Electropolished
Height:	50 mm.
Depth:	134 mm. (outer dimension)
Thread diameter:	8 mm.
Weight:	171 grams
MOQ:	5 pcs.

CONTACT INFO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

SHARPS CONTAINER HOLDER - ROUND

Other Assortment

JBMEDICO



JB 264-00-00

Sharps Container Holder for e.g. Mediq /USON
5 and 7 liters round container.



JB 266-00-00

Sharps Container Holder for e.g. USON 11 liters and
Jokey container from Mediq.



JB 267-00-00

Sharps Container Holder for e.g. USON 21 liters and
Jokey container from Mediq.



Item no.: JB 264-00-00
Inner diameter: Ø200 mm.
Material: Stainless steel (AISI 304)
Surface: Electropolished
Height: 50 mm.
Depth: 244 mm. (outer dimension)
Thread diameter: 8 mm.
Weight: 286 grams
MOQ: 5 pcs.

Item no.: JB 266-00-00
Inner diameter: Ø251 mm.
Material: Stainless steel (AISI 304)
Surface: Electropolished
Height: 50 mm.
Depth: 295 mm. (outer dimension)
Thread diameter: 8 mm.
Weight: 345 grams
MOQ: 5 pcs.

Item no.: JB 267-00-00
Inner diameter: Ø280 mm.
Material: Stainless steel (AISI 304)
Surface: Electropolished
Height: 50 mm.
Depth: 324 mm. (outer dimension)
Thread diameter: 8 mm.
Weight: 368 grams
MOQ: 5 pcs.

Accessories:
Wall bracket with T-slot
Rail clamp with T-slot (DK)
Rail clamp with T-slot (EU)
Multibracket with T-slot
GCX rail bracket with T-slot


CONTACT INFO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

WIRE BASKETS AND WIRE PRODUCTS

Stainless steel

JBMEDICO



JB 161-00-00

Wire basket 6 liters, conical
with hook for T-slot.



JB 161-01-01

Wire basket 7 liters, conical
with hook for T-slot.



JB 161-02-01

Wire basket 15 liters, conical
with hook for T-slot.

Item no.: JB 161-00-00
Dimensions: Width x Depth x Height
231mm. x 161mm. x 242mm.
Material: Stainless steel (AISI 304)
Surface: Electropolished
Weight: 853 grams
Volumen: 6 liters
MOQ: 5 pcs.

Accessories for all 3 types of wire baskets:
Wall bracket with T-slot
Rail clamp with T-slot (DK)
Rail clamp with T-slot (EU)
Multibracket with T-slot
GCX rail bracket with T-spur

Item no.: JB 161-01-01
Dimensions: Width x Depth x Height
623mm. x 217mm. x 132mm.
Material: Stainless steel (AISI 304)
Surface: Electropolished
Weight: 667 grams
Volumen: 7 liters
MOQ: 5 pcs.

Item no.: JB 161-02-01
Dimensions: Width x Depth x Height
337mm. x 233mm. x 267mm.
Material: Stainless steel (AISI 304)
Surface: Electropolished
Weight: 959 grams
Volumen: 15 liters
MOQ: 5 pcs.

CONTACT INFO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

WIRE BASKETS AND WIRE ACCESSORIES

JBMEDICO

WALL BRACKET Max. Ø38 mm. tube



JB 161-03-00

Wire basket 25 liters, conical.



JB 281-00-00

Wall bracket for hanging 25 liters wire basket.



JB 162-00-00

Mask and tubing hook, with a guide for T-slot.



JB 146-45-38

Wall bracket for hanging heavy equipment.

Enables safe suspension of various hospital and care items.

Item no.: JB 161-03-00
Dimensions: Width x Depth x Height
348 mm. x 248mm. x 400 mm.
Material: Stainless steel (AISI 304)
Surface: Electropolished
Weight: 1.241 grams
Volumen: 25 liters
MOQ: 5 pcs.

Item no.: JB 281-00-00
Dimensions: Width x Depth x Height
23,5 mm. (above) 13,8 mm.
(below) x 21 mm. x 53,5 mm.
Material: Stainless steel (AISI 304)
Surface: Electropolished
Weight: 27 grams
Hole size: Ø4,25 mm.
MOQ: 40 pcs. in a bag

Item no.: JB 162-00-00
Dimensions: Width x Depth x Height
30 mm. x 76 mm. x 100 mm.
Material: Stainless steel (AISI 304)
Surface: Electropolished
Weight: 80 grams
MOQ: 10 pcs.

Item no.: JB 146-45-38
Dimensions: Width (above) or Width
(below) x Depth x Height
40 mm. (above) 20 mm.
(below) x 50 mm. x 93 mm.
Material: Stainless steel (AISI 304)
Surface: Electropolished
Weight: 117 grams
Hole size: Ø6 mm.
MOQ: 10 pcs. in a bag

CONTACT INFO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

BRACKET FOR SUCTION CANISTER AND CATHETER CONTAINER



JB 100-00-00

Bracket Ø100, for 2 liters
suction canister, Ø100 mm. hole.



JB 239-00-00

Catheter container
with guide for T-slot.

Item no.: JB 100-00-00
Dimensions: Ø100 mm. / Ø 8 mm. Thread
Material: Stainless steel (AISI 304)
Surface: Electropolished
Weight: 189 grams
MOQ: 5 pcs.

Item no.: JB 239-00-00
Height: 400 mm.
Material: Stainless steel (AISI 304)
Surface: Electropolished
Weight: 1.000 grams
Volumen: 2 liters
MOQ: 5 pcs.



CONTACT INFO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

VARIOUS ACCESSORIES



JB 92-00-00

Hygienic universal urine bottle holder.

Fits all types of hospital beds and all common bottles.
Machine washable and autoclavable.
Optimised fit and easy to clean.

Item no.:	JB 92-00-00
Dimensions:	Width x Depth x Height 123,6 mm. x 145 mm. x 192 mm.
Material:	Stainless steel (AISI 304)
Surface:	Electropolished
Weight:	853 grams
MOQ:	5 pcs.



JB 15-100-240

Storage box in stainless steel
with T-track bracket.

Item no.:	JB IP 9025
Material:	Stainless steel (Sling cleaned)
Height:	100 mm.
Length:	240 mm. (without T-track)
Width:	100 mm.
Weight:	1.2 kg.
MOQ:	2 pcs.



JB 115-00-00

Hospital journals holder in hard frosted plastic.

With handle for hanging on beds.
Instrument washable up to 90° C.
and fits all beds.

Item no.:	JB 115-00-00
Material:	Plastic PP (Polypropylene)
Weight:	800 grams
Colour:	White (Matte)
MOQ:	4 pcs.



Depth:	56 mm. without handle
Depth:	98 mm. with handle
Weight:	282 mm.
Height:	372 mm.

CONTACT INFO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

VARIOUS ACCESSORIES
DANISH Hospital Power Cords



DK Hospital Power Cords

	Item no.:
Angled 0.5 meter, red DK hospital + straight C13 appliance connector and shielded 3G1.0 mm. ² cable	1210715
Angled 1 meter Red DK hospital + Straight C13 appliance connector and shielded 3G1.0 mm. ² cable	1190110
Angled 2 meters Red DK hospital + Straight C13 appliance connector and shielded 3G1,0 mm. ² cable	1190111
Angled 3 meters Red DK hospital + Straight C13 appliance connector and shielded 3G1.0 mm. ² cable	1190112
Angled 4.5 metres Red DK hospital + Straight C13 appliance connector and shielded 3G1.0 mm. ² cable	1190113
Angled 8 meters Red DK hospital + Straight C13 device connector and shielded 3G1,5 mm. ² cable	1211755

	Item no.:
Angled 3,0 meters Red DK hospital + Straight C13 appliance connector and shielded 3G1.0 mm. ² cable	1212273
Angled 4,5 meters Red DK hospital + Straight C13 appliance connector and shielded 3G1.0 mm. ² cable	1212274



CONTACT INFO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

VARIOUS ACCESSORIES

Charging Towers or Power Towers



JB 680-400-15

Charging Towers or Power Towers
on wheels for mobility.

Varenr.: JB 680-400-15
Dimensioner: 400 mm. (Bottom plate)
Wheels: Ø 75 mm. (2 wheels wit brakes)
Max. load: 10A / 2.300 W.
Height: 762 mm.
Width x Depth: 165 x 135 mm. (Tower)
Material: White paineted steedal, Ral 9016
Sockets w. ground: 15 plugs.
Weight: 13 kg.
Cord: 8 meters



JB 680-400-15-UH

Charging Towers or Power Towers
without wheels.

Varenr.: JB 680-400-15-UH
Dimensioner: 400 mm. (Bottom plate)
Max. load: 10A / 2.300 W.
Height: 625 mm.
Width x Depth: 165 x 135 mm. (Tower)
Material: White paineted steedal, Ral 9016
Sockets w. ground: 15 plugs
Weight: 13 kg.
Cord: 8 meters



CONTACT INFO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

VARIOUS ACCESSORIES

Toilet brush holder



JB 350-00-00

Toilet brush holder for most round and oval brushes in plastic.

Item no.: JB 350-00-00
 Dimensions: Depth x Height
 140 mm. x 153 mm.
 Thickness: 5 mm.
 Material: Aluminium
 Surface: Anodized
 Colour: Metal
 Weight: 173 grams
 MOQ: 5 pcs.



Accessories:
 Wall brackets
 Wall bracket with
 hex allen screw for
 locking

Item no.
 JB 47-00-00
 JB 47-01-01

VARIOUS ACCESSORIES

Tubing holder



JB 292-00-00

IV tubing holder with a guide for t-slot.

Item no.: JB 292-00-00
 Dimensions: Depth x Height x Width
 71 mm. x 71 mm. 29 mm.
 Material: Aluminium
 Surface: Anodized
 Colour: Metal
 Weight: 388 grams
 MOQ: 10 pcs.

Accessories:
 - Clamp with T-slot for
 DK 10 x 30 mm. standard
 rails.
 - Clamp with T-slot for
 EU-DIN 10 x 25 mm.
 standard rails
 - Wall bracket with
 T-slot

Item no.
 JB 120-00-00
 JB 146-00-00
 JB 47-00-00



CONTACT INFO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

VARIOUS ACCESSORIES
Tubing holders



JB 600-00-01



Tube holder in red plastic
To be placed directly on the rail.
Fits DK standard rails only.

Item no.: JB 600-00-01
Dimensions: Depth x Height
140 mm. x 153 mm.
Material: POM
Colour: Red
Weight: 8,7 grams
MOQ: 40 pcs.

JB 600-10-25

Tube holder in red plastic
To be placed directly on the rail.
Fits EU DIN standard rails only.

Item no.: JB 600-10-25
Dimensions: Depth x Height
40 mm. x 60 mm.
Material: POM
Colour: Red
Weight: 8,7 grams
MOQ: 40 pcs.



CONTACT INFO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361



JB IP 9025

Emergency pull cord 2,5 meters.



JB IP 24155

"Bell" for emergency pull cord in red PP-Plastic



JB IP 9004

Break joint, black, for emergency pull cord.



JB IP 3000-RED

Emergency pull cord on cable reel

Pull force: 60 N.m (newton meters)

Item no.:	JB IP 9025
Material:	Plastic PE LDPE soft Polyethylene
Weight:	7 grams (182 grams)
MOQ:	25 stk.

Item no.:	JB IP 4155
Dimensions:	Diameter x Height 15,5 mm. (bottom) x 24 mm.
Material:	PP, Plastic (Polypropylen))
Weight:	1,4 grams (140 grams)
MOQ:	100 stk. / pose

Item no.:	JB IP 9004
Dimensions:	Length: 22 mm. Diameter: 9 mm.
Material:	Plastic ABS
Colour:	Black
Weight:	0,5 grams (100 grams)
MOQ:	200 stk.

Item no.:	JB IP 3000-RED
Item no.:	JB IP 1000-RED
Item no.:	JB IP 500-RED
Item no.:	JB IP 100-RED

CONTACT INFO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

TERMS OF SALE AND DELIVERY

JB-Medico ApS supplies private and public companies

Delivery:

Items usually arrive 1-5 business days after placing the order.

Products that need to be manufactured will of course have a longer delivery time.

Prices:

Unless otherwise stated, prices are pr. pcs. and in EURO

All government taxes are included our prices, with the exception of VAT.

All prices are pr. pcs.

Orders are taken at the prices applicable at the time of delivery.

We reserve the right to make adjustments and changes in prices, product specifications, tax changes, printing errors and image errors.

Payment:

Bank: Spar Nord Bank A/S

Account no.: 9004 / 1120732826

IBAN no: DK2390041120732826

SWIFT code: DABADKK

Transport damage:

Must be reported immediately upon receipt of goods, both to the carrier and to JB-Medico ApS.

A note must be made on the consignment note.

Custom-made goods and custom-made products, that have been specially taken home/manufactured for a customer, can not be returned.

Complaints:

Any error in the delivery must be reported immediately to customer service +45 7020 5826

Or info@jbmedico.com

Returns:

Upon receipt of goods from JB-Medico ApS, the customer must check that what is received corresponds correctly, to what was ordered and according to the enclosed delivery note.

Incorrect delivery must be reported immediately to JB-Medico ApS so that an exchange can take place.

All returns must be agreed with JB-Medico's dispatch department. Otherwise JB-Medico ApS may refuse receipt.

Returned goods must be returned in unbroken and undamaged original packaging.

In addition, the original packaging must be protected by an outer packaging during return transport. The goods must be cleaned.

Correspondence must be enclosed stating what is being returned, the reason for the return and with whom the return has been agreed.

CONTACT INFO

info@jbmedico.com | +45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

FAVORITES

Item no.

[illegible][illegible]

NOTES

	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64	65	66	67	68	69	70	71	72	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96	97	98	99	100	101	102	103	104	105	106	107	108	109	110	111	112	113	114	115	116	117	118	119	120	121	122	123	124	125	126	127	128	129	130	131	132	133	134	135	136	137	138	139	140	141	142	143	144	145	146	147	148	149	150	151	152	153	154	155	156	157	158	159	160	161	162	163	164	165	166	167	168	169	170	171	172	173	174	175	176	177	178	179	180	181	182	183	184	185	186	187	188	189	190	191	192	193	194	195	196	197	198	199	200	201	202	203	204	205	206	207	208	209	210	211	212	213	214	215	216	217	218	219	220	221	222	223	224	225	226	227	228	229	230	231	232	233	234	235	236	237	238	239	240	241	242	243	244	245	246	247	248	249	250	251	252	253	254	255	256	257	258	259	260	261	262	263	264	265	266	267	268	269	270	271	272	273	274	275	276	277	278	279	280	281	282	283	284	285	286	287	288	289	290	291	292	293	294	295	296	297	298	299	300	301	302	303	304	305	306	307	308	309	310	311	312	313	314	315	316	317	318	319	320	321	322	323	324	325	326	327	328	329	330	331	332	333	334	335	336	337	338	339	340	341	342	343	344	345	346	347	348	349	350	351	352	353	354	355	356	357	358	359	360	361	362	363	364	365	366	367	368	369	370	371	372	373	374	375	376	377	378	379	380	381	382	383	384	385	386	387	388	389	390	391	392	393	394	395	396	397	398	399	400	401	402	403	404	405	406	407	408	409	410	411	412	413	414	415	416	417	418	419	420	421	422	423	424	425	426	427	428	429	430	431	432	433	434	435	436	437	438	439	440	441	442	443	444	445	446	447	448	449	450	451	452	453	454	455	456	457	458	459	460	461	462	463	464	465	466	467	468	469	470	471	472	473	474	475	476	477	478	479	480	481	482	483	484	485	486	487	488	489	490	491	492	493	494	495	496	497	498	499	500	501	502	503	504	505	506	507	508	509	510	511	512	513	514	515	516	517	518	519	520	521
--	---	---	---	---	---	---	---	---	---	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----