

Product Catalogue

2025

SUSPENSION SOLUTIONS
FOR THE HEALTHCARE SECTOR

MADE OF NECESSITY AND DESIGNED TO LAST

JB MEDICO

Highlights

- + Suspension solutions, hygiene, drip and infusion stands
- + Various stainless steel dispensers
- + Design, development and production of new products
- + Assembly

CONTACTS

Production & Design

Jan Haurum Kristensen
Phone: +45 7020 5826
E-mail: jan@jbmedico.com

Purchasing & Warehouse:

Jesper A. Berg
Phone: +45 7020 5826
Mobile: +45 2868 5826
Area: Denmark and Exports
E-mail: jb@jbmedico.com

Sales:

Ana Belén Martín
Mobile: +45 9386 3558
E-mail: abm@jbmedico.com

ENVIRONMENTAL OBJECTIVES

JB Medico develops and produces high quality products with strong materials to increase the lifetime of the products. In this way we spare the environment in the long term.

Having a Green Environmental Profile and constantly focusing on sustainability means a lot to us and is something we deal with on a daily basis and constantly try to improve. We think we have come a long way, but there are of course places where we can improve even more.

Packaging:

We pack a lot of parcels every day and therefore use a lot of packaging in the form of boxes and stuffing. So this is one of the areas where we focus a lot on making it as sustainable as possible. We like to recycle cardboard boxes and fillings that still have a certain quality. So when our customers receive a box that looks used or has a logo from someone other than us, it's simply an attempt to protect the environment as much as possible.

- FSC Eco-certified cardboard and paper from responsible sources. (www.fsc.org)
- Minimal use of Plastic for box filling. Otherwise it is recycling from received shipments.
- Recycling of cardboard boxes and filling.

Waste sorting:

We recycle for example:

- Paper and Cardboard.
- Packaging.
- Soft Plastics.
- Hard Plastic.
- Metal.
- Wiring.

Our hope for our customers:

We would be delighted if you would help us to take care of our planet by sorting and returning the packaging you receive from us to the recycling centre. Thank you!

Thanks in advance!

With best regards
Team JB Medico

JBMEDICO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

TABLE OF CONTENTS

Product Catalogue	1	SHARPS CONTAINER HOLDER - ROUND	62
CONTACTS	2	SHARPS CONTAINER HOLDER - OVAL	63
ENVIRONMENTAL OBJECTIVES	2	SHARPS CONTAINER HOLDER - SQUARED	64
MEDICAL RAIL CLAMPS and SLIDING CLAMPS	4	SHARPS CONTAINER HOLDER - ROUND	66
MEDICAL RAIL CLAMPS FOR EQUIPMENT RAILS	5	WIRE BASKETS AND WIRE PRODUCTS	68
RAIL CLAMPS	6	WIRE BASKETS AND WIRE ACCESSORIES	70
MEDICAL RAIL CLAMPS - SLIDING CLAMPS	10	WALL BRACKET	70
RAIL CLAMPS FOR DRIP STANDS, IT AND MEDICAL EQUIPMENT	13	BRACKET FOR SUCTION CANISTER AND CATHETER CONTAINER	71
RAIL CLAMPS - EXTRA SOLID - FOR LONG AND HEAVY EXAMINATION LIGHTS	15	VARIOUS ACCESSORIES	72
ACCESSORIES FOR VARIOUS RAIL CLAMPS	16	EMERGENCY PULL CORD + ACCESSORIES	77
MULTIBRACKET - ADJUSTABLE - 16-41 MM.	17	VARIOUS ACCESSORIES	78
SUSPENSION HOUSING FOR RAIL	19	TERMS OF SALE AND DELIVERY	79
CLAMPS AND MULTIBRACKET	19	FAVORITES	80
FITTINGS SUITABLE FOR GCX RAILS	20	NOTES	80
MEDICAL EQUIPMENT RAILS	24		
ACCESSORIES FOR MEDICAL RAILS	25		
WALL BRACKETS FOR BOTH RAIL STANDARDS	26		
WALL BRACKET FOR PROTECTION OF MEDICAL EQUIPMENT RAILS AND MORE	27		
MEDICAL LIGHTING FROM GLAMOX LUXO	28		
ELECTRIC PILL CRUSHER V3 - SEVERO	31		
ACCESSORIES FOR SEVERO ELECTRIC PILL CRUSHER V3	32		
MONITOR ARMS AND ACCESSORIES	33		
WALL BRACKETS	34		
ACCESSORIES FOR MONITOR ARMS	35		
MONITOR ARMS AND ACCESSORIES	36		
ROD SUSPENSION	38		
SOAP AND ALCOHOL DISPENSERS	39		
DRIP TRAYS FOR SOAP AND ALCOHOL DISPENSERS	41		
DISPENSERS WITH DRIP TRAY AND ADAPTER BRACKET	42		
TABLE AND WALL DISPENSER FOR SQUARE AND ROUND BOTTLES	44		
GLOVE, WET WIPE AND SANITARY FACEMASK HOLDER	45		
WALL DISPENSER FOR WET WIPES	46		
ACCESSORIES FOR GLOVE BOX HOLDERS	47		
ALUMINIUM SLIDING BRACKET FOR T-SLOT RAILS	48		
GLOVE BOX HOLDER	49		
IV STANDS / IV POLES	50		
IV POLE XL	51		
IV POLE XXL	52		
IV STAND / IV POLE	53		
IV STAND / INFUSION POLE - SUPPORT COLUMNS	54		
INFUSION AND IV POLE - ACCESSORIES	55		
WALL BRACKET FOR CO2 BOTTLES, Ø183 MM. HOLE	60		
OXYGEN AND GAS CYLINDERS HOLDERS	61		

MEDICAL RAIL CLAMPS and SLIDING CLAMPS

EU DIN (10x25 mm) & Scandinavian (10x30 mm) standard
for medical equipment rails



Medical rail clamps for 10x25 & 10x30 mm. rails

All our medical rail clamps fit EU DIN (10x25mm) & Scandinavian (10x30 mm) standard rails.

When planning and ordering, it is important for your to know which of our systems you would like to use.

JB-Medico offers several different rail variants. We have developed two own rails 10x25 mm. and 10x30 mm., which can handle most of the cleats available on the market.

JB-Medico customises solutions for all clinical departments.

JBMEDICO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361


MEDICAL RAIL CLAMPS FOR EQUIPMENT RAILS

EU DIN standard 10x25 mm.




JB 146-00-00

Rail clamp, narrow model, with one ball lock and T-slot.



JB 147-00-00

Rail clamp narrow model, with one ball lock.



JB 147-03-00

Rail clamp, narrow model, with one ball lock and three countersunk Ø6,6 holes.

"Prepared for suspension."



JB 143-00-00

Rail clamp, wide model, locked using two hex allen screws.

Item no.: **JB 146-00-00**
Material: **Aluminium (Ceramic Coating)**
Dimensions: **Height x Width x Depth: 61,5 mm. x 64 mm. x 28,6 mm.**
Weight: **132 grams**
MOQ: **5 pcs.**

Item no.: **JB 147-00-00**
Material: **Aluminium (Ceramic Coating)**
Dimensions: **Height x Width x Depth: 61,5 mm. x 64 mm. x 28,6 mm.**
Weight: **168 grams**
MOQ: **5 pcs.**

Item no.: **JB 147-03-00**
Material: **Aluminium (Ceramic Coating)**
Dimensions: **Height x Width x Depth: 61,5 mm. x 64 mm. x 28,6 mm.**
Weight: **161 grams**
MOQ: **5 pcs.**

Item no.: **JB 143-00-00**
Material: **Aluminium (Ceramic Coating)**
Dimensions: **Height x Width x Depth: 61,5 mm. x 90 mm. x 28,6 mm.**
Weight: **308 grams**
MOQ: **5 pcs.**

JBMEDICO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

RAIL CLAMPS

EU DIN standard for 10x25 mm. clamp rails



JB 143-03-00

Rail clamp, wide model, locked using two hex allen screws. It has three countersunk holes Ø6,6 mm.

"Prepared for suspension."

Item no.: **JB 143-03-00**
Material: **Aluminium (Ceramic Coating)**
Dimensions: **Height x Width x Depth:**
61,5 mm. x 90 mm. x 28,6 mm.
Weight: **303 grams**
MOQ: **5 pcs.**



JB 144-00-00

Rail clamp, wide model, locked with two ball locks.

Item no.: **JB 144-00-00**
Material: **Aluminium (Ceramic Coating)**
Dimensions: **Height x Width x Depth:**
61,5 mm. x 110 mm. x 28,6 mm.
Weight: **328 grams**
MOQ: **5 pcs.**



JB 144-03-00

Rail clamp, wide model, locked with two ball locks. It has three countersunk holes Ø6,6 mm.

"Prepared for suspension."

Item no.: **JB 144-03-00**
Material: **Aluminium (Ceramic Coating)**
Dimensions: **Height x Width x Depth:**
61,5 mm. x 110 mm. x 28,6 mm.
Weight: **318 grams**
MOQ: **5 pcs.**



JB 145-00-00

Rail clamp, wide model, locked with one ball lock.

Item no.: **JB 145-00-00**
Material: **Aluminium (Ceramic Coating)**
Dimensions: **Height x Width x Depth:**
61,5 mm. x 110 mm. x 28,6 mm.
Weight: **325 grams**
MOQ: **5 pcs.**

JBMEDICO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

RAIL CLAMPS

EU DIN standard for 10x25 mm. clamp rails



JB 145-03-00

Rail clamp, wide model, locked with one ball lock. It has three countersunk holes Ø6,6 mm.

"Prepared for suspension."

Item no.: **JB 145-03-00**
Material: **Aluminium (Ceramic Coating)**
Dimensions: **Height x Width x Depth:**
61,5 mm. x 110 mm. x 28,6 mm.
Weight: **317 grams**
MOQ: **5 pcs.**



JB 045-03-00

Rail clamp, narrow model, locked with two socket screws. It has three countersunk Ø6,6 mm. holes.

"Prepared for suspension."

Item no.: **JB 045-03-00**
Material: **Aluminium (Ceramic Coating)**
Dimensions: **Height x Width x Depth:**
61,5 mm. x 45 mm. x 28,6 mm.
Weight: **151 grams**
MOQ: **5 pcs.**



JB 143-03-18

Rail clamp, wide model, locked using two hex allen screws and adapter fitting Ø18 mm.

Three bush sizes are available.
Handles can be purchased.

Item no.: **JB 143-03-18**
Material: **Aluminium (Ceramic Coating)**
Dimensions: **Lenght x Width x Depth:**
61,5 mm. x 90 mm. x 64,2 mm.
Weight: **468 grams**
MOQ: **2 pcs.**



JB 145-03-18

Rail clamp, wide model, 1 pcs. ball lock and adapter fitting Ø18 mm.

Three bush sizes are available.
Handles can be purchased.

Item no.: **JB 145-03-18**
Material: **Aluminium (Ceramic Coating)**
Dimensions: **Lenght x Width x Depth:**
61,5 mm. x 110 mm. x 64,2 mm.
Weight: **501 grams**
MOQ: **2 pcs.**

JBMEDICO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

LAMP RAIL CLAMPS - EXTRA SOLID - FOR LONG AND HEAVY LAMPS

EU DIN Standard for 10x25 mm.



JB 144-03-18

Rail clamp, wide model, 2 pcs. ball locks and adapter fitting Ø18 mm.

Three bush sizes are available.
Handles can be purchased (JB 94-00-00)



JB 147-03-18

Rail clamp, narrow model, 1 pcs. ball locks and adapter fitting Ø18 mm.

Three bush sizes available.
Handles can be purchased (JB 94-00-00)



JB 143-03-20

Rail clamp, wide model, locked using two hex allen screws and adapter fitting Ø18 mm.

Handles can be purchased (JB 94-00-00)



JB 145-03-20

Rail clamp, wide model, with a one-ball lock, adapter fitting and brass bush Ø20 mm. hole.

Handles can be purchased (JB 94-00-00)

Item no.: **JB 144-03-18**
Material: **Aluminium (Ceramic Coating)**
Dimensions: **Lenght x Width x Depth:
61,5 mm. x 110 mm. x 64,2 mm.**
Weight: **502 grams**
MOQ: **2 pcs.**

Item no.: **JB 147-03-18**
Material: **Aluminium (Ceramic Coating)**
Dimensions: **Lenght x Width x Depth:
61,5 mm. x 64 mm. x 64,2 mm.**
Weight: **343 grams**
MOQ: **4 pcs.**

Item no.: **JB 143-03-20**
Material: **Aluminium (Ceramic Coating)**
Dimensions: **Lenght x Width x Depth:
61,5 mm. x 90 mm. x 60,6 mm.**
Weight: **536 grams**
MOQ: **2 pcs.**

Item no.: **JB 145-03-20**
Material: **Aluminium (Ceramic Coating)**
Dimensions: **Lenght x Width x Depth:
61,5 mm. x 110 mm. x 64,2 mm.**
Weight: **550 grams**
MOQ: **2 pcs.**

JBMEDICO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

SUSPENSION RAIL CLAMPS FOR DROP STANDS, IT AND MEDICAL EQUIPMENT

EU DIN Standard for 10x25 mm.



JB 144-03-20

Rail clamp, wide model, with a two-ball lock, adapter fitting and brass bush Ø20 mm. hole.

Handles can be purchased (JB 94-00-00)

Item no.: **JB 144-03-20**
Material: **Aluminium (Ceramic Coating)**
Dimensions: **Lenght x Width x Depth:**
61,5 mm. x 110 mm. x 64,2 mm.

Weight: **551 grams**
MOQ: **2 pcs.**



JB 147-03-20

Rail clamp, narrow model, with a one-ball lock, adapter fitting and brass bush Ø20 mm. hole.

Handles can be purchased (JB 94-00-00)

Item no.: **JB 147-03-20**
Material: **Aluminium (Ceramic Coating)**
Dimensions: **Lenght x Width x Depth:**
61,5 mm. x 64 mm. x 64,2 mm.

Weight: **550 grams**
MOQ: **2 pcs.**



JB 045-03-20

Rail clamp, narrow model, with two hex allen screws, adapter fitting and brass bush Ø20 mm. hole.

Handles can be purchased (JB 94-00-00)

Item no.: **JB 045-03-20**
Material: **Aluminium (Ceramic Coating)**
Dimensions: **Lenght x Width x Depth:**
61,5 mm. x 45 mm. x 64,2 mm.

Weight: **384 grams**
MOQ: **2 pcs.**



JBMEDICO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

MEDICAL RAIL CLAMPS - SLIDING CLAMPS

Scandinavian standard for 10x30 mm. medical equipment rails



JB 88-00-00

Rail clamp, wide model, locked using two hex allen screws.



JB 103-00-00

Rail clamp, wide model, locked with one ball lock.



JB 104-00-00

Rail clamp, wide model, locked with two ball locks.



JB 207-00-00

Rail clamp, wide model, locked with two ball locks. It has three countersunk holes Ø6,6 mm.

"Prepared for suspension."

Item no: **JB 88-00-00**
Material: **Aluminium (Ceramic coating)**
Dimensions: **Length x Width x Depth: 61,5 mm. x 90 mm. x 25 mm.**
Weight: **280 grams**
MOQ: **5 pcs.**

Item no: **JB 103-00-00**
Material: **Aluminium (Ceramic coating)**
Dimensions: **Length x Width x Depth: 61,5 mm. x 110 mm. x 25 mm.**
Weight: **295 grams**
MOQ: **5 pcs.**

Item no: **JB 104-00-00**
Material: **Aluminium (Ceramic coating)**
Dimensions: **Length x Width x Depth: 61,5 mm. x 110 mm. x 25 mm.**
Weight: **304 grams**
MOQ: **5 pcs.**

Item no: **JB 207-00-00**
Material: **Aluminium (Ceramic coating)**
Dimensions: **Length x Width x Depth: 61,5 mm. x 110 mm. x 25 mm.**
Weight: **299 grams**
MOQ: **5 pcs.**

JBMEDICO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

MEDICAL RAIL CLAMPS - SLIDING CLAMPS

Scandinavian standard for 10x30 mm. medical equipment rails



JB 206-00-00

Rail clamp, wide model, locked using two hex allen screws. It has three countersunk holes Ø6,6 mm.

"Prepared for suspension."

Item no: **JB 206-00-00**
Material: **Aluminium (Ceramic coating)**
Dimensions: **Lenght x Width x Depth: 61,5 mm. x 90 mm. x 25 mm.**
Weight: **275 grams**
MOQ: **5 pcs.**




JB 217-00-00

Rail clamp, wide model, locked with one ball lock. It has three countersunk holes Ø6,6 mm.

"Prepared for suspension."


Item no: **JB 217-00-00**
Material: **Aluminium (Ceramic coating)**
Dimensions: **Lenght x Width x Depth: 61,5 mm. x 110 mm. x 25 mm.**
Weight: **293 grams**
MOQ: **5 pcs.**



JB 120-00-00

Rail clamp, narrow model, locked with one ball lock, and a T-slot slot.

Item no: **JB 120-00-00**
Material: **Aluminium (Ceramic coating)**
Dimensions: **Lenght x Width x Depth: 61,5 mm. x 64 mm. x 25 mm.**
Weight: **118 grams**
MOQ: **5 pcs.**



JB 121-00-00

Rail clamp, narrow model, locked with one ball lock.

Item no: **JB 121-00-00**
Material: **Aluminium (Ceramic coating)**
Dimensions: **Lenght x Width x Depth: 61,5 mm. x 64 mm. x 25 mm.**
Weight: **155 grams**
MOQ: **5 pcs.**

JBMEDICO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

MEDICAL RAIL CLAMPS - SLIDING CLAMPS

Scandinavian standard for 10x30 mm. medical equipment rails



JB 121-03-00

Rail clamp, narrow model, locked with one ball lock. It has three countersunk Ø6,6 mm. holes.

"Prepared for suspension."



JB 216-03-00

Rail clamp, narrow model, locked with two socket screws. It has three countersunk Ø6,6 mm. holes.

"Prepared for suspension."



JB 176-00-180

Rail Clamp, an extra wide model, is locked using four hex allen screws.



JB 176-00-180 supporting a transformer

Item no.: **JB 121-03-00**
Material: **Aluminium (Ceramic coating)**
Dimensions: **Lenght x Width x Depth:
61,5 mm. x 64 mm. x 25 mm.**
Weight: **150 grams**
MOQ: **5 pcs.**

Item no.: **JB 216-03-00**
Material: **Aluminium (Ceramic coating)**
Dimensions: **Lenght x Width x Depth:
61,5 mm. x 45 mm. x 25 mm.**
Weight: **132 grams**
MOQ: **5 pcs.**

Item no.: **JB 176-00-180**
Material: **Aluminium (Ceramic coating)**
Dimensions: **Lenght x Width x Depth:
61,5 mm. x 180 mm. x 25 mm.**
Weight: **560 grams**
MOQ: **4 pcs.**

JBMEDICO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

RAIL CLAMPS FOR DRIP STANDS, IT AND MEDICAL EQUIPMENT

Scandinavian standard 10x30 mm.



JB 206-03-20

Rail clamp, wide model, locks with two hex allen screws and adapter bracket Ø20 mm.

Three bush sizes available.
Handles can be purchased (JB 94-00-00)

Item no.:	JB 206-03-20
Material:	Aluminium (Ceramic Coating)
Dimensions:	Length x Width x Depth: 61,5 mm. x 90 mm. x 60,5 mm.
Weight:	507 grams
MOQ:	2 pcs.

JBMEDICO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

RAIL CLAMPS FOR DRIP STANDS, IT AND MEDICAL EQUIPMENT

Scandinavian standard 10x30 mm.



JB 217-03-20

Rail clamp, wide model, with a one-ball lock, adapter fitting and brass bush Ø20 mm. hole.

Three bush sizes available.
Handles can be purchased (JB 94-00-00)

Item no.: **JB 217-03-20**
Material: **Aluminium (Ceramic Coating)**
Dimensions: **Length x Width x Depth:**
61,5 mm. x 110 mm. x 60,5 mm.
Weight: **530 grams**
MOQ: **2 pcs.**



JB 207-03-20

Rail clamp, wide model, with a two-ball lock, adapter fitting and brass bush Ø20 mm. hole.

Three bush sizes available.
Handles can be purchased (JB 94-00-00)

Item no.: **JB 207-03-20**
Material: **Aluminium (Ceramic Coating)**
Dimensions: **Length x Width x Depth:**
61,5 mm. x 110 mm. x 60,5 mm.
Weight: **537 grams**
MOQ: **2 pcs.**



JB 216-03-20

Rail clamp, narrow model, with two hex allen screws, adapter fitting and brass bush Ø20 mm. hole.

Three bush sizes available.
Handles can be purchased (JB 94-00-00)

Item no.: **JB 216-03-20**
Material: **Aluminium (Ceramic Coating)**
Dimensions: **Length x Width x Depth:**
61,5 mm. x 45 mm. x 60,5 mm.
Weight: **368 grams**
MOQ: **4 pcs.**



JB 121-03-20

Rail clamp, narrow model, with a one-ball lock, adapter fitting and brass bush Ø20 mm. hole.

Three bush sizes available.
Handles can be purchased (JB 94-00-00)

Item no.: **JB 121-03-20**
Material: **Aluminium (Ceramic Coating)**
Dimensions: **Length x Width x Depth:**
61,5 mm. x 64 mm. x 60,5 mm.
Weight: **383 grams**
MOQ: **4 pcs.**

JBMEDICO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

RAIL CLAMPS - EXTRA SOLID - FOR LONG AND HEAVY EXAMINATION LIGHTS

Scandinavian standard for 10x30 mm. medical equipment rails



JB 206-03-18

Rail clamp, wide model, locked using two hex allen screws and adapter fitting Ø18 mm.

Three bush sizes are available.
Handles can be purchased.

Item no.: **JB 206-03-18**
Material: **Aluminium (Ceramic Coating)**
Dimensions: **Length x Width x Depth:**
61,5 mm. x 90 mm. x 60,6 mm.
Weight: **468 grams**
MOQ: **2 pcs.**



JB 217-03-18

Rail clamp, wide model, 1 pcs. ball lock and adapter fitting Ø18 mm.

Three bush sizes are available.
Handles can be purchased.

Item no.: **JB 217-03-18**
Material: **Aluminium (Ceramic Coating)**
Dimensions: **Length x Width x Depth:**
61,5 mm. x 110 mm. x 60,6 mm.
Weight: **475 grams**
MOQ: **2 pcs.**



JB 207-03-18

Rail clamp, wide model, 2 pcs. ball locks and adapter fitting Ø18 mm.

Three bush sizes are available.
Handles can be purchased.

Item no.: **JB 207-03-18**
Material: **Aluminium (Ceramic Coating)**
Dimensions: **Length x Width x Depth:**
61,5 mm. x 110 mm. x 60,6 mm.
Weight: **482 grams**
MOQ: **2 pcs.**



JB 216-03-18

Rail clamp, narrow model, locks with 2 pcs. hex allen screws and adapter fitting Ø18 mm.

Three bush sizes are available.
Handles can be purchased.

Item no.: **JB 216-03-18**
Material: **Aluminium (Ceramic Coating)**
Dimensions: **Length x Width x Depth:**
61,5 mm. x 45 mm. x 60,6 mm.
Weight: **316 grams**
MOQ: **4 pcs.**

JBMEDICO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

ACCESSORIES FOR VARIOUS RAIL CLAMPS

JB 128-00-09

POM bushing with bushing hole Ø 9.6 mm.

Item no.: JB 128-00-09
Material: Plastic (White)
Weight: 12 grams
MOQ: 10 pcs.



JB 22-00-05

Brass bushing with bushing hole Ø 20 mm.

Item no.: JB 22-00-05
Material: Brass
Weight: 88 grams
MOQ: 10 pcs.



JB 128-00-01

POM bushing with bushing hole Ø 13 mm.

Item no.: JB 128-00-01
Material: Plastic (White)
Weight: 8 grams
MOQ: 10 pcs.



JB 94-00-00

Handle with M6 x 20 mm. thread

Item no.: JB 94-00-00
Material: Stainless steel (AISI 304)
Weight: 42 grams
MOQ: 5 pcs.



JB 127-00-01

POM bushing with bushing hole Ø 16 mm.

Item no.: JB 127-00-01
Material: Plastic (White)
Weight: 3,6 grams
MOQ: 10 pcs.



JB 10-32-00

Handle with M10 x 40 mm. thread

Item no.: JB 10-32-00
Material: Stainless steel (AISI 304)
Weight: 170 grams
MOQ: 2 pcs.



JBMEDICO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

MULTIBRACKET - ADJUSTABLE - 16-41 MM.




JB 158-00-00

Multibracket - The fitting can be connected anywhere. Fits Ø16-41 mm.

This multibracket makes it possible to assemble various equipment together.

Item no.: **JB 158-00-00**
 Material: **Aluminium (Ceramic coating), SS 304, POM og Nylon (PGA)**
 Dimensions: **Lenght x Width x Depth: 36 mm. x 51 mm. x 134/96 mm.**
 Weight: **209 grams**
 MOQ: **4 pcs.**




JB 158-00-05

Multibracket with 2xØ6 clearance holes. Fits up to. Ø16-41 mm.

This multibracket makes it possible to assemble various equipment together.

Item no.: **JB 158-00-05**
 Material: **Aluminium (Ceramic coating), SS 304, POM og Nylon (PGA)**
 Dimensions: **Lenght x Width x Depth: 36 mm. x 51 mm. x 134/96 mm.**
 Weight: **204 grams**
 MOQ: **4 pcs.**



JB 158-06-05

Double multibracket for mounting of for example bed and drop stand.

Includes 2 pcs. Unbrako bolt screws.

Item no.: **JB 158-06-05**
 Material: **Aluminium (Ceramic coating), SS 304, POM og Nylon (PGA)**
 Weight: **421 grams**
 MOQ: **2 pcs.**



JB 158-00-163

Multibracket with T-slot fitting fits from Ø16-41 mm. Fitting can be placed on table, trolley, drop stand, DIN rail and columns.

Includes 2 pcs. hex allen screws.

Item no.: **JB 158-00-163**
 Material: **Aluminium (Ceramic coating), SS 304, POM og Nylon (PGA)**
 Dimensions: **Lenght x Width x Depth: 61,5 mm. x 64 mm. x 60,5 mm.**
 Weight: **245 grams**
 MOQ: **4 pcs.**

JBMEDICO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

MULTIBRACKET - ADJUSTABLE - 16-41 MM.



JB 158-03-18

Multibracket with a fixing device, Ø18mm. hole.
Fits Ø16-41 mm.

Three bush sizes available.
Handle can be purchased separately.

Item no.: **JB 158-03-18**
Material: **Aluminium (Ceramic coating),
SS 304, POM, Nylon (PA6)**
Dimensions: **Length x Width x Depth:
36 mm. x 51 mm. x 134/96 mm.**
Weight: **288 grams**
MOQ: **2 pcs.**

ACCESSORIES: Handle can be purchased.



JB 158-03-20

Multibracket with a fixing
device and brass bushing.
Fits up to. Ø16-41 mm.

Handle can be purchased separately.

Item no.: **JB 158-03-20**
Material: **Aluminium (Ceramic coating),
SS 304, Brass, POM, Nylon (PA6)**
Dimensions: **Length x Width x Depth:
36 mm. x 51 mm. x 134/96 mm.**
Weight: **440 grams**
MOQ: **2 pcs.**

ACCESSORIES: Handle and T-slot fitting can be
purchased.



JB 163-00-01

T-slot fitting for multi-bracket.

The holes are only suitable for multibracket.

Item no.: **JB 163-00-01**
Material: **Aluminium (Ceramic coating)**
Dimensions: **Height x Width x Depth:
60 mm. x 36,2 mm. x 12,5 mm.**
Weight: **36 grams**
MOQ: **1 bag of 20 pcs.**

JBMEDICO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

SUSPENSION HOUSING FOR RAIL CLAMPS AND MULTIBRACKET



JB 205-00-00

Fixing device for rail clamps and multi-brackets. Hole dimension: Ø25 mm.

Prepared for suspension.

Item no.: JB 205-00-00
Material: Aluminium (Ceramic coating)
Dimensions: Length x Width x Depth:
36 mm. x 35,6 mm. x 61 mm.
Weight: 133 grams
MOQ: 2 pcs.

ACCESSORIES: Handle and bushing is purchased separatley.



JB 205-00-18

Fixing device for rail clamps and multi brackets. Hole dimension: Ø18 mm.

Prepared for suspension.

Item no.: JB 205-00-18
Material: Aluminium (Ceramic coating)
Dimensions: Length x Width x Depth:
36 mm. x 35,6 mm. x 61 mm.
Weight: 171 grams
MOQ: 2 pcs.

ACCESSORIES: Handle and bushing is purchased separatley.



JBMEDICO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

FITTINGS SUITABLE FOR GCX RAILS



JB 820-625-06

GCX rail fitting, prepared with mounting holes.

Item no.: **JB 820-625-06**
 Material: **Aluminium (Ceramic coating)**
 Colour: **Metal**
 Dimensions: **Height x Width x Depth: 82 mm. x 63 mm. x 6 mm.**
 Weight: **83 grams**
 MOQ: **4 pcs.**

ACCESSORIES:

JB 163-00-01 **T-slot fitting**
 JB 205-00-00 **Adapter fitting/ Suspension housing Ø20**
 JB 205-00-18 **Adapter fitting/ suspension housing Ø18**



JB 820-625-163

GCX rail fitting with a T-slot.

Item no.: **JB 820-625-163**
 Material: **Aluminium (Ceramic coating)**
 Colour: **Metal**
 Dimensions: **Height x Width x Depth: 82 mm. x 36 mm. x 18,5 mm.**
 Weight: **123 grams**
 MOQ: **4 pcs.**



JB 820-625-20

GCX rail fitting with suspension housing. Ø20 mm. bushing hole.

Item no.: **JB 820-625-20**
 Material: **Aluminium (Ceramic coating) and Brass**
 Colour: **Metal**
 Dimensions: **Height x Width x Depth: 82 mm. x 63 mm. x 42 mm.**
 Weight: **313 grams**
 MOQ: **2 pcs.**

ACCESSORIES: Handle can be purchased.



POM bushing not included as standard

JB 820-625-18

GCX rail fitting with suspension housing Ø18 mm. hole.

Optional plastic bushing Ø9,6/13 or 16 mm.

Item no.: **JB 820-625-18**
 Material: **Aluminium (Ceramic coating) and Plastic**
 Colour: **Metal**
 Dimensions: **Height x Width x Depth: 82 mm. x 63 mm. x 42 mm.**
 Weight: **270 grams**
 MOQ: **2 pcs.**

ACCESSORIES: Handle can be purchased.

JBMEDICO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

FITTINGS SUITABLE FOR GCX RAILS




JB 820-625-20
For flat screen / monitor holder.



JB 820-625-163
For catheter holder.



JB 820-625-18
For light fixture.



JB 820-625-163
For sharpsafe holder.

JBMEDICO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

FITTINGS SUITABLE FOR GCX RAILS



JB 820-625-20
For monitor holder.



JB 820-625-163
For wire basket.



JB 820-625-20
Fitting with suspension housing.




JB 820-625-163
For catheter holder.

JBMEDICO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

MEDICAL EQUIPMENT RAILS AND ACCESSORIES

Rail Scandinavian 10x30 mm. standard



JBM 100-01-06

JBM medical equipment rail in hygienic design.

With increased Depth to avoid collision with cable trays. This allows suspension of medical equipment.

Item no.:	JBM 100-01-06
Material:	Aluminium (Ceramic coating)
Dimensions:	Length - Standard 6 meters
	Depth: 81 mm.
	Height: 93 mm.
Weight:	3,8 kg pr. meter
Standard length:	6 meters
End caps:	Plastic, easily mounted with wetroom silicone or screws
Screw holes:	Unplugged with plugs

Benefits include:

- 1. Very easy to clean**
- 2. Optimal hygiene, the rail has smooth surfaces, no cracks or transitions that can be dirt collectors.**
- 3. A steep slope, 25 degrees, so the rail cannot be used as a shelf.**
- 4. Mounting holes are pre-drilled pr. 100 mm.**




JBMEDICO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

MEDICAL EQUIPMENT RAILS

EU 10x25mm & Scandinavian standard 10x30 mm



JB 252-00-00

10x25 mm. medical equipment rail.

EU DIN standard.



JB 249-00-00

End cap.

For EU DIN 10x25 mm. rail profile.



JB 234-00-01

10x30 mm. medical equipment rail.

Scandinavian standard.



JB 248-00-00

End cap.

For 10x30 mm. Scandinavian standard rail profile.

Item no.: **JB 252-00-00**
Material: **Aluminium**
Surface: **Natural Anodized Metal**
Colour: **Metal**
Dimensions: **Height x Depth: 25 mm. x 10 mm.**
Weight: **700 grams pr. m.**
Standard Længde: **6 meters**

Cutting:
Drilling Holes:

Item no.: **JB 249-00-00**
Material: **Plastic (PE)**
Colour: **Grey**
Dimensions: **Height x Width 25 mm. x 10 mm.**
MOQ: **50 pcs.**

Item no.: **JB 234-00-01**
Material: **Aluminium**
Surface: **Natural Anodized Metal**
Colour: **Metal**
Dimensions: **Height x Depth: 30 mm. x 10 mm.**
Weight: **800 grams pr. meter**
Standard Længde: **6 meters**

Cutting:
Drilling Holes:

Item no.: **JB 248-00-00**
Material: **Plastic (PE)**
Colour: **Grey**
Dimensions: **Height x Width 30 mm. x 10 mm.**
MOQ: **50 pcs.**

JBMEDICO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

ACCESSORIES FOR MEDICAL RAILS

Scandinavian standard



JB M 100-01-07

End cap for JBM rail, right side.

Item no.: JB 100-01-07
Material: ABS Plastic
Colour: Grey
Dimensions: Height x Width
93 mm. x 81 mm.
Weight: 6 grams
MOQ: 50 pcs.

JB M 100-01-08

End cap for JBM rail, left side.

Item no.: JB 100-01-08
Material: ABS Plastic
Colour: Grey
Dimensions: Height x Width
93 mm. x 81 mm.
Weight: 6 grams
MOQ: 50 pcs.



JB M 100-01-10

Screw plugs for JBM rail.

Item no.: JB 100-01-10
Material: ABS Plastic
Colour: Grey
Dimensions: Diameter
Ø12 mm.
MOQ: 50 pcs.



JBMEDICO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

WALL BRACKETS FOR BOTH RAIL STANDARDS



JB 286-00-00

Wall or equipment fitting for two types of medical equipment rails.

Item no.: **JB 286-00-00**
Material: **Stainless steel (AISI 304)**
Colour: **Metal**
Dimensions: **Diameter x Depth**
Ø50 mm. x 10 mm.
Weight: **137 grams**
MOQ: **10 pcs.**



JB 400-00-00

Wall or equipment fitting for two types of medical equipment rails.

Item no.: **JB 400-00-00**
Material: **Aluminium**
Colour: **Metal**
Dimensions: **Diameter x Depth**
Ø50 mm. x 10 mm.
Weight: **48 grams**
MOQ: **10 pcs.**



JB 286-00-18 (small)

JB 286-00-36 (large)

Spacer to two types of wall rails.

Item no.: **JB 286-00-18 (small)**
Material: **Stainless steel (AISI 304)**
Colour: **Metal**
Dimensions: **Height / Outer diameter /**
Hole diameter
18 mm. / Ø15 mm. / Ø6,5 mm.
Weight: **20 grams**
MOQ: **40 pcs.**

Item no.: **JB 286-00-36 (large)**
Material: **Stainless steel (AISI 304)**
Colour: **Metal**
Dimensions: **Height / Outer diameter /**
Hole diameter
36 mm. / Ø15 mm. / Ø6,5 mm.
Weight: **40 grams**
MOQ: **20 pcs.**



JB 400-00-18 (small)

JB 400-00-36 (large)

Spacer to two types of wall rails.

Item no.: **JB 400-00-18 (small)**
Material: **Aluminium**
Colour: **Metal**
Dimensions: **Height / Outer diameter /**
Hole diameter
18 mm. / Ø15 mm. / Ø6,5 mm.
Weight: **7 grams**
MOQ: **40 pcs.**

Item no.: **JB 400-00-36 (large)**
Material: **Aluminium**
Colour: **Metal**
Dimensions: **Height / Outer diameter /**
Hole diameter
36 mm. / Ø15 mm. / Ø6,5 mm.
Weight: **14 grams**
MOQ: **20 pcs.**

JBMEDICO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

WALL BRACKET FOR PROTECTION OF MEDICAL EQUIPMENT RAILS AND MORE



JB 180-02-02

**Bracket that protects
medical equipment rails etc.**

It absorbs shocks from beds that accidentally
collide with the rails.

Item no.:	JB 180-02-02
Material:	5 mm. Stainless steel (AISI 304)
Colour:	Metal
Dimensions:	Width x Height x Depth 400 mm. x 3.380 mm. x 1.450 mm.
Weight:	1.000 grams
MOQ:	4 pcs.



JBMEDICO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

MEDICAL LIGHTING FROM GLAMOX LUXO

Examination lamps



Carelite LED

Carelite's LEDs are extremely energy efficient and have a long lifetime of at least 50,000 hours.

The powerful LED module has an excellent light output providing approximately 2,000 lux at 0.6 meters distance.

Item nr.

**CAG026487
CAG026488
CAG026497
CAG026498
CAG026499**



LHH LED

Modern medical examination lamp LED - the light source of the future.

LHH's LEDs are extremely energy efficient and have a long lifetime of at least 50,000 hours. The powerful LED module has an excellent light output, providing more than 24,000 lux at 1 meter.

Item nr.

**LHG026534
LHG026533
LHG026535
LHG026547
LHG026545
LHG026546**

JBMEDICO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

MEDICAL LIGHTING FROM GLAMOX LUXO

Magnifiers



Circus LED

Circus LED is suitable for demanding work in beauty salons, clinics and laboratories, or as a tool for people with impaired vision.

Circus LED has a large lens, aluminium head and a very flexible self-balancing arm and screen mount. The large lens allows comfortable viewing with both eyes. Standard version: 3.5/5 diopter glass lens.

Item nr.

CIL026695
CIL026694
CIL026692
CIL026693



KFM LED

Advanced LED magnifier for many purposes
The KFM LED has a robust metal construction with a slim design, a flexible arm and excellent light.

It is suitable for many industrial applications, clinics and laboratories, and as an aid for people with impaired vision. Standard version comes with a 3 or 5 diopter glass lens. An included fabric bag protects the lens from dust and prevents the lamp from becoming a burning glass via sunlight. Colours: White or light grey.

Item nr.

KFL026034
KFL026035
KFL026036
KFL026037



JBMEDICO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

MEDICAL LIGHTING FROM GLAMOX LUXO

Magnifiers and work lamps



KFM LED ESD

Developed for the electronics industry.

Static discharges can be devastating to electronic components. It is designed to dissipate static discharges in a controlled manner.

Item nr.

KFL026050
KFL026099



LFM LED G2

LFM LED is a lightweight magnifier with a metal construction, LED light sources, and a flexible spring balanced arm in the Luxo tradition.

Voltage from (V): 100, Lumen out: 700.

Item nr.

LFG028194
LFG028195



Wave LED

Wave LED provides a shadow-free magnification with light coming from both sides and which three-dimensional magnification with light from the left or right.

The 3D feature is particularly suitable for working with printed circuit boards and similar tasks.

Item nr.

WAL025948
WAL025949
WAL025950
WAL025951



JBMEDICO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

ELECTRIC PILL CRUSHER V3 - SEVERO



JB 0830-2017

Crushes several standard pills at a time.



Item no.:	JB 0830-2017	Weight:	
Brand:	Severo		
Guarantee:	2 years on the device (1 year on battery - but does not include user errors)		
Cup volume:	30 ml.	Accessories:	Type:
Milling time:	8 seconds for standard pills	JB MED-020	Cups
Volume:	52 dB (low volume)	JB MED-002	Lids
Fully charged:	2 hours	JB MED-003	Rondels
		JB 3000001	Metal stamps
Contacts:	On/Off		
Indicator:	Battery level		
Dimensions:	Height x Width x Depth:		
	240 mm. x 70 mm. x 110 mm.		
	1.5 kg.		

JBMEDICO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

ACCESSORIES FOR SEVERO ELECTRIC PILL CRUSHER V3



JB MED-020

Severo cups 2.100 pcs.

Disposable cup in hard plastic. The bottom has a serrated edge, so the cup can stand firmly in the Severo pellet crusher.

Item no.: **JB MED-020**
Material: **PP**
Dimensions: **Height x Upper diameter**
29 mm. x 47mm.
Volume: **30 ml.**
MOQ: **2.100 pcs.**



JB MED-002

Severo lid 2.000 pcs. for Severo cups.

Item no.: **JB MED-002**
Material: **PE (soft polyethylene)**
Colour: **White**
MOQ: **2.000 pcs.**



JB MED-003

Rondel/Pill crusher 2.100 pcs. pr. box.

Keeps the pill crusher clean and ready to use immediately.

Item no.: **JB MED-003**
Material: **PE**
Diameter: **8 mm.**
Thickness: **50 My**
Certificates: **CE**
MOQ: **2 pcs. boxes with 2.100 pieces in each.**



JB 3000001

Severo, extra metal stamp for tablet crushing.


Item no.: **JB 3000001**
Material: **Stainless steel (AISI 304)**
Colour: **Metal**
Diameter: **Ø55 mm.**
MOQ: **2 pcs.**

JBMEDICO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

MONITOR ARMS AND ACCESSORIES

Adjustable arm with gas spring - Junior Armstrong



JB 21-00-00

Gas arm 2-8 kg. with handle for external locking of height position.

Mounting on wall, table, ceiling, equipment rails, equipment etc. All arms can be potential equalised. Supplied with all necessary accessories.

Item no.: **JB 21-00-00**
Material: **Stainless steel (AISI 304)**
Lifting height: **382 mm.**
Range: **334 mm.**
Weight: **2300 grams**

JB 21-00-15

Gas arm 8-18 kg. with handle for external locking of height position.

Mounting on wall, table, ceiling, equipment rails, equipment etc. All arms can be potential equalised. Supplied with all necessary accessories.

Item no.: **JB 21-00-15**
Material: **Stainless steel (AISI 304)**
Lifting height: **382 mm.**
Range: **334 mm.**
Weight: **2300 grams**




JBMEDICO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

MONITOR ARMS AND ACCESSORIES

Junior Armstrong continued



JB 22-00-00

Manual extension arm with Ø20 mm. shaft, concealed cable routing.

The extension arm can be mounted on e.g. Wall bracket, table-clamp, ceiling-tube, rail clamps and IV-Poles.

Item no.: **JB 22-00-00**
Material: **Stainless steel (AISI 304)**
Total length: **255 mm.**
Height: **150 mm.**
Range: **220 mm.**
Weight: **1.100 grams**

Accessories:

Item no.: **JB 94-00-00**



JB 27-00-00

Monitor bracket with Ø20 mm. shaft.

Can be potential equalised.

Item no.: **JB 27-00-00**
Material: **Stainless steel (AISI 304)**
VESA: **75 mm. x 75 mm.**
VESA: **100 mm. x 100 mm.**
Dimensions: **Depth x Height x Width**
95 mm. x 159 mm. x 114 mm.
Weight: **1000 grams**

Accessories:

Item no.: **JB 10-32-00**



JB 31-00-00

Circular wall bracket fits Ø20 mm. shaft.

Unbrako bolt is standard.

Item no.: **JB 31-00-00**
Material: **Stainless steel (AISI 304)**
Dimensions: **Depth x Height x Width**
58 mm. x 123 mm. x 132 mm.
Weight: **688 grams**

Accessories:

Item no.: **JB 94-00-00**



JB 241-00-00

Extended wall bracket fits Ø20 mm. shaft.

Unbrako bolt is standard.

Item no.: **JB 241-00-00**
Material: **Stainless steel (AISI 304)**
Dimensions: **Depth x Height x Width**
120 mm. x 123 mm. x 130 mm.
Weight: **1100 grams**

Accessories:

Item no.: **JB 94-00-00**

JBMEDICO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

ACCESSORIES FOR MONITOR ARMS



JB 43-01-00

Keyboard Tray with hand support

The plate can be tilted -25°/ +115°

Item no.: **JB 43-01-00**
 Material: **Stainless steel (AISI 304)**
 Dimensions: **Depth (without handle)
 x Height x Width
 277 mm. x 119 mm. x 579 mm.**
 Weight: **3000 grams**

Accessories:
 Item no.: **JB 10-32-00**



JB 209-00-00

Suspension Tray for laptops with Ø20mm shaft.

Item no.: **JB 209-00-00**
 Material: **Stainless steel (AISI 304)**
 Dimensions: **Depth (without handle)
 x Height x Width
 44 mm. x 62 mm. x 135 mm.**
 Weight: **3800 grams**

Accessories:
 Item no.: **JB 10-32-00**



JB 253-00-00

Suspension Tray for medical equipment.

The plate can be tilted -25°/ +115°

Item no.: **JB 253-00-00**
 Material: **Stainless steel (AISI 304)**
 Dimensions: **Depth (without handle)
 x Height x Width
 274 mm x 142 mm. x 370 mm.**
 Weight: **2400 grams**

Accessories:
 Item no.: **JB 10-32-00**



JBMEDICO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

MONITOR ARMS AND ACCESSORIES



JB 90-00-00

Standard clamp desk mount for Ø20 mm. shaft.

Mounting on the edge of the table top.
Table thickness 19-51 mm.
Unbrako bolt is standard.

Item no.: **JB 90-00-00**
Material: **Stainless steel (AISI 304)**
Dimensions: **Depth x Height x Width**
80 mm. x 186,5 mm. x 60 mm.
Weight: **2.500 grams**
Surface: **Sand blasted**
Bushing: **Ø20 mm.**
Cone height: **94 mm.**

Accessories:
Locking handle can be purchased
Item no.: **JB 94-00-00**



JB 91-00-00

Through hole table mount for Ø20 mm. shaft.

Mounting in the table top.
Table thickness 0-41 mm.
Unbrako bolt is standard.

Item no.: **JB 91-00-00**
Material: **Stainless steel (AISI 304)**
Dimensions: **Depth x Height x Width**
80 mm. x 135.5 mm. x 60 mm.
Weight: **2.000 grams**
Surface: **Sand blasted**
Bushing: **Ø20 mm.**
Cone height: **95 mm.**

Accessories:
Locking handle can be purchased
Item no.: **JB 94-00-00**



JB 10-32-00

Locking handle M10x32 mm.

For e.g. screen laptops and keyboard brackets
for tilt/rotate function.

Item no.: **JB 10-32-00**
Material: **Stainless steel (AISI 304)**
MOQ: **2 pcs.**



JB 49-00-00

Bracket for mounting IT equipment on anaesthesia equipment or as wall bracket for Ø20 mm. tube.

Different models for different brands of

Item no.: **JB 49-00-00**
Material: **Stainless steel (AISI 304)**
Dimensions: **Depth x Height x Width**
122 mm. x 162 mm. x 73 mm.
Weight: **1.000 grams**
Surface: **Sand blasted**
Materiale Bushing: **Brass**

Accessories:
Locking handle can be purchased
Item no.: **JB 94-00-00**



JB 94-00-00

Locking handle M6x20 mm.

For e.g. clamping and manual extension arm.

Item no.: **JB 94-00-00**
Material: **Stainless steel (AISI 304)**
MOQ: **5 pcs.**

JBMEDICO

MONITOR ARMS AND ACCESSORIES



JB 62-00-00

Pole clamp mount Ø20/Ø30/Ø20 mm.
for mounting IT equipment.



JB 63-00-00

Pole clamp mount Ø20/Ø20 mm.
for mounting IT equipment.



JB 64-00-00

Pole clamp mount Ø20/Ø20 mm.
for mounting IT equipment.



JB 65-00-00

Pole clamp mount Ø30/Ø20 mm.
for mounting IT equipment.

Item no.: **JB 62-00-00**
Material: **Stainless steel (AISI 304)**
Dimensions: **Depth (without handle) x
Height x Width
44 mm. x 62 mm. x 235 mm.**
Weight: **1400 grams**

Item no.: **JB 63-00-00**
Material: **Stainless steel (AISI 304)**
Dimensions: **Depth (without handle) x
Height x Width
44 mm. x 62 mm. x 135 mm.**
Weight: **800 grams**

Item no.: **JB 63-00-00**
Material: **Stainless steel (AISI 304)**
Dimensions: **Depth (without handle) x
Height x Width
44 mm. x 62 mm. x 135 mm.**
Weight: **800 grams**

Item no.: **JB 65-00-00**
Material: **Stainless steel (AISI 304)**
Dimensions: **Depth (without handle) x
Height x Width
44 mm. x 62 mm. x 235 mm.**
Weight: **1300 grams**

Accessories:
Item no.: **Locking handle available
JB 94-00-00**

Accessories:
Item no.: **Locking handle available
JB 94-00-00**

Accessories:
Item no.: **Locking handle available
JB 94-00-00**

Accessories:
Item no.: **Locking handle available
JB 94-00-00**

JBMEDICO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

ROD SUSPENSION



JB 29-00-04

Solid rod
Ø20 x 600 mm.
2 end plugs.

Item no.: JB 29-00-04
Material: Stainless steel (AISI 304)
Dimensions: Diameter x Length
Ø20 mm. x 600 mm.
Surface: Glassblown
Weight: 1700 grams



JB 29-00-10

Solid rod
Ø20 x 500 mm.
2 end plugs.

Item no.: JB 29-00-10
Material: Stainless steel (AISI 304)
Dimensions: Diameter x Length
Ø20 mm. x 500 mm.
Surface: Glassblown
Weight: 1400 grams



JBMEDICO

SOAP AND ALCOHOL DISPENSERS

For bottles



JB 06-21-30

Wire dispenser with 14 cm. arm
for 1.000 ml. bottles.

Item no.: **JB 06-21-30**
Material: **Stainless steel (AISI 304)**
Surface: **Electropolished**
Dimensions: **Depth (incl. arm operation) x
Height x Width
149 mm. x 175 mm. x 90 mm.**
Weight: **165 grams**
MOQ: **12 pcs.**



JB 06-21-40

Wire dispenser with 10 cm. arm
for 1.000 ml. bottles.

Item no.: **JB 06-21-40**
Material: **Stainless steel (AISI 304)**
Surface: **Electropolished**
Dimensions: **Depth (incl. arm operation) x
Height x Width
145 mm. x 175 mm. x 90 mm.**
Weight: **160 grams**
MOQ: **12 pcs.**



JB 06-21-50

Wire dispenser with 6 cm. arm
for 1.000 ml. bottles.

Item no.: **JB 06-21-50**
Material: **Stainless steel (AISI 304)**
Surface: **Electropolished**
Dimensions: **Depth (incl. arm operation) x
Height x Width
141 mm. x 175 mm. x 90 mm.**
Weight: **150 grams**
MOQ: **12 pcs.**

JBMEDICO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

SOAP AND ALCOHOL DISPENSERS

For bags



JB 42-79-00

Wire dispenser with 14 cm. arm
for 1.000 ml. bags.



JB 42-90-03

Wire dispenser with 10 cm. arm
for 1.000 ml. bags.



JB 42-98-03

Wire dispenser with 6 cm. arm
for 1.000 ml. bags.



JB 728-87-67

Spacer, 85 mm. for all wire dispensers.
Used to mount the dispenser a distance from the
wall, e.g. over a table or sink.

Item no.: **JB 42-79-00**
Material: **Stainless steel (AISI 304)**
Surface: **Electropolished**
Dimensions: **Depth (incl. arm operation) x
Height x Width
165 mm. x 240 mm. x 100 mm.**
Weight: **165 grams**
MOQ: **12 pcs.**

Item no.: **JB 42-90-03**
Material: **Stainless steel (AISI 304)**
Surface: **Electropolished**
Dimensions: **Depth (incl. arm operation) x
Height x Width
161 mm. x 240 mm. x 100 mm.**
Weight: **145 grams**
MOQ: **12 pcs.**

Item no.: **JB 42-98-03**
Material: **Stainless steel (AISI 304)**
Surface: **Electropolished**
Dimensions: **Depth (incl. arm operation) x
Height x Width
157 mm. x 240 mm. x 100 mm.**
Weight: **144 grams**
MOQ: **12 pcs.**

Item no.: **JB 728-87-67**
Material: **Stainless steel (AISI 304)**
Surface: **Electropolished**
Dimensions: **Depth x Height x Width
86 mm. x 58 mm. x 90 mm.**
Weight: **132 grams**
MOQ: **12 pcs.**

JBMEDICO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

DRIP TRAYS FOR SOAP AND ALCOHOL DISPENSERS



JB 102-00-01

**Drip tray to suit everyone
soap and alcohol dispensers.**

For hanging on equipment or directly on the wall.
Protects floor and wall.

Item no.: **JB 102-00-01**
Material: **Stainless steel (AISI 304)**
Surface: **Electropolished**
Dimensions: **Depth x Height x Width**
84,5 mm. x 189,5 mm. x 91,5 mm.
Weight: **209 grams**
MOQ: **12 pcs.**



JB 71-72-51

**Drip tray, drip catcher Ø60 mm. for
mounting on wire soap and alcohol dispensers.**

Protects the floor.

Item no.: **JB 71-72-51**
Material: **Stainless steel (AISI 304)**
Surface: **Electropolished**
Dimensions: **Depth x Height x Width**
80 mm. x 180 mm. x 70 mm.
Weight: **67 grams**
MOQ: **12 pcs.**



JB 213-00-00

**Adapter T-slot fitting for soap and alcohol
detergent Wet Wipe dispenser and drip trays.**

Enables removal and cleaning of dispensers
without the use of too.

Item no.: **JB 213-00-00**
Material: **Stainless steel (AISI 304)**
Surface: **Electropolished**
Dimensions: **Depth x Height x Width**
14 mm. x 50 mm. x 92 mm.
Weight: **174 grams**
MOQ: **10 pcs.**

JBMEDICO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

DISPENSERS WITH DRIP TRAY AND ADAPTER BRACKET

For bottles




JB 50-213-102

Steel dispenser with 6 cm. long arm for 1.000 ml. bottles with drip tray and adapter fitting.

Item no.: JB 50-213-102
Material: Stainless steel (AISI 304)
Colour: Metal
Surface: Polished
Weight: 537 grams
MOQ: 2 pcs.

Accessories:
Item no.: JB 47-00-00

Item no.: JB 120-00-00



JB 40-213-102

Steel dispenser with 10 cm. long arm for 1.000 ml. bottles with drip tray and adapter fitting.

Item no.: JB 40-213-102
Material: Stainless steel (AISI 304)
Colour: Metal
Surface: Polished
Weight: 547 grams
MOQ: 2 pcs.

Accessories:
Item no.: JB 47-00-00

Item no.: JB 120-00-00



JB 30-213-102

Steel dispenser with 14 cm. long arm for 1.000 ml. bottles with drip tray and adapter fitting.

Item no.: JB 30-213-102
Material: Stainless steel (AISI 304)
Colour: Metal
Surface: Polished
Weight: 551 grams
MOQ: 2 pcs.

Accessories:
Item no.: JB 47-00-00

Item no.: JB 120-00-00



JBMEDICO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

DISPENSERS WITH DRIP TRAY AND ADAPTER BRACKET

For bags



JB 98-213-102

Steel dispenser with 6 cm. long arm
for 1.000 ml. bags with drip tray
and adapter fitting.

Item no.: JB 98-213-102
Material: Stainless steel (AISI 304)
Colour: Metal
Surface: Polished
Weight: 537 grams
MOQ: 2 pcs.

Accessories:
Item no.: JB 47-00-00

Item no.: JB 120-00-00

JB 90-213-102

Steel dispenser with 10 cm. long arm
for 1.000 ml. bags with drip tray
and adapter fitting.

Item no.: JB 90-213-102
Material: Stainless steel (AISI 304)
Colour: Metal
Surface: Polished
Weight: 532 grams
MOQ: 2 pcs.

Accessories:
Item no.: JB 47-00-00

Item no.: JB 120-00-00

JB 79-213-102

Steel dispenser with 14 cm. long arm for 1.000 ml.
bags with drip tray and adapter fitting.

Item no.: JB 79-213-102
Material: Stainless steel (AISI 304)
Colour: Metal
Surface: Polished
Weight: 553 grams
MOQ: 2 pcs.

Accessories:
Item no.: JB 47-00-00

Item no.: JB 120-00-00



TABLE AND WALL DISPENSER FOR SQUARE AND ROUND BOTTLES



JB 162-67-87

Table and wall dispenser for round and square 500- 600 ml. bottles

With a handle for T-slot, stainless steel

It consists of 2 pieces:

1. Wire dispensers in 5 mm. wire
2. Table clamp

Fits most square alcohol and round soap bottles. (500 and 600 ml.)
Square bottles 67x67 mm.
Round bottles Ø70 mm.

Item no.:
Material:
Colour:
Surface:
Dimensions:

Weight:
MOQ:

Accessories:

Wall fitting with T-slot
Rail clamp with T-slot (DK)
Rail clamp with T-slot (EU)
Multibracket with T-slot
GCX rail fitting with T-slot

JB 162-67-87
Stainless steel (AISI 304)
Metal
Electropolished
Depth x Height x Width
87 mm. x 162,5 mm. x 77 mm.
100 grams
5 pcs. pr. pcs.

JB 47-00-00
JB 120-00-00
JB 146-00-00
JB 158-00-163
JB 820-625-163





JBMEDICO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

GLOVE, WET WIPE AND SANITARY FACEMASK HOLDER

New HandsUp



JB 800-00-00

With “New HandsUp”, staff can finally see the glove sizes, which are always written on the sides of the boxes.

Fits 99 percent of people glove boxes.

Fits Wet Wipes (the large MAXI)

Fits sanitary towels and mask boxes

“New HandsUp” can be mounted in both horizontal and vertical position.

Mounts to wall fitting, T-slot or wall rail

Autoclavable

Item no.:
Material:

Surface:
Dimensions:

Weight:
MOQ:

Accessories:

Wall bracket with T-slot
Rail clamp with T-slot (DK)
Rail clamp with T-slot (EU)
Multibracket with T-slot
GCX rail bracket with T-slot

JB 800-00-00
Stainless steel (AISI 304) -
4mm. spring steel
Electropolished
Depth x Height x Width
105 mm. x 148 mm. x 260 mm.
362 grams
4 pcs.

JB 47-00-00
JB 120-00-00
JB 146-00-00
JB 158-00-163
JB 820-625-163



JBMEDICO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

WALL DISPENSER FOR WET WIPES



JB 16-05-41-02

Steel dispenser, Wet Wipe MAXI

Stainless steel wall holder for wall mounting of Wet Wipe products.

Item no.:
Material:
Brand:
Colour:
Surface:
Weight:
MOQ:

JB 16-05-41-02
Stainless steel (AISI 316)
Wet Wipe
Metal
Electropolished
230 grams
10 pcs.

Accessories:
JB 213-00-00



JB 16-05-40-02

Steel dispenser, Wet Wipe MINI

Stainless steel wall holder for wall mounting of Wet Wipe products.

Item no.:
Material:
Brand:
Colour:
Surface:
Weight:
MOQ:

JB 16-05-40-02
Stainless steel (AISI 316)
Wet Wipe
Metal
Electropolished
180 grams
10 pcs.

Accessories:
JB 213-00-00

JBMEDICO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

ACCESSORIES FOR GLOVE BOX HOLDERS



JB 47-00-00

T-slot wall bracket for 1 pcs. glove box holder.

C-profile with 2 pcs. countersunk mounting holes and stop screw.

Item no.: **JB 47-00-00**
 Material: **Aluminium**
 Surface: **Anodized**
 Dimensions: **Depth x Height x Width**
12,5 mm. x 60 mm. x 36 mm.
 Hole Size: **Ø 6,4 mm.**
 Weight: **36 grams**
 MOQ: **10 pcs.**



JB 47-00-02

T-slot wall bracket for 2 pcs. glove box holder.

C-profile with 4 pcs. countersunk mounting holes and stop screw.

Item no.: **JB 47-00-02**
 Material: **Aluminium**
 Surface: **Anodized**
 Dimensions: **Depth x Height x Width**
12,5 mm. x 210 mm. x 36 mm.
 Hole Size: **Ø 6,4 mm.**
 Weight: **129 grams**
 MOQ: **5 pcs.**



JB 47-00-03

T-slot wall bracket for 3 pcs. glove box holder.

C-profile with 6 pcs. countersunk mounting holes and stop screw.

Item no.: **JB 47-00-03**
 Material: **Aluminium**
 Surface: **Anodized**
 Dimensions: **Depth x Height x Width**
12,5 mm. x 360 mm. x 36 mm.
 Hole Size: **Ø 6,4 mm.**
 Weight: **222 grams**
 MOQ: **5 pcs.**



JB 47-00-04

T-slot wall bracket for 4 pcs. glove box holder.

C-profile with 6 pcs. countersunk mounting holes and stop screw.

Item no.: **JB 47-00-04**
 Material: **Aluminium**
 Surface: **Anodized**
 Dimensions: **Depth x Height x Width**
12,5 mm. x 520 mm. x 36 mm.
 Hole Size: **Ø 6,4 mm.**
 Weight: **320 grams**
 MOQ: **5 pcs.**

JBMEDICO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

ALUMINIUM SLIDING BRACKET FOR T-SLOT RAILS



JB 90-92-05-30

Aluminium sliding bracket, lockable

Hole patterns in plate part fit various accessories and dispensers

Item no.:
Material:
Colour:
Dimensions:

Weight:
MOQ:

JB 90-92-05-30
Aluminium
Metal
Depth x Height x Width
87 mm. x 162,5 mm. x 77 mm.
150 grams
4 pcs.


Accessories:

Steel dispenser, Wet Wipe MAXI
Steel dispenser, Wet Wipe MINI
T-slot bracket length 210 mm.
T-slot bracket length 360 mm.
T-slot bracket length 520 mm.
CombiPlum wire dispenser 14 cm.
CombiPlum wire dispenser 10 cm.
CombiPlum thread dispenser 6 cm.
Dispenser for 1.000 ml. bottles 14 cm.
Dispenser for 1.000 ml. bottles 10 cm.
Dispenser for 1.000 ml. bottles 6 cm.
Drip tray for wire dispensers
Fixing device Ø25 bushing hole
Fixing device Ø18 socket hole

JB 16-05-41-02
JB 16-05-40-02
JB 47-00-02
JB 47-00-03
JB 47-00-04
JB 42-79-00
JB 42-90-03
JB 42-98-03
JB 06-21-50
JB 06-21-40
JB 06-21-30
JB 102-00-01
JB 205-00-00
JB 205-00-18



Accessories proposal No. 1



Accessories proposal No. 2



Accessories proposal No. 3



Accessories proposal No. 4

JBMEDICO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

GLOVE BOX HOLDER



JB 110-00-00

Glove holder for 3 glove boxes wall mounted.

Suitable for glove boxes, face wipes, aprons and handkerchiefs, etc.

Item no.: JB 110-00-00
Material: Stainless steel (AISI 304)
Dimensions: Depth x Height x Width (outside)
93 mm. x 406 mm. x 285 mm.
Opening: 183 mm.
Weight: 1854 grams
MOQ: 4 pcs.



JB 110-00-04

Glove holder for 4 glove boxes wall mounted.

Suitable for glove boxes, face wipes, aprons and handkerchiefs, etc.

Item no.: JB 110-00-04
Material: Stainless steel (AISI 304)
Dimensions: Depth x Height x Width (outside)
93 mm. x 520 mm. x 285 mm.
Opening: 183 mm.
Weight: 2346 grams
MOQ: 4 pcs.

JBMEDICO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

IV STANDS / IV POLES



JB 306-2-317-211

**"One Hand" "Small" straight legs, large wheels,
long straight rod, 2 strong angled hooks**

The flexible IV stand can be modified as needed. The solid and rust-proof specially designed base allows the attachment of various items in a safe way without the risk of the stand tipping over.

Item no.:

Material / Surface
Color of Plastic Parts
Hooks
Load capacity pr. hook
Weight without equipment
Max. load capacity
Height with/without wheels (telescopic max.)
Height with/without wheels (telescopic min.)
Outer tube diameter and length
Inner tube diameter and length
Setting
Height adjustment
Base diameter, straight leg
Wheel diameter according to DIN EN 60601

JB 306-2-317-211

**Stainless steel (AISI 304) / Electropolished
Black
2 strong angled hooks
5 kg.
15 kg.
40 kg.
2.369 mm. / 2.232 mm.
1.605 mm. / 1.470 mm.
Ø32/1.281 mm.
Ø20 x 940 mm.
"One Hand" adjustable
765 mm.
515 mm.
Ø100 mm.
1 x antistatic wheels
2 x lockable wheels
2 x standard wheels**

Safety distance between both tubes when outer tube reaches its lowest height.

Integrated stop prevents separation of the adjustable tubes.

**Stainless steel ball bearing: Yes
An integrated base-weight facilitates a safe
of various pumps and other medical
equipment.
CE Approved**



Security

- Columns are equipped with plastic stops (no risk of damage)
- Sufficient safety distance in the collapsed state as a precaution against risk of injury (e.g. pinching of fingers)
- Increased stability due to conical connection between stand and foot section
- Safe control of the lowering of the inner tube
- No risk of injury due to rounded design of hooks
- Due to the geometry/angle of the hooks, the infusion bottles cannot fall out
- Flexible handle for drip stand can be adapted to the individual user.

Technology

- All plastic parts are made of lightweight materials
- Pipes in Stainless steel according to DIN EN 10217-7
- Material ultrasonic and eddy current tested (according to EN 10246-9 /3).
- Integrated pull-out limiter prevents unintentional separation of the adjustable tubes
- Safety bottle cross according to DIN ISO 15375
- Good running characteristics through high quality 100 mm. mono Rolls, 1 Rolls are antistatic apply (One Hand stands).
- The introduction of hygienically sensitive liquids is prevented by the sealing lid.
- Reinforced bottle cross allows connection of 5 liters bottles.

JB MEDICO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

IV POLE XL

One Hand - Four Hooks



JB 332-2-320-190

Drop stand "One Hand" with curved legs, large wheels, curved long rod, 4 strong angled hooks.

The flexible IV stand can be modified as needed. The solid and rust-proof specially designed base allows the attachment of various items in a safe way without the risk of the stand tipping over.

Item no.:

Material / Surface
Color of Plastic Parts
Hooks
Load capacity pr. hook
Weight without equipment
Max. load capacity
Height with/without wheels (telescopic max.)
Height with/without wheels (telescopic min.)
Outer tube diameter
Inner tube diameter
Setting
Height adjustment
Base diameter, straight leg
Wheel diameter according to DIN EN 60601

JB 332-2-320-190

Electropolished Stainless steel (AISI 304)
Black
Four strong hooks
5 kg.
15 kg.
50 kg.
2.281 mm. / 2.144 mm.
1.516 mm. / 1.379 mm.
Ø32 mm.
Ø20 mm.
"One Hand" adjustable
765 mm.
545 mm.
Ø100 mm.
1 x antistatic wheels
2 x lockable wheels
2 x standard wheels

Safety distance between both tubes when outer tube reaches its lowest height

Integrated stop prevents separation of the adjustable tubes.

Stainless steel ball bearing: Yes
An integrated base-weight facilitates a secure attachment of various pumps and other medical equipment.
CE Approved



Security

- Pillars are equipped with plastic stops (no risk of injury)
- Sufficient safety clearance in collapsed condition as a precaution against risk of injury (e.g. pinching of fingers)
- Increased stability due to conical connection between stand and foot section
- Safe control of the lowering of the inner tube
- No risk of injury due to rounded design of hooks
- Due to the geometry/angle of the hooks, the infusion bottles cannot fall out
- Flexible handle for drip stand can be adapted to the individual user.

Technology

- All plastic parts are made of lightweight materials
- Pipes in Stainless steel according to DIN EN 10217-7
- Material ultrasonic and eddy current tested (according to EN 10246-9 /3).
- Integrated pull-out limiter prevents unintentional separation of the adjustable tubes
- Safety bottle cross according to DIN ISO 15375
- Good running characteristics through high quality 100 mm. mono Rolls, 1 Rolls are antistatic apply (One Hand stands).
- The introduction of hygienically sensitive liquids is prevented by the sealing lid.
- Reinforced bottle cross allows connection of 5 liters bottles.

JB MEDICO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

IV POLE XXL

One Hand - Six Hooks



JB 332-2-321-216

Drop stand "One Hand" with curved legs, large wheels, curved long rod, 6 strong straight hooks

The flexible IV stand can be modified as needed. The solid and stainless steel frame, allows the attachment of various items in a safe way without the risk of the stand tipping over.

Item no.:

Material / Surface
Color of Plastic Parts
Hooks
Load pr. hook
Weight without equipment
Max. load
Height with/without wheels (telescopic max.)
Height with/without wheels (telescopic min.)
Outer tube diameter
Inner tube diameter
Setting
Height adjustment
Base diameter, straight leg
Wheel diameter according to DIN EN 60601

JB 332-2-321-216

Electropolished Stainless steel (AISI 304)
Black
Six strong straight hooks
5 kg.
15 kg.
50 kg.
2.281 mm. / 2.144 mm.
1.516 mm. / 1.379 mm.
Ø32 mm.
Ø20 mm.
"One Hand" adjustable
765 mm.
545 mm.
Ø100 mm.
1 x antistatic wheels
2 x lockable wheels
2 x standard wheels

Safety distance between both tubes when outer tube reaches its lowest height.

Integrated stop prevents separation of the adjustable tubes.

Stainless steel ball bearing: Yes
An integrated base-weight facilitates a secure attachment of various pumps and other medical equipment
CE Approved



Security

- Columns are equipped with plastic stops (no risk of damage)
- Sufficient safety distance in the collapsed state as a precaution against risk of injury (e.g. pinching of fingers)
- Increased stability due to conical connection between stand and foot section
- Safe control of the lowering of the inner tube
- No risk of injury due to rounded design of hooks
- Due to the geometry/angle of the hooks, the infusion bottles cannot fall out
- Flexible handle for drip stand can be adapted to the individual user.

Technology

- All plastic parts are made of lightweight materials
- Pipes in Stainless steel according to DIN EN 10217-7
- Material ultrasonic and eddy current tested (according to EN 10246-9 /3).
- Integrated pull-out limiter prevents unintentional separation of the adjustable tubes
- Safety bottle cross according to DIN ISO 15375
- Good running characteristics through high quality 100 mm. mono Rolls, 1 Rolls are antistatic apply (One Hand stands).
- The introduction of hygienically sensitive liquids is prevented by the sealing lid.
- Reinforced bottle cross allows connection of 5 liters bottles.

JB MEDICO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

IV STAND / IV POLE



JB 306-2-317-190

Drop stand "One Hand" with straight legs, large wheels, straight long rod, 4 strong angled hooks.

The flexible IV stand can be modified as needed. The solid and stainless steel frame, allows the attachment of various items in a safe way without the risk of the stand tipping over.

Item no.:

Material / Surface
Color of Plastic Parts
Hooks
Load pr. hook
Weight without equipment
Max. load
Height with/without wheels (telescopic max.)
Height with/without wheels (telescopic min.)
Outer tube diameter
Inner tube diameter
Setting
Height adjustment
Base diameter, straight leg
Wheel diameter according to DIN EN 60601

JB 306-2-317-190

Electropolished Stainless steel (AISI 304)
Black
Four strong straight hooks
5 kg.
15 kg.
40 kg.
2.369 mm. / 2.232 mm.
1.605 mm. / 1.470 mm.
Ø32/1.281 mm.
Ø20 x 940 mm.
"One Hand" adjustable
765 mm.
515 mm.
Ø100 mm.
1 x antistatic wheels
2 x lockable wheels
2 x standard wheels

Safety distance between both tubes when outer tube reaches its lowest height.

Integrated stop prevents separation of the adjustable tubes.

Stainless steel ball bearing: Yes
An integrated base-weight facilitates a secure attachment of various pumps and other medical equipment
CE Approved



Security

- Columns are equipped with plastic stops (no risk of damage)
- Sufficient safety distance in the collapsed state as a precaution against risk of injury (e.g. pinching of fingers)
- Increased stability due to conical connection between stand and foot section
- Safe control of the lowering of the inner tube
- No risk of injury due to rounded design of hooks
- Due to the geometry/angle of the hooks, the infusion bottles cannot fall out
- Flexible handle for drip stand can be adapted to the individual user.


Technology

- All plastic parts are made of lightweight materials
- Pipes in Stainless steel according to DIN EN 10217-7
- Material ultrasonic and eddy current tested (according to EN 10246-9 /3).
- Integrated pull-out limiter prevents unintentional separation of the adjustable tubes
- Safety bottle cross according to DIN ISO 15375
- Good running characteristics through high quality 100 mm. mono Rolls, 1 Rolls are antistatic apply (One Hand stands).
- The introduction of hygienically sensitive liquids is prevented by the sealing lid.
- Reinforced bottle cross allows connection of 5 liters bottles.

JB MEDICO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361


IV STAND / INFUSION POLE - SUPPORT COLUMNS



JB 320-00-02

Rod long, curved.


Item no.: **JB 320-00-02**
 Dimensions: **Length x Depth:**
 1.281 mm. x 32 mm.
 Weight: **2000 grams**



JB 317-00-06

Rod long, straight.

Item no.: **JB 317-00-06**
 Dimensions: **Length x Width**
 1.281 mm. x Ø32 mm.
 Weight: **1800 grams**



JB 317-00-13

Telescopic pole
 "One Hand".

Item no.: **JB 317-00-13**
 Dimensions: **Length x Depth:**
 940 mm. x Ø20 mm.
 Weight: **1000 grams**



INFUSION AND IV POLE - ACCESSORIES

Hooks



JB 204-00-00

Hook top - straight 2 strong hooks including bolt.

Item no.: **JB 204-00-00**
Material: **Stainless steel (AISI 304)**
Surface: **Electropolished**
Dimensions: **Height x Width**
406 mm. x 285 mm.
Hook distance: **252 mm.**
Load per hook: **5 kg.**
Weight: **248 grams**



JB 211-00-00

Hook top - angled 2 strong hooks including bolt.

Item no.: **JB 211-00-00**
Material: **Stainless steel (AISI 304)**
Surface: **Electropolished**
Dimensions: **Height x Width**
406 mm. x 285 mm.
Hook distance: **204 mm.**
Load per hook: **5 kg.**
Weight: **248 grams**



JB 212-00-00

Hook top - angled 3 strong hooks including bolt.

Item no.: **JB 212-00-00**
Material: **Stainless steel (AISI 304)**
Surface: **Electropolished**
Dimensions: **Height x Width**
84 mm. x 234 mm.
Hook distance: **120 mm.**
Load per hook: **5 kg.**
Weight: **363 grams**



JB 190-00-08

Hook top - 4 strong hooks including bolt.

Item no.: **JB 190-00-08**
Material: **Stainless steel (AISI 304)**
Surface: **Electropolished**
Dimensions: **Height x Width**
406 mm. x 285 mm.
Hook distance: **126 mm.**
Load per hook: **5 kg.**
Weight: **448 grams**

JBMEDICO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

IV STAND / IV POLE - ACCESSORIES

Hooks



JB 216-00-00

Hook top - 6 strong hooks including bolt.


Item no.: **JB 216-00-00**
 Material: **Stainless steel (AISI 304)**
 Surface: **Electropolished**
 Dimensions: **Height x Width**
406 mm. x 285 mm.
 Hook distance: **120 mm.**
 Load per hook: **5 kg.**
 Weight: **596 grams**



JB 170-02-10

IV 2-Hooks for infusion bags, AISI 304


Item no.: **JB 170-02-10**
 Material: **10 mm. stainless steel (AISI 304)**
 Surface: **Electrically polished**
 Dimensioner: **Height x Width**
105mm x 175 mm
 Total load capacity: **4 kg.**
 Distance between hooks: **110 mm.**
 Weight: **229 grams**
 MOQ: **2 pcs.**



JB 280-03-10

IV 3-Hooks for infusion bags, AISI 304

Item no.: **JB 280-03-10**
 Material: **10 mm. stainless steel (AISI 304)**
 Surface: **Electrically polished**
 Dimensioner: **Height x Width**
44,88 mm. x 280 mm.
 Total load capacity: **6 kg.**
 Distance between hooks: **110 mm.**
 Weight: **332 grams**
 MOQ: **2 pcs.**



JB 390-04-10

IV 4-Hooks for infusion bags, AISI 304

Item no.: **JB 390-04-10**
 Material: **10 mm. stainless steel (AISI 304)**
 Surface: **Electrically polished**
 Dimensioner: **Height x Width**
105mm x 390 mm
 Total load capacity: **8 kg.**
 Distance between hooks: **110 mm.**
 Weight: **433 grams**
 MOQ: **2 pcs.**

JBMEDICO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

IV STAND / IV POLE - ACCESSORIES

Wheels and bases



5387PJC100P30-13

Swivel castor, 100 mm. total brake.
Stainless (Levina), Grey.

Item no.:	5387PJC100P30-13
Wheel Ø	Ø13 mm.
Centre hole Ø:	Ø13 mm.
Overall height:	137 mm.
Load:	100 kg. (dynamic)
Standard:	EN 12530
Weight:	430 grams
Stainless steel ball bearing:	Yes
Electrical conductivity:	No
MOQ:	18 pcs.



5380PJC100P30-13

Swivel castor without brake 100 mm.
Stainless (Levina), Grey.

Item no.:	5380PJC100P30-13
Wheel Ø	Ø13 mm.
Centre hole Ø:	Ø13 mm.
Overall height:	137 mm.
Load:	100 kg (dynamisk)
Standard:	EN 12530
Weight:	430 grams
Stainless steel ball bearing:	Yes
Electrical conductivity:	No
MOQ:	20 pcs.



5380XSC100P30-13

Swivel castor without brake 100 mm.
Stainless (Levina) Grey, Conductive.

Item no.:	5380XSC100P30-13
Wheel Ø	Ø13 mm.
Centre hole Ø:	Ø13 mm.
Overall height:	137 mm.
Load:	100 kg (dynamisk)
Standard:	EN 12530
Weight:	430 grams
Stainless steel ball bearing:	Yes
Electrical conductivity:	Yes
MOQ:	20 pcs.



JB 332-00-00

Main base, bent legs for drop stand.
(Tente wheels are included)

Item no.:	JB 332-00-00
Base diameter:	Ø545 mm.
Leg Profile:	Ø25 mm.
Weight:	12,3 kg.



JB 306-00-00

Main base, straight leg for drop stand.
(Tente wheels are included)

Item no.:	JB 306-00-00
Base diameter:	Ø545 mm.
Leg Profile:	Ø25 mm.
Weight:	11,6 kg.

JBMEDICO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

IV STAND / IV POLES - ACCESSORIES



JB 294-00-05

Stainless steel handle with an ergonomic grip.

The handle is easy to mount and adjustable, can easily be set in relation to other equipment/user's height.

Item no.: **JB 294-00-05**
 Material: **Stainless steel (AISI 304)**
 Surface: **Electropolished**
 Dimensions: **Height x Width x Length**
36 mm. x 194 mm. x 280 mm.
 Handle diameter: **Ø20 mm.**
 Weight: **733 grams**

Accessories: **JB 163-00-01**
 Dimensions: **H: 60 mm. x W: 36,2 mm.**
x D: 12,5 mm.
 Weight: **36 grams**



JB 190-00-05

Telescopic rod 940 mm. on rail or

tube- mounted.

Several combinations (Picture is with hook top - to

Item no.: **JB 190-00-05**
 Material: **Stainless steel (AISI 304)**
 Surface: **Electropolished**
 Length: **940 mm.**
 Diameter: **Ø20 mm.**
 Hook distance: **204 mm.**
 Weight: **795 grams**



JB 289-00-00

Transducer holder for 4 pcs. A needles.

A-cannulas made of aluminium.
 Can be washed, fits most transducers.

Item no.: **JB 289-00-00**
 Material: **Stainless steel (AISI 304)**
 Surface: **Electropolished**
 Dimensions: **Height x Width x Depth**
90 mm. x 185 mm. x 19,5 mm.
241 grams
 Weight: **241 grams**
 MOQ: **2 pcs.**

Accessories: **JB 158-00-163**
(Multibracket med T-spors beslag)



JB 205-08-05

Drip pole for infusion bags for IV fluids.


Fits all our products with T-slot handle.

Item no.: **JB 205-08-05**
 Material: **Stainless steel (AISI 304)**
 Surface: **Electropolished**
 Dimensions: **Height x Width**
50 mm. x 205 mm.
5-6 drip bags
 Capacity: **5-6 drip bags**
 Weight: **135 grams**
 MOQ: **4 pcs.**

JBMEDICO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

IV STAND / IV POLE - ACCESSORIES



JB 25-38-BB

Lock for B. Braun Pump Rack

Securing against rotation of the "Docking Station" on the stand so that the weight of the pumps does not cause the Infusion Stand to be unstable. The lock consists of a front plate that slides into the "rack" and a rear plate that is clamped around the infusion stand with four bolts and four nuts.

+ an extra bolt to prevent unwanted removal without the use of tools

Item no.:
Material:
Surface:
Dimensions:

Weight:
MOQ:

Note:

JB 25-38-BB
Stainless steel (AISI 304)
Glassblown
Height x Width
130 mm. x 75 mm.
387 grams
2 pcs.

Suitable for pipe diameter
from 25-38 mm.



JBMEDICO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

WALL BRACKET FOR CO2 BOTTLES, Ø183 MM. HOLE



JB 283-01-00

Wall bracket for CO2 bottles fits for example Air Liquide 6 liters.

Bottle holder for wall for up to Ø183 mm. (Inside dimensions)

Item no.:

Material:

Surface:

Dimensions:

Weight:

Note:

JB 283-01-00

Stainless steel (AISI 304) - 4 mm.

Electropolished

Height x Width

306 mm. x 215 mm.

2.340 grams

Fits up to Ø183 mm. bottles



JBMEDICO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

OXYGEN AND GAS CYLINDERS HOLDERS



Picture is with our multi fitting

JB 290-01-00

Oxygen and gas cylinder holder top
for Ø117 mm. bottles.

Can be ordered with a bottom ring JB 290-02-00

Item no.: **JB 290-01-00**
Material: **Ø8 mm. SS, POM,
Nylon(PA6)**
Width: **175 mm. Ext. w. handle**
Weight: **255 grams**
MOQ: **2 pcs.**

Accessories: **Multibracket**
Item no.: **JB 158-00-163**



Picture is with our multi fitting

JB 290-02-00

Oxygen and gas cylinder holder-bottom
for Ø117 mm. cylinders.

Can be ordered with a bottom ring JB 290-02-00

Item no.: **JB 290-02-00**
Material: **Ø8 mm. SS, POM,
Nylon(PA6)**
Width: **175 mm. Ext. w. handle**
Height: **112 mm.**
Weight: **315 grams**
MOQ: **2 pcs.**

Accessories: **Multibracket**
Item no.: **JB 158-00-163**



Picture is with our multi fitting

JB 277-01-00

Oxygen and gas cylinder holder top
for Ø140 mm. cylinders.

Can be ordered with a bottom ring JB 290-02-00

Item no.: **JB 277-01-00**
Material: **Ø8 mm. SS, POM,
Nylon(PA6)**
Width: **189 mm. Ext. w. handle**
Weight: **278 grams**
MOQ: **2 pcs.**

Accessories: **Multibracket**
Item no.: **JB 158-00-163**



Picture is with our multi fitting

JB 277-02-00

Oxygen and gas cylinder holder-bottom
for Ø140 mm. cylinders.

Can be ordered with a bottom ring JB 290-02-00

Item no.: **JB 277-02-00**
Material: **Ø8 mm. SS, POM,
Nylon(PA6)**
Width: **198 mm. Ext. w. handle**
Height: **112 mm.**
Weight: **341 grams**
MOQ: **2 pcs.**

Accessories: **Multibracket**
Item no.: **JB 158-00-163**

JBMEDICO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

SHARPS CONTAINER HOLDER - ROUND

e.g. "Grathwol"




Sharps Container 0,7 and 1,5 liters

Sharps container 0,7 and 1,5 liters
AP Medical DISPO Line.

Yellow or red lid.

Accessories:

Wall bracket with T-slot
Rail clamp with T-slot (DK)
Rail clamp with T-slot (EU)
Multibracket with T-slot
GCX rail bracket with T-slot



JB 255-00-00

Sharps container holder til e.g.
AP Medical 0,7 and 1,5 liters Sharps Containers.

Item no.:	JB 255-00-00
Inner diameter:	Ca. Ø100 mm.
Material:	Stainless steel (AISI 304)
Surface:	Electropolished
Height:	50 mm.
Depth:	143 mm. (outer dimension)
Thread diameter:	8 mm.
Weight:	189 grams
MOQ:	5 pcs.

JBMEDICO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

SHARPS CONTAINER HOLDER - ROUND

e.g. "AP MEDICAL"



Sharps Container 2 liters

Sharps container 2 liters
AP Medical DISPO Line.

Yellow or red lid.



JB 149-00-00

Sharps container holder for
1.0 - 1,3 - 1,5 - 2 and 3 liters Sharps Containers.

Fits perfectly to Klinion EasyCare "Mediq /
UDSON" 1,3 and 2 liters AP Medical DISPO Line
boxes.

Item no.:	JB 149-00-00
Inner diameter:	Ø 134 mm.
Material:	Stainless steel (AISI 304)
Surface:	Electropolished
Height:	50 mm.
Depth:	172 mm. (outer dimension)
Thread diameter:	8 mm.
Weight:	220 grams
MOQ:	5 pcs.

SHARPS CONTAINER HOLDER - OVAL

e.g. "AP MEDICAL"



Sharps Container 2, 4 and 7 liters

Sharps container 2 liters
AP Medical DISPO Line.

Yellow or red lid.



JB 159-00-00

Sharps container holder for e.g. AP Medical 2, 4
and 7 liters oval boxes.

Fits perfectly to Klinion EasyCare "Mediq /UD-
SON" AP Medical DISPO Line boxes.
Yellow or red lid.

Item no.:	JB 159-00-00
Inner diameter:	133 mm. x 133 mm.
Material:	Stainless steel (AISI 304)
Surface:	Electropolished
Height:	50 mm.
Depth:	175 mm. (outer dimension)
Thread diameter:	8 mm.
Weight:	248 grams
MOQ:	5 pcs.

Accessories:
Wall bracket with T-slot
Rail clamp with T-slot (DK)
Rail clamp with T-slot (EU)
Multibracket with T-slot
GCX rail bracket with T-slot

JBMEDICO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

SHARPS CONTAINER HOLDER - OVAL

For e.g. "AP Medical"



Sharps Container 3 liters

Sharps container 3 liters
AP Medical DISPO Line.

Yellow or red lid.

JB 169-00-00

Sharps container holder for e.g. AP Medical
3 liters, oval container.

Accessories:

Wall bracket with T-slot
Rail clamp with T-slot (DK)
Rail clamp with T-slot (EU)
Multibracket with T-slot
GCX rail bracket with T-slot

Item no.:	JB 169-00-00
Inner diameter:	143 mm. x 143 mm.
Material:	Stainless steel (AISI 304)
Surface:	Electropolished
Height:	50 mm.
Depth:	175 mm. (outer dimension)
Thread diameter:	8 mm.
Weight:	319 grams
MOQ:	5 pcs.

SHARPS CONTAINER HOLDER - SQUARED

For e.g. "Sharpsafe"



Sharps Container 4 and 7 liters

Sharps containers 4 / 7 liters
"Sharpsafe" container.

Yellow or red lid

JB 154-00-00

Sharps container holder for e.g. "Sharpsafe"
4 and 7 liters squared container.

Item no.:	JB 154-00-00
Inner diameter:	156 mm. x 156 mm.
Material:	Stainless steel (AISI 304)
Surface:	Electropolished
Height:	50 mm.
Depth:	208 mm. (outer dimension)
Thread diameter:	8 mm.
Weight:	280 grams
MOQ:	5 pcs.



JBMEDICO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

SHARPS CONTAINER HOLDER - SQUARED

For e.g. "Sharpsafe"



JB 166-00-00

Sharps container holder for e.g.
"Sharpsafe Quiver"

7,5 liters, squared container (long model)



Sharps Container 7,5 Liters

Sharps container "Sharpsafe"
7.5 liters QUIVER.

Yellow or red lid.

Item no.: **JB 166-00-00**
Inner diameter: **101 mm. x 147 mm.**
Material: **Stainless steel (AISI 304)**
Surface: **Electropolished**
Height: **50 mm.**
Depth: **153 mm. (outer dimension)**
Thread diameter: **8 mm.**
Weight: **247 grams**
MOQ: **5 pcs.**

Accessories:
Wall bracket with T-slot
Rail clamp with T-slot (DK)
Rail clamp with T-slot (EU)
Multibracket with T-slot
GCX rail bracket with T-spur

SHARPS CONTAINER HOLDER - SQUARED

For e.g. "Europlásticos"



JB 170-00-00

Sharps container holder for e.g. "Sharpsafe"
1 liter squared container.



Sharps Container 1 Liters

Sharps container for 1 liter boks
from Europlásticos.

Yellow or red lid.

Item no.: **JB 170-00-00**
Inner diameter: **96 mm. x 96 mm.**
Material: **Stainless steel (AISI 304)**
Surface: **Electropolished**
Height: **50 mm.**
Depth: **208 mm. (outer dimension)**
Thread diameter: **8 mm.**
Weight: **200 grams**
MOQ: **5 pcs.**



JBMEDICO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

SHARPS CONTAINER HOLDER - SQUARED

Other Assortment - e.g. "Safebox"



Sharps Container 4L, 6L and 11 Liters

Sharps Container "Safebox" 4, 6 and 11 liters

Yellow or red lid.

Accessories:

Wall bracket with T-slot
Rail clamp with T-slot (DK)
Rail clamp with T-slot (EU)
Multibracket with T-slot
GCX rail bracket with T-slot



JB 166-00-00

Sharps Container Holder for e.g.
"Sharsafe Quiver"

7,5 liters, squared container (long model)

Item no.: **JB 166-00-00**
Inner diameter: **101 mm. x 147 mm.**
Material: **Stainless steel (AISI 304)**
Surface: **Electropolished**
Height: **50 mm.**
Depth: **153 mm. (outer dimension)**
Thread diameter: **8 mm.**
Weight: **247 grams**
MOQ: **5 pcs.**

SHARPS CONTAINER HOLDER - ROUND

Other Assortment



JB 150-00-00

Sharps Container Holder for
0,5- and 1 liter round container.

Fits e.g. Dispo 0,5 liter box from Grathwol.

Item no.: **JB 150-00-00**
Inner diameter: **Ø73 mm.**
Material: **Stainless steel (AISI 304)**
Surface: **Electropolished**
Height: **50 mm.**
Depth: **112 mm. (outer dimension)**
Thread diameter: **8 mm.**
Weight: **150 grams**
MOQ: **5 pcs.**

JB 151-00-00

Sharps Container Holder for
0.6 liter - round container.

Fits e.g. Multisafe 1 liter and USON 500 ml.

Item no.: **JB 151-00-00**
Inner diameter: **Ø89 mm.**
Material: **Stainless steel (AISI 304)**
Surface: **Electropolished**
Height: **50 mm.**
Depth: **134 mm. (outer dimension)**
Thread diameter: **8 mm.**
Weight: **171 grams**
MOQ: **5 pcs.**

JBMEDICO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

SHARPS CONTAINER HOLDER - ROUND

Other Assortment



JB 264-00-00

Sharps Container Holder for e.g. Mediq /USON
5 and 7 liters round container.

Item no.: **JB 264-00-00**
Inner diameter: **Ø200 mm.**
Material: **Stainless steel (AISI 304)**
Surface: **Electropolished**
Height: **50 mm.**
Depth: **244 mm. (outer dimension)**
Thread diameter: **8 mm.**
Weight: **286 grams**
MOQ: **5 pcs.**



JB 266-00-00

Sharps Container Holder for e.g. USON 11 liters
and Jokey container from Mediq.

Item no.: **JB 266-00-00**
Inner diameter: **Ø251 mm.**
Material: **Stainless steel (AISI 304)**
Surface: **Electropolished**
Height: **50 mm.**
Depth: **295 mm. (outer dimension)**
Thread diameter: **8 mm.**
Weight: **345 grams**
MOQ: **5 pcs.**



JB 267-00-00

Sharps Container Holder for e.g. USON 21 liters
and Jokey container from Mediq.

Item no.: **JB 267-00-00**
Inner diameter: **Ø280 mm.**
Material: **Stainless steel (AISI 304)**
Surface: **Electropolished**
Height: **50 mm.**
Depth: **324 mm. (outer dimension)**
Thread diameter: **8 mm.**
Weight: **368 grams**
MOQ: **5 pcs.**



Accessories:
Wall bracket with T-slot
Rail clamp with T-slot (DK)
Rail clamp with T-slot (EU)
Multibracket with T-slot
GCX rail bracket with T-slot

JBMEDICO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

WIRE BASKETS AND WIRE PRODUCTS

Stainless steel



JB 161-150-100

Wire basket 1 liters, conical with hook for T-slot

Item no. **JB 161-150-100**
Dimensions: **Width x Depth x Height**
150mm. x 100mm. x 100mm.
Material: **Stainless steel (AISI 304)**
Surface: **Electropolished**
Weight: **325 grams**
Volume: **1 liters**
MOQ: **5 pcs.**



JB 161-235-163

Wire basket 3,5 liters, conical with hook for T-slot

Item no. **JB 161-235-163**
Dimensions: **Width x Depth x Height**
235mm. x 163mm. x 146mm.
Material: **Stainless steel (AISI 304)**
Surface: **Electropolished**
Weight: **800 grams**
Volume: **3,5 liters**
MOQ: **5 pcs.**



JB 161-350-150

Wire basket 6,9 liters, conical with hook for T-slot

Item no. **JB 161-350-150**
Dimensions: **Width x Depth x Height**
350mm. x 150mm. x 200mm.
Material: **Stainless steel (AISI 304)**
Surface: **Electropolished**
Weight: **800 grams**
Volume: **6,9 liters**
MOQ: **5 pcs.**

JBMEDICO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

WIRE BASKETS AND WIRE PRODUCTS

Stainless steel



JB 161-00-00

Wire basket 6 liters, conical
with hook for T-slot.



JB 161-01-01

Wire basket 7 liters, conical
with hook for T-slot.



JB 161-02-01

Wire basket 15 liters, conical
with hook for T-slot.

Item no.: JB 161-00-00
Dimensions: Width x Depth x Height
231mm. x 161mm. x 242mm.
Material: Stainless steel (AISI 304)
Surface: Electropolished
Weight: 853 grams
Volumen: 6 liters
MOQ: 5 pcs.

Accessories for all 3 types of wire baskets:
Wall bracket with T-slot
Rail clamp with T-slot (DK)
Rail clamp with T-slot (EU)
Multibracket with T-slot
GCX rail bracket with T-spur

Item no.: JB 161-01-01
Dimensions: Width x Depth x Height
326mm. x 217mm. x 132mm.
Material: Stainless steel (AISI 304)
Surface: Electropolished
Weight: 667 grams
Volumen: 7 liters
MOQ: 5 pcs.

Item no.: JB 161-02-01
Dimensions: Width x Depth x Height
337mm. x 233mm. x 267mm.
Material: Stainless steel (AISI 304)
Surface: Electropolished
Weight: 959 grams
Volumen: 15 liters
MOQ: 5 pcs.

JBMEDICO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

WIRE BASKETS AND WIRE ACCESSORIES

Stainless steel



JB 161-03-00

Wire basket 25 liters, conical.

Item no.: JB 161-03-00
 Dimensions: Width x Depth x Height
 348 mm. x 248mm. x 400 mm.
 Material: Stainless steel (AISI 304)
 Surface: Electropolished
 Weight: 1.241 grams
 Volumen: 25 liters
 MOQ: 5 pcs.



JB 281-00-00

Wall bracket for hanging 25 liters wire basket.

Item no.: JB 281-00-00
 Dimensions: Width x Depth x Height
 23,5 mm. (above) 13,8 mm. (below) x 21 mm. x 53,5 mm.
 Material: Stainless steel (AISI 304)
 Surface: Electropolished
 Weight: 27 grams
 Hole size: Ø4,25 mm.
 MOQ: 10 pcs. in a bag



JB 162-00-00

Mask and tubing hook, with a guide for T-slot.

Item no.: JB 162-00-00
 Dimensions: Width x Depth x Height
 30 mm. x 76 mm. x 100 mm.
 Material: Stainless steel (AISI 304)
 Surface: Electropolished
 Weight: 80 grams
 MOQ: 10 pcs.



JB 146-45-38

Wall bracket for hanging heavy equipment.

Enables safe suspension of various hospital and care items.


Item no.: JB 146-45-38
 Dimensions: Width (above) or Width (below) x Depth x Height
 40 mm. (above) 20 mm. (below) x 50 mm. x 93 mm.
 Material: Stainless steel (AISI 304)
 Surface: Electropolished
 Weight: 117 grams
 Hole size: Ø6,4 mm.
 MOQ: 10 pcs. in a bag

JBMEDICO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

BRACKET FOR SUCTION CANISTER AND CATHETER CONTAINER

Stainless steel



JB 100-00-00

Bracket Ø100, for 2 liters
suction canister, Ø100 mm. hole.

Item no.: JB 100-00-00
Dimensions: Ø100 mm. / Ø 8 mm.
Material: Thread
Surface: Stainless steel (AISI 304)
Electropolished
Weight: 189 grams
MOQ: 5 pcs.



JB 239-00-00

2 liter catheter container
with guide for T-slot.

Item no.: JB 239-00-00
Height: 400 mm.
Material: Stainless steel (AISI 304)
Surface: Electropolished
Weight: 1.000 grams
Volumen: 2 liters
MOQ: 5 pcs.



JB 500-115-115

6,6 liter catheter container
with guide for T-slot.

Item no.: JB 500-115-115
Height: 200 mm
Material: Stainless steel (AISI 304)
Surface: Electropolished
Weight: 700 grams
Volumen: 6,6 liters
MOQ: 5 pcs.



JBMEDICO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

VARIOUS ACCESSORIES



JB 92-00-00

Hygienic universal urine bottle holder.

Fits all types of hospital beds and all common bottles. Machine washable and autoclavable. Optimised fit and easy to clean.

Item no.: **JB 92-00-00**
 Dimensions: **Width x Depth x Height**
123,6 mm. x 145 mm. x 192 mm.
 Material: **Stainless steel (AISI 304)**
 Surface: **Electropolished**
 Weight: **853 grams**
 MOQ: **5 pcs.**



JB 15-100-240

Storage box in stainless steel with T-track bracket.

Item no.: **JB 15-100-240**
 Material: **Stainless steel (AISI 304)**
 (Sling cleaned)
 Height: **100 mm.**
 Length: **240 mm. (without T-track)**
 Width: **100 mm.**
 Weight: **1.2 kg.**
 MOQ: **2 pcs.**



JB 115-00-00

Hospital journals holder in hard frosted plastic.

With handle for hanging on beds. Instrument washable up to 90° C. and fits all beds.

Item no.: **JB 115-00-00**
 Material: **Plastic PP (Polypropylene)**
 Weight: **800 grams**
 Colour: **White (Matte)**
 MOQ: **4 pcs.**



Depth: **56 mm. without handle**
 Depth: **98 mm. with handle**
 Weight: **282 mm.**
 Height: **372 mm.**

JBMEDICO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

VARIOUS ACCESSORIES **DANISH Hospital Power Cords**



DK Hospital Power Cords

Angled 0.5 meter, red DK hospital + straight C13 appliance connector and shielded 3G1.0 mm. ² cable	Item no.: 1210715
Angled 1 meter Red DK hospital + Straight C13 appliance connector and shielded 3G1.0 mm. ² cable	1190110
Angled 2 meters Red DK hospital + Straight C13 appliance connector and shielded 3G1,0 mm. ² cable	1190111
Angled 3 meters Red DK hospital + Straight C13 appliance connector and shielded 3G1.0 mm. ² cable	1190112
Angled 4.5 metres Red DK hospital + Straight C13 appliance connector and shielded 3G1.0 mm. ² cable	1190113
Angled 8 meters Red DK hospital + Straight C13 device connector and shielded 3G1,5 mm. ² cable	1211755

Angled 3,0 meters Red DK hospital + Straight C13 appliance connector and shielded 3G1.0 mm. ² cable	Item no.: 1212273
Angled 4,5 meters Red DK hospital + Straight C13 appliance connector and shielded 3G1.0 mm. ² cable	1212274



JBMEDICO

VARIOUS ACCESSORIES


Charging Towers or Power Towers



JB 680-400-15-S

Charging Towers /
Power Towers on wheels with Schuko plugs

Item no.: JB 680-400-15-S
Dimensions: 400 mm. (Bottom plate)
Wheels: 2 with brake/ 2 without brake
Max. load: 10A / 240V
Height: 762 mm
Width x Depth: 165 x 135 mm. (Tower)
Material: Steel
Surface: White painted / RAL 9016
Plug w/ground: 15 pcs. Schuko plugs
Weight: 13 kg.
Cord: 8 m.



JB 25-26-2025

Support tray for up to 12 phones

Item no.: JB 25-26-2025
Weight: 1.2 kg.
Height: 370 mm
Material: Steel
Color: RAL 9016
Surface: Powder Coated
Capacity: Up to 12 phones
MOQ: 2 pcs.



JBMEDICO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

VARIOUS ACCESSORIES

Toilet brush holder



JB 350-00-00

Toilet brush holder for most round and oval brushes in plastic.

Item no.: **JB 350-00-00**
 Dimensions: **Depth x Height 140 mm. x 153 mm.**
 Thickness: **5 mm.**
 Material: **Aluminium**
 Surface: **Anodized**
 Colour: **Metal**
 Weight: **173 grams**
 MOQ: **5 pcs.**



Accessories:
 Wall brackets
 Wall bracket with hex allen screw for locking

Item no. **JB 47-00-00**
JB 47-01-01

VARIOUS ACCESSORIES

Tubing holder



JB 292-00-00

IV tubing holder with a guide for t-slot.

Item no.: **JB 292-00-00**
 Dimensions: **Depth x Height x Width 71 mm. x 71 mm. 29 mm.**
 Material: **Aluminium**
 Surface: **Anodized**
 Colour: **Metal**
 Weight: **388 grams**
 MOQ: **10 pcs.**



Accessories:
 - Clamp with T-slot for DK 10 x 30 mm. standard rails.
 - Clamp with T-slot for EU-DIN 10 x 25 mm. standard rails
 - Wall bracket with T-slot

Item no. **JB 120-00-00**
JB 146-00-00
JB 47-00-00

JBMEDICO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

VARIOUS ACCESSORIES

Tubing holders



JB 600-00-01

Tube holder in red plastic

To be placed directly on the rail.
Fits DK standard rails only.

Item no.: **JB 600-00-01**
Dimensions: **Depth x Height**
140 mm. x 153 mm.
Material: **POM**
Colour: **Red**
Weight: **8,7 grams**
MOQ: **40 pcs.**





JB 600-10-25

Tube holder in red plastic

To be placed directly on the rail.
Fits EU DIN standard rails only.

Item no.: **JB 600-10-25**
Dimensions: **Depth x Height**
40 mm. x 60 mm.
Material: **POM**
Colour: **Red**
Weight: **8,7 grams**
MOQ: **40 pcs.**



JBMEDICO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

EMERGENCY PULL CORD + ACCESSORIES



JB IP 9025

Emergency pull cord 2,5 meters.



JB IP 24155

"Bell" for emergency pull cord in red PP-Plastic



JB IP 9004

Break joint, black, for emergency pull cord.



JB IP 3000-RED

Emergency pull cord on cable reel

Pull force: 60 N.m (newton meters)

Item no.: **JB IP 9025**
Material: **Plastic PE
LDPE soft Polyethylene**
Tensile force before
breaking: **2.2 kg.**
Weight: **9 grams (225 grams)**
Max. Tensile
Strength: **2,2 kg.**
MOQ: **25 stk.**

Item no.: **JB IP 4155**
Dimensions: **Diameter x Height
15,5 mm. (bottom) x 24 mm.**
Material: **PP, Plastic (Polypropylen))**
Weight: **2,5 grams (250 grams)**
MOQ: **100 stk. / pose**

Item no.: **JB IP 9004**
Dimensions: **Length: 22 mm.
Diameter: 9 mm.**
Material: **Plastic ABS**
Colour: **Black**
Tensile force before
breaking: **2.2 kg.**
Weight: **0,5 grams (100 grams)**
MOQ: **200 stk.**

Item no.: **JB IP 3000-RED**
Item no.: **JB IP 1000-RED**
Item no.: **JB IP 500-RED**
Item no.: **JB IP 100-RED**

JBMEDICO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

VARIOUS ACCESSORIES

Suction cups all 3 sizes fits our Sharps Container Brackets: JB 150-00-00



JB 06-08-230

Suction cup 230 ml conical
in stainless steel AISI 201

Item no. **JB 06-08-230**
Dimensions: **Height 80 mm.
Diameter bottom 60 mm.
Diameter top 75 mm.**
Material: **Stainless steel (AISI 201)**
Surface: **Polished**
Weight: **66 grams**
Volumen: **230 ml.**
MOQ:: **5 pcs.**



JB 06-08-330

Suction cup 330 ml conical
in stainless steel AISI 201

Item no. **JB 06-08-330**
Dimensions: **Height 110 mm.
Diameter bottom 59 mm.
Diameter top 79 mm.**
Material: **Stainless steel (AISI 201)**
Surface: **Polished**
Weight: **74 grams**
Volumen: **330 ml.**
MOQ:: **5 pcs.**



JB 06-08-500

Suction cup 500 ml conical
in stainless steel AISI 201

Item no. **JB 06-08-500**
Dimensions: **Height 120 mm.
Diameter bottom 61 mm.
Diameter top 79 mm.**
Material: **Stainless steel (AISI 201)**
Surface: **Polished**
Weight: **103 grams**
Volumen: **500 ml.**
MOQ:: **5 pcs.**



JBMEDICO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

TERMS OF SALE AND DELIVERY

JB-Medico ApS supplies private and public companies

Delivery:

Items usually arrive 1-5 business days after placing the order.

Products that need to be manufactured will of course have a longer delivery time.

Prices:

Unless otherwise stated, prices are pr. pcs. and in EURO

All government taxes are included our prices, with the exception of VAT.

All prices are pr. pcs.

Orders are taken at the prices applicable at the time of delivery.

We reserve the right to make adjustments and changes in prices, product specifications, tax changes, printing errors and image errors.

Payment:

Bank: Spar Nord Bank A/S

Account no.: 9004 / 1120732826

IBAN no: DK2390041120732826

SWIFT code: DABADKK

Transport damage:

Must be reported immediately upon receipt of goods, both to the carrier and to JB-Medico ApS.

A note must be made on the consignment note.

Custom-made goods and custom-made products, that have been specially taken home/manufactured for a customer, can not be returned.

Complaints:

Any error in the delivery must be reported immediately to customer service +45 7020 5826

Or info@jbmedico.com

Returns:

Upon receipt of goods from JB-Medico ApS, the customer must check that what is received corresponds correctly, to what was ordered and according to the enclosed delivery note.

Incorrect delivery must be reported immediately to JB-Medico ApS so that an exchange can take place.

All returns must be agreed with JB-Medico's dispatch department. Otherwise JB-Medico ApS may refuse receipt.

Returned goods must be returned in unbroken and undamaged original packaging.

In addition, the original packaging must be protected by an outer packaging during return transport. The goods must be cleaned.

Correspondence must be enclosed stating what is being returned, the reason for the return and with whom the return has been agreed.

JBMEDICO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361

FAVORITES

Item no.

[illegible][illegible]

NOTES

[illegible]



info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361