

MEDICAL RAIL CLAMPS - SLIDING CLAMPS
 Scandinavian standard for 10x30 mm. medical equipment rails



JB 88-00-00

Rail clamp, wide model, locked using two hex allen screws.



JB 103-00-00

Rail clamp, wide model, locked with one ball lock.



JB 104-00-00

Rail clamp, wide model, locked with two ball locks.



JB 207-00-00

Rail clamp, wide model, locked with two ball locks. It has three countersunk holes Ø6,6 mm.

“Prepared for suspension.”

Item no: JB 88-00-00
 Material: Aluminium (Natural Anodized)
 Dimensions: Length x Width x Depth:
 61.5 mm. x 90 mm. x 25 mm.
 Weight: 280 grams
 MOQ: 5 pcs.

Item no: JB 103-00-00
 Material: Aluminium (Natural Anodized)
 Dimensions: Length x Width x Depth:
 61.5 mm. x 110 mm. x 25 mm.
 Weight: 295 grams
 MOQ: 5 pcs.

Item no: JB 104-00-00
 Material: Aluminium (Natural Anodized)
 Dimensions: Length x Width x Depth:
 61,5 mm. x 110 mm. x 25 mm.
 Weight: 304 grams
 MOQ: 5 pcs.

Item no: JB 207-00-00
 Material: Aluminium (Natural Anodized)
 Dimensions: Length x Width x Depth:
 61,5 mm. x 110 mm. x 25 mm.
 Weight: 299 grams
 MOQ: 5 pcs.

CONTACT INFO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361