

**LAMP RAIL CLAMPS - EXTRA SOLID -  
FOR LONG AND HEAVY LAMPS  
EU DIN STANDARD FOR 10X25**

**SUSPENSION RAIL CLAMPS FOR  
DROP STANDS, IT AND MEDICAL EQUIPMENT  
EU DIN STANDARD FOR 10X25 MM.**



**JB 144-03-18**

Rail clamp, wide model, 2 pcs. ball locks and adapter fitting Ø18 mm.

Three bush sizes are available.  
Handles can be purchased (JB 94-00-00)

**JB 147-03-18**

Rail clamp, narrow model, 1 pcs. ball locks and adapter fitting Ø18 mm.

Three bush sizes available.  
Handles can be purchased (JB 94-00-00)

**JB 143-03-20**

Rail clamp, wide model, locked using two hex allen screws and adapter fitting Ø18 mm.

Handles can be purchased (JB 94-00-00)

**JB 145-03-20**

Rail clamp, wide model, with a one-ball lock, adapter fitting and brass bush Ø20 mm. hole.

Handles can be purchased (JB 94-00-00)

Item no.: JB 144-03-18  
Material: Aluminium (Ceramic Coating)  
Dimensions: Length x Width x Depth:  
61,5 mm. x 110 mm. x 64,2 mm.  
Weight: 502 grams  
MOQ: 2 pcs.

Item no.: JB 147-03-18  
Material: Aluminium (Ceramic Coating)  
Dimensions: Length x Width x Depth:  
61,5 mm. x 64 mm. x 64,2 mm.  
Weight: 343 grams  
MOQ: 4 pcs.

Item no.: JB 143-03-20  
Material: Aluminium (Ceramic Coating)  
Dimensions: Length x Width x Depth:  
61,5 mm. x 90 mm. x 60,6 mm.  
Weight: 536 grams  
MOQ: 2 pcs.

Item no.: JB 145-03-20  
Material: Aluminium (Ceramic Coating)  
Dimensions: Length x Width x Depth:  
61,5 mm. x 110 mm. x 64,2 mm.  
Weight: 550 grams  
MOQ: 2 pcs.

**CONTACT INFO**

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361