



JB 286-00-00

Wall or equipment fitting for two types of medical equipment rails.

Item no.: JB 286-00-00
 Material: Stainless steel (AISI 304)
 Colour: Metal
 Dimensions: Diameter x Depth
 Ø50 mm. x 10 mm.
 Weight: 137 grams
 MOQ: 10 pcs.



JB 400-00-00

Wall or equipment fitting for two types of medical equipment rails.

Item no.: JB 400-00-00
 Material: Aluminium
 Colour: Metal
 Dimensions: Diameter x Depth
 Ø50 mm. x 10 mm.
 Weight: 48 grams
 MOQ: 10 pcs.



JB 286-00-18 (small)
JB 286-00-36 (large)

Spacer to two types of wall rails.

Item no.: JB 286-00-18 (small)
 Material: Stainless steel (AISI 304)
 Colour: Metal
 Dimensions: Height / Outer diameter /
 Hole diameter
 18 mm. / Ø15 mm. / Ø6,5 mm.
 Weight: 20 grams
 MOQ: 40 pcs.

Item no.: JB 286-00-36 (large)
 Material: Stainless steel (AISI 304)
 Colour: Metal
 Dimensions: Height / Outer diameter /
 Hole diameter
 36 mm. / Ø15 mm. / Ø6,5 mm.
 Weight: 40 grams
 MOQ: 20 pcs.



JB 400-00-18 (small)
JB 400-00-36 (large)

Spacer to two types of wall rails.

Item no.: JB 400-00-18 (small)
 Material: Aluminium
 Colour: Metal
 Dimensions: Height / Outer diameter /
 Hole diameter
 18 mm. / Ø15 mm. / Ø6,5 mm.
 Weight: 7 grams
 MOQ: 40 pcs.

Item no.: JB 400-00-36 (large)
 Material: Aluminium
 Colour: Metal
 Dimensions: Height / Outer diameter /
 Hole diameter
 36 mm. / Ø15 mm. / Ø6,5 mm.
 Weight: 14 grams
 MOQ: 20 pcs.

CONTACT INFO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361