

## MONITOR ARMS AND ACCESSORIES

Junior Armstrong continued

## WALL BRACKETS

**JB**MEDICO



**JB 22-00-00**

Manual extension arm with Ø20 mm. shaft, concealed cable routing.

The extension arm can be mounted on e.g. Wall bracket, table-clamp, ceiling-tube, rail clamps and IV-Poles.

Item no.: JB 22-00-00  
 Material: Stainless steel (AISI 304)  
 Total length: 255 mm.  
 Height: 150 mm.  
 Range: 220 mm.  
 Weight: 1.100 grams

Accessories:  
 Item no.: JB 94-00-00



**JB 27-00-00**

Monitor bracket with Ø20 mm. shaft.

Can be potential equalised.

Item no.: JB 27-00-00  
 Material: Stainless steel (AISI 304)  
 VESA: 75 mm. x 75 mm.  
 VESA: 100 mm. x 100 mm.  
 Dimensions: Depth x Height x Width  
 95 mm. x 159 mm. x 114 mm.  
 Weight: 1.000 grams

Accessories:  
 Item no.: JB 46-00-00



**JB 31-00-00**

Circular wall bracket fits Ø20 mm. shaft.

Unbrako bolt is standard.

Item no.: JB 31-00-00  
 Material: Stainless steel (AISI 304)  
 Dimensions: Depth x Height x Width  
 58 mm. x 123 mm. x 132 mm.  
 Weight: 688 grams

Accessories:  
 Item no.: JB 94-00-00



**JB 241-00-00**

Extended wall bracket fits Ø20 mm. shaft.

Unbrako bolt is standard.

Item no.: JB 241-00-00  
 Material: Stainless steel (AISI 304)  
 Dimensions: Depth x Height x Width  
 120 mm. x 123 mm. x 130 mm.  
 Weight: 1.100 grams

Accessories:  
 Item no.: JB 94-00-00

## CONTACT INFO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361