



JB 90-00-00

Standard clamp desk mount for Ø20 mm. shaft.

Mounting on the edge of the table top.
Table thickness 19-51 mm.
Unbrako bolt is standard.

Item no.: JB 90-00-00
Material: Stainless steel (AISI 304)
Dimensions: Depth x Height x Width
80 mm. x 186,5 mm. x 60 mm.
Weight: 2.500 grams
Surface: Sand blasted
Bushing: Ø20 mm.
Cone height: 94 mm.

Accessories: Locking handle can
be purchased

Item no.: JB 94-00-00



JB 91-00-00

Through hole table mount for Ø20 mm. shaft.

Mounting in the table top.
Table thickness 0-41 mm.
Unbrako bolt is standard.

Item no.: JB 91-00-00
Material: Stainless steel (AISI 304)
Dimensions: Depth x Height x Width
80 mm. x 135.5 mm. x 60 mm.
Weight: 2.000 grams
Surface: Sand blasted
Bushing: Ø20 mm.
Cone height: 95 mm.

Accessories: Locking handle can
be purchased

Item no.: JB 94-00-00



JB 46-00-00

Locking handle M10x40 mm.

For e.g. screen laptops and keyboard brackets for
tilt/rotate function.

Item no.: JB 46-00-00
Material: Stainless steel (AISI 304)
MOQ: 2 pcs.



JB 49-00-00

Bracket for mounting IT equipment on anaesthesia
equipment or as wall bracket for Ø20 mm. tube.

Different models for different brands of
equipment, e.g. Dräger and Dameca.

Item no.: JB 49-00-00
Material: Stainless steel (AISI 304)
Dimensions: Depth x Height x Width
122 mm. x 162 mm. x 73 mm.
Weight: 1.000 grams
Surface: Sand blasted
Materiale Bushing: Brass

Accessories: Locking handle can be
purchased

Item no.: JB 94-00-00



JB 94-00-00

Locking handle M6x20 mm.

For e.g. clamping and manual extension arm.

Item no.: JB 94-00-00
Material: Stainless steel (AISI 304)
MOQ: 5 pcs.

CONTACT INFO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361