

ACCESSORIES FOR GLOVE BOX HOLDERS



JB 47-00-00

T-slot wall bracket for 1 pcs. glove box holder.

C-profile with 2 pcs. countersunk mounting holes and stop screw.

Item no.: JB 47-00-00
 Material: Aluminium
 Surface: Anodized
 Dimensions: Depth x Height x Width
 12,5 mm. x 60 mm. x 36 mm.
 Hole Size: Ø 6,4 mm.
 Weight: 36 grams
 MOQ: 10 pcs.



JB 47-00-02

T-slot wall bracket for 2 pcs. glove box holder.

C-profile with 4 pcs. countersunk mounting holes and stop screw.

Item no.: JB 47-00-02
 Material: Aluminium
 Surface: Anodized
 Dimensions: Depth x Height x Width
 12,5 mm. x 210 mm. x 36 mm.
 Hole Size: Ø 6,4 mm.
 Weight: 129 grams
 MOQ: 5 pcs.



JB 47-00-03

T-slot wall bracket for 3 pcs. glove box holder.

C-profile with 6 pcs. countersunk mounting holes and stop screw.

Item no.: JB 47-00-03
 Material: Aluminium
 Surface: Anodized
 Dimensions: Depth x Height x Width
 12,5 mm. x 360 mm. x 36 mm.
 Hole Size: Ø 6,4 mm.
 Weight: 222 grams
 MOQ: 5 pcs.



JB 47-00-04

T-slot wall bracket for 4 pcs. glove box holder.

C-profile with 6 pcs. countersunk mounting holes and stop screw.

Item no.: JB 47-00-04
 Material: Aluminium
 Surface: Anodized
 Dimensions: Depth x Height x Width
 12,5 mm. x 520 mm. x 36 mm.
 Hole Size: Ø 6,4 mm.
 Weight: 320 grams
 MOQ: 5 pcs.

CONTACT INFO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361