


# WIRE BASKETS AND WIRE PRODUCTS

Stainless steel



## JB 161-00-00

Wire basket 6 liters, conical with hook for T-slot.

## JB 161-01-01

Wire basket 7 liters, conical with hook for T-slot.

## JB 161-02-01

Wire basket 15 liters, conical with hook for T-slot.

Item no.: JB 161-00-00  
Dimensions: Width x Depth x Height  
231mm. x 161mm. x 242mm.  
Material: Stainless steel (AISI 304)  
Surface: Electropolished  
Weight: 853 grams  
Volumen: 6 liters  
MOQ: 5 pcs.

Accessories for all 3 types of wire baskets:  
Wall bracket with T-slot  
Rail clamp with T-slot (DK)  
Rail clamp with T-slot (EU)  
Multibracket with T-slot  
GCX rail bracket with T-spur

Item no.: JB 161-01-01  
Dimensions: Width x Depth x Height  
623mm. x 217mm. x 132mm.  
Material: Stainless steel (AISI 304)  
Surface: Electropolished  
Weight: 667 grams  
Volumen: 7 liters  
MOQ: 5 pcs.

Item no.: JB 161-02-01  
Dimensions: Width x Depth x Height  
337mm. x 233mm. x 267mm.  
Material: Stainless steel (AISI 304)  
Surface: Electropolished  
Weight: 959 grams  
Volumen: 15 liters  
MOQ: 5 pcs.

### CONTACT INFO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361