

MEDICAL RAIL CLAMPS FOR EQUIPMENT RAILS
EU DIN standard 10x25 mm.



JB 146-00-00

Rail clamp, narrow model, with one ball lock and T-slot.



JB 147-00-00

Rail clamp narrow model, with one ball lock.



JB 147-03-00

Rail clamp, narrow model, with one ball lock and three countersunk Ø6,6 holes.

“Prepared for suspension.”



JB 143-00-00

Rail clamp, wide model, locked using two hex allen screws.

Item no.: JB 146-00-00
Material: Aluminium (Ceramic Coating)
Dimensions: Height x Width x Depth:
61,5 mm. x 64 mm. x 28,6 mm.
Weight: 132 grams
MOQ: 5 pcs.

Item no.: JB 147-00-00
Material: Aluminium (Ceramic Coating)
Dimensions: Height x Width x Depth:
61,5 mm. x 64 mm. x 28,6 mm.
Weight: 168 grams
MOQ: 5 pcs.

Item no.: JB 147-03-00
Material: Aluminium (Ceramic Coating)
Dimensions: Height x Width x Depth:
61,5 mm. x 64 mm. x 28,6 mm.
Weight: 161 grams
MOQ: 5 pcs.

Item no.: JB 143-00-00
Material: Aluminium (Ceramic Coating)
Dimensions: Height x Width x Depth:
61,5 mm. x 90 mm. x 28,6 mm.
Weight: 308 grams
MOQ: 5 pcs.

CONTACT INFO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361