



**JB 158-03-18**

Multibracket with a fixing device, Ø18mm. hole.  
Fits Ø16-41 mm.

Three bush sizes available.  
Handle can be purchased separately.

Item no.: JB 158-03-18  
Material: Aluminium (Natural Anodized),  
SS 304, POM, Nylon (PA6)  
Dimensions: Length x Width x Depth:  
36 mm. x 51 mm. x 134/96 mm.  
Weight: 288 grams  
MOQ: 2 pcs.

ACCESSORIES: Handle can be purchased.



**JB 158-03-20**

Multibracket with a fixing device and brass bushing.  
Fits up to. Ø16-41 mm.

Handle can be purchased separately.

Item no.: JB 158-03-20  
Material: Aluminium (Natural Anodized),  
SS 304, Brass, POM, Nylon (PA6)  
Dimensions: Length x Width x Depth:  
36 mm. x 51 mm. x 134/96 mm.  
Weight: 440 grams  
MOQ: 2 pcs.

ACCESSORIES: Handle and T-slot fitting can be purchased.



**JB 163-00-01**

T-slot fitting for multi-bracket.

The holes are only suitable for multibracket.

Item no.: JB 163-00-01  
Material: Aluminium (Natural Anodized)  
Dimensions: Height x Width x Depth:  
60 mm. x 36,2 mm. x 12,5 mm.  
Weight: 36 grams  
MOQ: 1 bag of 20 pcs.

**CONTACT INFO**

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361