

# SUSPENSION HOUSING FOR RAIL CLAMPS AND MULTIBRACKET



**JB 205-00-00**

Fixing device for rail clamps and multi-brackets. Hole dimension: Ø25 mm.

Prepared for suspension.

Item no.: JB 205-00-00  
Material: Aluminium (Natural Anodized)  
Dimensions: Length x Width x Depth:  
36 mm. x 35,6 mm. x 61 mm.  
Weight: 133 grams  
MOQ: 2 pcs.

ACCESSORIES: Handle and bushing is purchased separatley.



**JB 205-00-18**

Fixing device for rail clamps and multi brackets. Hole dimension: Ø18 mm.

Prepared for suspension.

Item no.: JB 205-00-18  
Material: Aluminium (Natural Anodized)  
Dimensions: Length x Width x Depth:  
36 mm. x 35,6 mm. x 61 mm.  
Weight: 171 grams  
MOQ: 2 pcs.

ACCESSORIES: Handle and bushing is purchased separatley.



## CONTACT INFO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361