



JB 820-625-06

GCX rail fitting, prepared with mounting holes.

Item no.: JB 820-625-06
Material: Aluminium (Natural Anodized)
Colour: Metal
Dimensions: Height x Width x Depth:
82 mm. x 63 mm. x 6 mm.
Weight: 83 grams
MOQ: 4 pcs.

ACCESSORIES:
JB 163-00-01 T-slot fitting
JB 205-00-00 Adapter fitting/ Suspension housing Ø20
JB 205-00-18 Adapter fitting/ suspension housing Ø18



JB 820-625-163

GCX rail fitting with a T-slot.

Item no.: JB 820-625-163
Material: Aluminium (Natural Anodized)
Colour: Metal
Dimensions: Height x Width x Depth:
82 mm. x 36 mm. x 18,5 mm.
Weight: 123 grams
MOQ: 4 pcs.



JB 820-625-20

GCX rail fitting with suspension housing. Ø20 mm. bushing hole.

Item no.: JB 820-625-20
Material: Aluminium (Natural Anodized) and Brass
Colour: Metal
Dimensions: Height x Width x Depth:
82 mm. x 63 mm. x 42 mm.
Weight: 313 grams
MOQ: 2 pcs.

ACCESSORIES: Handle can be purchased.



POM bushing not included as standard

JB 820-625-18

GCX rail fitting with suspension housing Ø18 mm. hole.

Optional plastic bushing Ø9,6/13 or 16 mm.

Item no.: JB 820-625-18
Material: Aluminium (Natural Anodized) and Plastic
Colour: Metal
Dimensions: Height x Width x Depth:
82 mm. x 63 mm. x 42 mm.
Weight: 270 grams
MOQ: 2 pcs.

ACCESSORIES: Handle can be purchased.

CONTACT INFO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361