

INFUSION AND IV POLE - ACCESSORIES

Hooks



JB 204-00-00

Hook top - straight 2 strong hooks including bolt.



JB 211-00-00

Hook top - angled 2 strong hooks including bolt.



JB 212-00-00

Hook top - angled 3 strong hooks including bolt.



JB 190-00-08

Hook top - 4 strong hooks including bolt.

Item no.: JB 204-00-00
Material: Stainless steel (AISI 304)
Surface: Electropolished
Dimensions: Height x Width
406 mm. x 285 mm.
Hook distance: 252 mm.
Weight: 248 grams

Item no.: JB 211-00-00
Material: Stainless steel (AISI 304)
Surface: Electropolished
Dimensions: Height x Width
406 mm. x 285 mm.
Hook distance: 204 mm.
Weight: 248 grams

Item no.: JB 212-00-00
Material: Stainless steel (AISI 304)
Surface: Electropolished
Dimensions: Height x Width
406 mm. x 285 mm.
Hook distance: 120 mm.
Weight: 363 grams

Item no.: JB 190-00-08
Material: Stainless steel (AISI 304)
Surface: Electropolished
Dimensions: Height x Width
406 mm. x 285 mm.
Hook distance: 126 mm.
Weight: 448 grams

CONTACT INFO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361