

IV STAND / IV POLE - ACCESSORIES

Hooks and Wheels



JB 216-00-00

Hook top - 6 strong hooks including bolt.

Item no.: JB 216-00-00
Material: Stainless steel (AISI 304)
Surface: Electropolished
Dimensions: Height x Width
406 mm. x 285 mm.
Hook distance: 120 mm.
Weight: 596 grams



5387PJC100P30-13

Swivel castor, 100 mm. total brake.
Stainless (Levina), Grey.

Item no.: 5387PJC100P30-13
Wheel Ø: Ø13 mm.
Centre hole Ø: Ø13 mm.
Overall height: 137 mm.
Load: 100 kg. (dynamic)
Standard: EN 12530
Weight: 430 grams
Stainless steel ball bearing: Yes
Electrical conductivity: No
MOQ: 18 pcs.



5380PJC100P30-13

Swivel castor without brake 100 mm.
Stainless (Levina), Grey.

Item no.: 5380PJC100P30-13
Wheel Ø: Ø13 mm.
Centre hole Ø: Ø13 mm.
Overall height: 137 mm.
Load: 100 kg (dynamisk)
Standard: EN 12530
Weight: 430 grams
Stainless steel ball bearing: Yes
Electrical conductivity: No
MOQ: 20 pcs.



5380XSC100P30-13

Swivel castor without brake 100 mm.
Stainless (Levina) Grey, Conductive.

Item no.: 5380XSC100P30-13
Wheel Ø: Ø13 mm.
Centre hole Ø: Ø13 mm.
Overall height: 137 mm.
Load: 100 kg (dynamisk)
Standard: EN 12530
Weight: 430 grams
Stainless steel ball bearing: Yes
Electrical conductivity: Yes
MOQ: 20 pcs.

CONTACT INFO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361