

WIRE BASKETS AND WIRE ACCESSORIES

WALL BRACKET Max Ø38 mm. tube



JB 161-03-00

Wire basket 25 liters, conical.



JB 281-00-00

Wall bracket for hanging 25 liters wire basket.



JB 162-00-00

Mask and tubing hook, with a guide for T-slot.



JB 146-45-38

Wall bracket for hanging heavy equipment.

Enables safe suspension of various hospital and care items.

Item no.: JB 161-03-00
 Dimensions: Width x Depth x Height
 348 mm. x 248mm. x 400 mm.
 Material: Stainless steel (AISI 304)
 Surface: Electropolished
 Weight: 1.241 grams
 Volumen: 25 liters
 MOQ: 5 pcs.

Item no.: JB 281-00-00
 Dimensions: Width x Depth x Height
 23,5 mm. (above) 13,8 mm.
 (below) x 21 mm. x 53,5 mm.
 Material: Stainless steel (AISI 304)
 Surface: Electropolished
 Weight: 27 grams
 Hole size: Ø4,25 mm.
 MOQ: 40 pcs. in a bag

Item no.: JB 162-00-00
 Dimensions: Width x Depth x Height
 30 mm. x 76 mm. x 100 mm.
 Material: Stainless steel (AISI 304)
 Surface: Electropolished
 Weight: 80 grams
 MOQ: 10 pcs.

Item no.: JB 146-45-38
 Dimensions: Width (above) or Width
 (below) x Depth x Height
 40 mm. (above) 20 mm.
 (below) x 50 mm. x 93 mm.
 Material: Stainless steel (AISI 304)
 Surface: Electropolished
 Weight: 117 grams
 Hole size: Ø6 mm.
 MOQ: 10 pcs. in a bag

CONTACT INFO

info@jbmedico.com | +(45) 7020 5826 | www.jbmedico.com | CVR: 34212414 | JB-Medico ApS | Birkemosevej 3A, Hasselager, DK-8361